



## Ductless mobile fume hoods with modular filtration column

Designed to protect users during applications emitting vapors and/or chemical particles, the Captair® Flex® line offers a level of performance that ensures your safety while offering an environmentally-friendly alternative to traditional ducted systems.

Based on the Flex® technology -a flexible, adaptable modular filtration column- this line of chemical protection enclosures offers a wide range of possibilities and allows you to carry out a variety of applications in your laboratory.

The high containment and filtration performance of this technology offer users a high degree of protection, in accordance with the AFNOR NF X 15-211: 2009 standard, class 1 and class 2.

This technology is suited for many different industries, such as:  
**chemistry, pharmaceuticals, cosmetics, biochemistry,  
academics, petrochemistry, forensics, manufacturing,  
agro-food, hospitals, etc.**

*Distributed and Supported in the UK by:*



24 Norman Way Industrial Estate, Over, Cambridge, CB24 5WE  
Tel: +44(0)1954 233 100 Fax: +44(0)1954 233 101  
Email: [sales@camlab.co.uk](mailto:sales@camlab.co.uk) Web: [www.camlab.co.uk](http://www.camlab.co.uk)

Modular filtration technology  
customized for your applications

Air face velocity  
monitoring system

Electrical and fluid lines  
can be run into the enclosure

Vibration-absorbent  
work surface, high chemical  
and mechanical resistance




Automatic filter  
saturation detection



Bright, energy-efficient lighting

Ergonomically designed  
slanted front shield

Ergonomic openings



## Flex® Technology



**Commissioning,  
training, and  
maintenance**  
offered by our Asura  
department and its  
network of licensed technicians.



**M 32I**



Tested according to the  
ASHRAE 110: 1995 standard  
& compliant with the BS7989 standard  
Tests and markings **CE**

**Dimensions (mm)** M 32I & Midcap

	L	D	H mini/max
Interior	764	543	866
Exterior	800	630	1160/1345

### Technical specifications

	M 32I
Number of filtration columns	1
Number of fans (IP44)	1
Processed air flow	230 m³/h
Air velocity at openings (in on-position)	0,4 to 0,6 m/s
Voltage/frequency	90 - 264 V / 50-60 Hz

	M 32I
Including power used for lighting	70 W
Type of opening	Oblong
Structure	Anti-corrosion steel coated with 100% polyester
Panels	8 mm synthetic glass
Filtration module	Polypropylene

**M 39I**



**M 48I**



Flex<sup>®</sup> Technology



Except M 48I



Tested according to the  
ASHRAE 110: 1995 standard  
& compliant with the BS7989 standard  
Tests and markings **CE**

**Dimensions (mm)**

	M 39I			M 48I		
	L	D	H min x max	L	D	H min/max
Interior	965	522	860	1240	522	860
Exterior	1000	630	1160/1345	1275	630	1160/1345

### Technical specifications

	M 39I	M 48I
Number of filtration columns	1	
Number of fans (IP44)	1	
Processed air flow	230 m <sup>3</sup> /h	
Air velocity at openings (in on-position)	0,4 to 0,6 m/s	
Voltage/frequency	90 - 264 V / 50-60 Hz	

	M 39I	M 48I
Including power used for lighting (max)		70 W
Type of opening		Oblong
Structure	Anti-corrosion steel coated with 100% polyester	
Panels	6 mm synthetic glass	
Filtration module	Polypropylene	

## Flex® Technology



**XLS 392**



**XLS 483**



More info p.32



Tested according to the  
ASHRAE 110: 1995 standard  
& compliant with the BS7989 standard  
Tests and markings 

### Dimensions (mm)

	XLS 392			XLS 483		
	L	D	H mini/max	L	D	H mini/max
Interior	965	695	1040	1173	695	1040
Exterior	1000	800	1315/1495	1275	800	1315/1495

### Technical specifications

	XLS 392	XLS 483
Number of filtration columns	2	3
Number of fans (IP44)	2	3
Processed air flow	460 m <sup>3</sup> /h	690 m <sup>3</sup> /h
Air velocity at openings (in on-position)	0,4 to 0,6 m/s	
Voltage/frequency	90 - 264 V / 50-60 Hz	

### XLS 392 XLS 483

Including power used for lighting	121 W	191 W
Type of opening	Total	
Structure	Anti-corrosion steel coated with 100% polyester	
Panels	6 mm synthetic glass	
Filtration module	Polypropylene	





**XLS 633**



**XLS 714**



## Flex<sup>®</sup> Technology



Except XLS 633



Tested according to the  
ASHRAE 110: 1995 standard  
& compliant with the BS7989 standard  
Tests and markings **CE**

### Dimensions (mm)

	XLS 633			XLS 714		
	L	D	H mini/max	L	D	H mini/max
Interior	1566	695	1040	1765	695	1040
Exterior	1600	800	1315/1495	1800	800	1315/1495

### Technical specifications

	XLS 633	XLS 714
Number of filtration columns	3	4
Number of fans (IP44)	3	4
Processed air flow	690 m <sup>3</sup> /h	920 m <sup>3</sup> /h
Air velocity at openings (in on-position)	0,4 to 0,6 m/s	
Voltage/frequency	90 - 264 V / 50-60 Hz	

### XLS 633 XLS 714

Including power used for lighting	191 W	261 W
Type of opening	Trapezoid	Total
Structure	Anti-corrosion steel coated with 100% polyester	
Panels	6 mm synthetic glass	
Filtration module	Polypropylene	

### Control panel

#### Flow monitor

This device allows for continuous monitoring of the ventilation flow rate and alerts the user via visible and audible alarm in the event of a ventilation system failure.

#### Adjustable timer

This timer records the number of hours that the device has been in operation and, every 60 hours, notifies the user the need to test the saturation level of the molecular filter. (In accordance with the requirements of the AFNOR NF X 15-211: 2009 standard).

Digital display for optimal data read-out



### Sampling port



(Equipment not included on devices set with the Molecodet S automatic saturation detection sensor)

This port allows the user to sample the air within the detection chamber of the filtration module in order to evaluate the saturation level of the molecular filter, using color changing reagent tubes (not included).

### Anemometer



This system continuously monitors the face velocity, which must fall between 0.4 and 0.6 m/s. (in accordance with the requirements of the AFNOR NF X 15-211: 2009 standard).

### Chemical Listing

#### A guide of Erlab-approved chemicals

This guide includes a comprehensive list of chemicals that Erlab certifies as tested and authorized for use within the hood, under the conditions set forth by the AFNOR NF X 15-211: 2009 standard.

The guide includes almost 700 chemicals and lists the following for each of these chemicals: name of the chemical, its formula, its CAS number, its boiling point, its molecular mass, its saturation vapor pressure, the filter designed to trap this chemical and the retention capacity of this filter, the type of filter saturation detection system, the maximum mass of the chemical that may be introduced within the enclosure, and the name of the testing laboratory that performed the test related to this chemical handling.

### Energy-efficient lighting



**Internal Lighting**  
18 W - 500 lux  
- IP67.  
Compact fluorescent tube lights. One

to three tubes, depending on the model. Dust and vapor-tight. Even, bright lighting of the work surface.

### Side panel utility ports



Located on the enclosure sides, these ports allow electrical cables and/or fluid lines to enter the enclosure with ease.



The product of 40 years of R&D!

### Work surfaces

#### Glass work surface

- Tempered glass work surface with framing
- Ergonomic arm rest to work comfortably.



#### Phenolic resin work surface

- Work surface with built-in spill tray, made of phenolic resin, with an ergonomic arm rest to work comfortably.
- High chemical and mechanical resistance.
- Ideal for precision weighing operations.



#### Work surface in stainless steel 304 L

- High chemical and mechanical resistance. Rounded corners to facilitate cleaning operations.
- Built in spill tray.

(Only available for the models : M 321, M 391, XLS 483, XLS 714)



### Work benches and shelves

#### Mobicap™\*

- Metal rolling cart, equipped with 4 wheels (2 locking wheels).
- Allows the device to be moved safely.




\*Only available for the Captair® Flex® M 321 and Captair® Flex® M 391 models

#### Benchcap™

- Fixed metal work bench.
- Equipped with 4 height adjustment jacks.



#### Internal metal sliding shelf for Mobicap™ and Benchcap™.



### Molecode™ S



#### Large-spectrum filter saturation alarm.

(Equipment required by class I of the AFNOR NF X 15-211:2009 standard)

- 1 sensor is located in the detection chamber and automatically detects when the filter has become saturated by solvents.

- 1 sensor is in contact with the laboratory air and indicates an eventual pollution rise with solvents

### Particle pre-filter

Eliminates particles > 0.3 µm to optimize the performance of the HEPA H14 filter.



### Transparent back panel

- Made of synthetic glass.
- Offers 360° visibility of handlings performed into the enclosure
- Optimizes lighting conditions.



### Rear access panel

- Made of steel.
  - Located on the back side of the enclosure, this door provides easy access for large, heavy instruments.
  - Ideal for maintenance operations.
- (Except on the Captair® Flex® M 321 model)

