T1

**500 mg/l:** K<sup>+</sup>, Na<sup>+</sup>, Cl<sup>-</sup>

**100 mg/l:** Ag<sup>+</sup>

**50 mg/l:** Pb<sup>2+</sup>, Zn<sup>2+</sup>, Ni<sup>2+</sup>, Fe<sup>3+</sup>, Cd<sup>2+</sup>, Sn<sup>2+</sup>, Ca<sup>2+</sup>, Cu<sup>2+</sup>

**10 mg/l:** Co<sup>2+</sup>, Fe<sup>2+</sup>

**5 mg/l:** Cr<sup>6+</sup>

**Datatablel / Data table**

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GB

**LCK 339 Nitrate**

**Please check the "Edition Date" (see data table).**

**Safety advice and expiry date on package.**

**Principle**

Nitrate ions in solutions containing sulphuric and phosphoric acids react with 2.6-dimethylphenol to form 4-nitro-2.6-dimethylphenol.

**Range of Application**

Waste water (beware of interferences!), drinking water, raw water, surface water, soils, substrates, nutrient solutions

**Interferences**

The ions listed in **T1** have been individually checked up to the given concentrations and do not cause interference. We have not determined cumulative effects and the influence of other ions. High loads of oxidizable organic substances (COD) cause the reagent to change colour and to give high-bias results. The test can thus only be used for waste water analyses if the COD is less than 200 mg/l.

The measurement results must be subjected to plausibility checks (dilute and/or spike the sample).

**Removal of Interferences**

Nitrite concentrations of more than 2.0 mg/l interfere (high-bias results) and can be removed by the addition of a spatula-tipfull of amidosulphonic acid. The chloride can be precipitated out as silver chloride by adding silver sulphate.

High calcium concentrations cause turbidity. This interferes with the determination but can be prevented by adding a spatula-tipfull of EDTA to the sample.

**pH** sample .....3 – 10  
**Temperature** sample/reagents .....20 – 24°C

**In case of not working at the right recommended temperature an incorrect result may be obtained.**

Not more than 3 hours should elapse between sampling and analysis. **Store in a cool place!**

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