

GW4090 GLASSWARE WASHER



TECHNICAL MANUAL

CAREFULLY READ THE INSTRUCTION MANUAL

Failure to read or fully understand the instruction manual, or incorrect interpretation of the instructions herein may cause damage to the appliance as well as being a source of danger for the operator and lowering the performances provided by the machine to a considerable extent.

The manufacturer declines all liability for uses differing from those listed below.



The appliance must only be installed, serviced and repaired by authorized personnel.



The warranty provided will immediately become void if the machine is used in ways that FAIL TO CONFORM to the instructions given by the manufacturer.

This manual is for informational purposes only. The contents of this manual and the appliance described herein may be liable to modification without prior notice. In no case may SMEG be held liable for any direct or accidental damages deriving from or concerning the use of this manual.

MANUAL N°				
REV	ECR/ECN	DATE	DESCRIPTION	BY

Do you need information or Assistance ?

Please contact us from 8:00 to 18:00 hrs at the following numbers and addresses:

Camlab Ltd, Norman Way, Over, Cambridge CB24 5WE

Tel.: +44 (0)1954 233100 Fax: +44 (0)1954 233101

web: www.camlab.co.uk email: sales@camlab.co.uk

Our Sales Department staff will provide information about prices and offers. You can examine our entire product range in our Internet web site along with our innovations.

Our Technical Assistance Department can tell you anything you need to know about how to use your appliance in the proper way and can also put you in contact with your nearest authorized Assistance Centre if necessary.

International customers, please contact your local distributor.

TABLE OF CONTENTS

1. KEY TO THE SYMBOLS USED IN THE MANUAL AND ON THE MACHINE	7
2. INTRODUCTION	8
3. TECHNICAL SPECIFICATIONS	9
4. INSTALLATION	10
4.1 LIFTING AND TRANSPORT	10
4.2 POSITIONING	10
4.3 LEVELLING	10
4.4 WATER COCK CONNECTIONS	10
4.5 DEMINERALIZED WATER CONNECTION	13
4.6 OPERATION WITHOUT THE DEMINERALIZED WATER CONNECTION	13
4.7 DRAIN PIPE CONNECTION	14
4.8 ELECTRICAL CONNECTIONS	16
4.9 AFTER INSTALLATION	16
5. ADJUSTMENT OF THE WASHING PARAMETERS	17
5.1 CONNECTION TO THE WATER SUPPLY	19
5.2 EXTRA FILLING TIME	19
5.3 HOW TO FLUSH OUT THE WATER TRAP	20
5.4 REGENERATION	20
5.5 COOLER DRAINING	22
5.6 DRAINED WATER RECOVERY	22
5.7 CONDENSER ACTIVATION	22
5.8 PERISTALTIC DETERGENT PUMPS	23
5.9 DETERGENT FLOWMETERS	23
5.10 CONDUCTIVITY SENSOR	23
5.11 DWP ACTIVATION	23
5.12 TEMPERATURE PROBE SETTINGS	24

6.	UTILITIES	25
6.1.1	ABOUT	26
6.2	HOW TO CHANGE THE PASSWORD/S	26
6.3	SELECTION OF THE MENU LANGUAGES	26
6.4	CALENDAR	27
6.5	ARCHIVE	27
6.6	CONTINUOUS PRINTOUT OF THE CYCLES IN PROGRESS IN REAL TIME	27
6.7	ACTIVATION OF THE USER CONTROL OPTION	27
6.8	INSTALLATION OF THE MACHINE'S DEFAULT PARAMETERS	28
	The machine must be turned off when the physical connection via USB or RS232 is being carried out.	28
6.9	DOWNLOAD PROCEDURE OF THE UPDATED FIRMWARE VIA FRONT END	29
6.10	COMMUNICATION CONFIGURATION VIA LAN	30
6.11	THE MACHINE'S DEFAULT PARAMETERS	30
6.12	SERIAL NUMBER OF THE MACHINE	30
6.13	PRINTER SETTINGS	30
6.14	MAINTENANCE	31
6.15	STATUS OF THE OUTPUTS	32
6.16	STATUS OF THE INPUTS	33
7.	PROGRAMMING AND PROGRAM EDITING	35
7.1	PROGRAMMING MENU	36
7.2	MAIN PROGRAMMING MENU	36
7.3	HOW TO CREATE A NEW PROGRAM	37
7.3.1	How to select the program phases	38
7.3.2	How to edit each individual phase	39
7.4	HOW TO EDIT AN EXISTING PROGRAM	40
7.4.1	How to modify an existing program	40
7.4.2	How to copy an existing program	41
7.4.3	How to delete an existing program	41
8.	SUPPLEMENTARY MAINTENANCE	42
8.1	PUMP BLOCKED AFTER A LONG IDLE PERIOD	42

9.	TABLE OF PROGRAMS	43
10.	HYDRAULIC CIRCUIT	44
10.1	HYDRAULIC LAYOUT	
10.2	16.2 MAIN HYDRAULIC CIRCUIT (OVERALL VIEW)	44
10.2	16.2 MAIN HYDRAULIC CIRCUIT (OVERALL VIEW)	45
10.3	DETERGENT DISPENSING SYSTEM	46
10.4	DRYING SYSTEM	47
11.	WIRING DIAGRAMS	49
11.1	LOWER RACK IN DETERGENT COMPARTMENT (PLAN VIEW)	50
11.2	TERMINAL BOARD ASSEMBLY (DETAIL)	51
11.3	UPPER RACK IN DETERGENT COMPARTMENT	
11.4	ELECTRONIC BOARD IN THE UPPER RACK	53
12.	TROUBLESHOOTING – FAILURE SEARCH AND SOLUTIONS	57
12.1	PROBLEMS or FAILURES NOT ASSOCIATED WITH A FAULT	57
12.1.1	The machine does not turn on	57
12.1.2	Water leaks	57
12.1.3	Problems regarding drying	57
12.1.4	Problems with the quality of washing results	57
12.1.4.1	Arrangement of the instruments in the washing trolley	58
12.1.4.2	Quality and quantity of water	58
12.1.4.3	Water pressure	58
12.1.4.4	Washing temperature	58
12.1.4.5	Detergents used	58
12.1.4.6	Spots on the instruments	59
12.2	FAILURES ASSOCIATED WITH A FAULT	59
13.	INSTALLATION LAYOUTS	67
13.1	MODEL GW40900	67

1. KEY TO THE SYMBOLS USED IN THE MANUAL AND ON THE MACHINE



Read with the utmost care



Warning, danger



Warning, hot surfaces

2. INTRODUCTION

This manual is an integral part of the machine

It must be kept in a good condition and ready to hand for the entire life cycle of the machine.

You are advised to carefully read this manual and all the instructions it contains before using the appliance.

This appliance conforms to directives currently in force and to the applicable reference standards.

This appliance has been built for the following function:

- **Washing various kinds of glassware by means of Chemical Disinfection or Thermal disinfection;**
- **The appliance cannot be used for sterilizing the instruments or any other device.**

All other use is considered improper

The manufacturer declines all liability for uses differing from those indicated.



SMEG declines all liability for damage caused by washing, in the machine, instruments for which the relative manufacturers have not explicitly authorized automatic decontamination.

3. TECHNICAL SPECIFICATIONS

Electronic control	3 microprocessors (+1 on the optional LAN comm. Module)
Standard programs memorized	20
Programs that can be entered	10 (expandable up to 50)
Backlight LCD graphic display	128 x 64 pixels
Clock and date indicator	YES
Programmable phases	10
Phase parameters	Type of water, amount of detergent, target temperature, extension time in minutes, drying temperature and duration
Temperature inside tub	From 5°C to 95°C
Precision	0.1°C
Temperature sensors in tub	1 PT 1000 CLASS B IEC 60751
Time display	5 digits
Peristaltic pumps (0-220 ml)	4 (max)
Security lock	YES, with electromagnetic release
Alarm indications	Approx. 80
Troubleshooting menu	YES
Programme editing	YES (by means of password)
Password	4 levels
Languages	4: 4: Italian, English, French, German (on request: Spanish, Polish, Swedish, Russian, Japanese)
AUXILIARY FUNCTIONS	
Duct for inserting external probes	YES
Drain separating solenoid valve control	YES
RS232 printer output	YES
RS232 serial port for PC	YES
USB serial port (on demand)	YES
LAN connection (on demand)	YES
Cycle filing	YES
Cycle archive downloading	YES
DRYING SYSTEM	
Drying fan	0.4 kW
Heating elements for drying	2.5 kW
Class C 98% prefilter	YES
Class S 99.999% HEPA filter	YES
Temperature of the drying air	75°C-100°C
WATER SUPPLY (PRESSURE 1.5-5 BAR)	
Cold/hot mains water hardness	Max 42° F
Demineralized water conductivity	<20µS/cm
Pump for demineralized water	Optional
Built-in softener	YES
Recirculation pump	400L/MIN
WATER HEATING	
Electrical	6.3 kW MAX
Water preheating by boiler	Optional
STEAM CONDENSER	
	YES
DIMENSIONS LxDxH mm	
External (with built-in top)	900x630x850 (830)
Internal	520x515x545
Net weight (kg)	100
STEEL	
Washing tub	AISI 316L
External cladding	AISI 304
ELECTRIC POWER SUPPLY	
Max. power/voltage rating	1/N/PE 230V ~ 50Hz 2.8 kW 3/N/PE 400V ~ 50Hz 7.0 kW
NOISE LEVEL	
	50 dB
CONFORMITY	
	2006/95/EEC (low voltage equipment), 2004/108/EEC (emc) EN61010-1, EN61326

4. INSTALLATION



The appliance should be positioned against the wall (10 cm away at most) and must be installed by one of SMEG's authorized technicians, in compliance with the laws in force. The installer is responsible for ensuring that the appliance operates correctly after it has been installed. Moreover, he must provide all the information required for correct use of the appliance. Remove the scratchproof film from the external steel surfaces when the appliance is installed.



All adjustments, servicing work and so forth must be carried out with the appliance disconnected from the electricity main.

4.1 LIFTING AND TRANSPORT

Before leaving the factory, the base of the appliance is fixed to a pallet, which is also used for lifting and transporting the machine itself.

The appliance must be handled with a lift truck or transpallet.

4.2 POSITIONING

The accessory kit for assembly purposes (containing filters, seals and feet) is packed inside the washing tub.

The appliance can be fixed to the sides of the adjacent furniture. In this case, make sure that the steam vent at the rear is not obstructed in any way. It is therefore advisable for the wall at the back to be made of brickwork or impermeable material.

The heat must also be prevented from reaching any electrical circuits or sockets at the rear.

The appliance is equipped with pipes for supplying the water and draining it away. They can be positioned towards the right or left to allow the appliance to be installed properly.

The appliance can also be installed under a work top: this operation must be carried out by a specialized technician.

4.3 LEVELLING

Once the machine has been set in position, screw in the feet or unscrew them so as to regulate the height of the machine and level it until it is horizontal (the tolerance permitted is 2 degrees).

Good levelling will ensure that the machine operates in the correct way and, in models with dryers, will keep the doors flush.


4.4 WATER COCK CONNECTIONS

The machine is supplied with four water supply pipes that join to a non-return device inside the machine, in accordance with CENELEC HD27451 directives.

The pipes are designed for connection to cocks with $\frac{3}{4}$ " gas threaded bushings.

The following pipes must be connected to the water main:

- . cold water pipe;
- . hot water pipe;
- . pressurized demineralized water pipe (min. 1 bar max. 9 bar)
- . cold water pipe for the condenser.



If the connections are made with new pipes, it is advisable to allow the water to flow through in order to eliminate any rusty deposits or sludge.




fig. 2.3.2

WARNING

If the double hot / cold water supply is not available, the two black supply pipes (cold and hot) must be connected together by means of the relative Y connection (see fig. 2.3.2 circled figure).

Make sure that the pressure at which the mains water is supplied is within the operating limits: min. 1.5 bar max 5 bar.

The temperature of the hot water must not exceed 50°C.

Higher temperatures could impair the efficiency of the built-in water softener and damage the resins it contains.

The water supply cocks must be accessible.

Always shut off the water supply cocks when the machine is not operating.

WARNING

Chemical characteristics of the mains water incompatible with good washing actions

Washing glassware and stainless steel materials

If the supply water contains $\text{Fe}^{2+}/\text{Fe}^{3+}$ ions in an amount greater than 2 ppm and/or the supply water is harder than 45°F, the water must be pre-treated by installing a deferrization and/or water softening system to which the appliance must then be connected.

4.5 *DEMINERALIZED WATER CONNECTION*

The appliances are pre-engineered for connection to a pressurized demineralized water supply line at a minimum pressure of 1 bar and a maximum pressure of 5 bar.

Here again, install the relative cock in an accessible position.

Always shut off the water supply cocks when the machine is not operating.

4.6 *OPERATION WITHOUT THE DEMINERALIZED WATER CONNECTION*

Final rinsing with demineralized water is recommended as it completely eliminates the saline residues in the water and ensures optimum results.

However, good washing results can also be obtained without the use of demineralized water.

WARNING - IMPORTANT

If a demineralized water connection is not available, the settings of the water connections must be modified in the memory (see point 3.5 How to memorize the water connections).

WARNING

IF THE WATER SUPPLY FOR THE FINAL RINSES IS NOT MODIFIED, THE APPLIANCE WILL STOP AT THE FIRST DEMINERALIZED WATER RINSE AND THE FOLLOWING ERROR MESSAGE WILL APPEAR ON THE DISPLAY: NO WATER INLET

IF THIS HAPPENS, FIRST MAKE A RESET AND THEN MODIFY THE SETTINGS OF THE WATER CONNECTIONS. AFTER THIS, YOU WILL BE ABLE TO RE-ENTER THE PROGRAM.

4.7 DRAIN PIPE CONNECTION


Models with steam condensers

These models have 2 drain pipes 1/2 an inch in diameter. We strongly recommend to have two separate drains, in order to avoid interferences and backflows between the 2 drain hoses.

General instructions for installing the drain pipe

Use a drain pipe with trap. Take the following precautions when installing:

- Since the temperature of the drained water can reach as much as 93°C, the end of the drain pipe must be fixed in the hose nipple in a permanent way with a screw or spring clamp.



- There must be no sharp bends and consequent throttles in the drain pipe.
- The end of the drain pipe must be neither more than 80 cm higher nor 65 cm lower than the bearing surface of the appliance. In no case must the end of the pipe be immersed in water
- The inner diameter of the drain pipe must be at least 40 mm.
- Install a 50 mm diameter drain pipe if possible.
- The models with steam condensers have two 1/2" drain pipes: an adapter union is available in order to connect the 2 drain pipes straight to a 50 mm diameter drain pipe.
- Extensions to the drain pipe/pipes must be no longer than 1 meter and must possess the same inner diameter, while the maximum height at which the free end must be positioned must be reduced from 80 to 50 cm.

WARNING: The draining system must comply with the international standard. Our company declines all liability if improper use of the machine gives rise to pollution.

Models equipped with the drain separating device

Models fitted with valves to separate the polluted drained water (conveyed to a storage tank) from the clean drained water (conveyed to the normal sewer) must be installed in compliance with the instructions (supplied with the valves) for assembly of the separation valves.

4.8 ELECTRICAL CONNECTIONS

The appliances have been designed for connection to electricity mains with the following voltage ratings:

- 3 / N / PE ~ 400/230V 50Hz
- 1 / N / PE ~ 230V 50Hz

The appliances must be PERMANENTLY connected to the electricity main. There is a four-core disconnecter in the detergent compartment at the side.



All the metal surfaces of the machine must be grounded: during the installation phase, check to make sure that the connection to the ground terminal that protects the system is securely made.

The appliances leave the factory already pre-engineered for 3 / N / PE ~ 400/230V 50Hz voltage and are supplied with a 5x2.5 mm² power flex plus relative locking system or for 1 / N / PE ~ 230V 50Hz voltage for single-phase connection with a 1x2.5 mm² power flex and relative locking system.

Although it is possible, 230V single-phase connection is inadvisable since it considerably lowers the performances of the appliance when it comes to operating speed since the time it takes to heat the water is triplicated and the overall cycle time (considering the same cycle) is more than double that obtained with a threephase connection.

Power flex specifications:

- . FROR 5 x 2.5 mm², 450/750 V (threephase version)
- . FROR 3 x 2.5 mm², 450/750 V (single-phase version)



The electricity main to which the appliance is connected must comply with the standards in force. Smeg declines all liability for damage to persons or property caused by a defective electricity main.

4.9 AFTER INSTALLATION

After the machine has been connected to the electricity and water mains by the Smeg technician, certain operating parameters must be entered via the keyboard (or computer).

The parameters listed below must be entered and will be described in detail in the following chapters:

- Salt addition to the softener/1st regeneration start
- Language of the dialog menus
- Hardness of the water supply
- Type of water effectively connected
- Volume of the liquid detergents to dispense (if Smeg's default programs are changed).

Before using the machine, make sure that it has been correctly connected to the electricity main, that the water fill pipes have been connected to the cocks and that these are open, and that the drain pipe has been positioned as indicated.

5. ADJUSTMENT OF THE WASHING PARAMETERS

The basic machine parameters can be entered via the console (the operator must have a superuser password level) or through the WDTRACE management/control software.

Access the initial mask

```
PR: PLASTIC 70°C
      (000003)
N° Phases: 5
Duration: 1:10:00
      <Dryer ON>

05/02/06      08:30:17
```

then press ESC to access the program selection menu:

```
<PREWASH>
<DEMI PREWASH>
<PLASTIC 70°C>
<BASES 70°C>
<BASES 80°C>
<NURS.BOT. 93°C>
```

press ESC again to reach the menu where the various sub-menus can be selected through the superuser password:

```
ENTER PASSWORD
-----
```

Then, by pressing ENTER, you reach the following submenu (use arrows to scroll):

```
<PROGRAMS>
<PROGRAMMING>
<UTILITIES>
<WASH PARAMETERS>
```

Press ENTER on the “WASHING PARAM” item to access the menu with the basic settings, which include the 12 items listed below:

```
<Water Connection>
<Extra Filling Time>
<Water trap Washing>
<Regeneration>
```

<Drain mixture>
<Drained Water Recovery>
<Activate Condenser>
<Detergent Dispensers>

<Deterg Flowmeter>
<Conduct. Sensor>
<DWP Activation>
<Temp Calibration>

The meanings of these settings will now be described one by one.

Saving the changes



Every time parameters are changed, the modifications are saved in two cases:

- 1) if the confirmation key (ENTER) can be pressed on the modification menu, the changes become permanent immediately. Use the WDTRACE® management software to verify this instantly.
- 2) If the changes cannot be confirmed by pressing the ENTER key on the menu, the changes will become permanent only after the machine is turned off and turned on again.

5.1 CONNECTION TO THE WATER SUPPLY

The washing programs have been formulated so as to ensure that the most suitable type of water (cold, hot or demineralized) is used for the required function.

If all three types of water are not available, the washing program can be corrected automatically so as to use the most appropriate alternative.

Select this item to access the menu where the type of water is selected.

Cold Water	<Yes>
Hot Water	<Yes>
Demi Water	<Yes>
Softened Water	<Yes>

Use the **+** and **-** keys to select or deselect the type of water. Once the selection has been made, press ENTER to confirm it.

If softened mains water is available, the built-in softener must be inhibited from the standard operating mode and salt must not be added to the softener's reservoir (see 3.2.1).

5.2 EXTRA FILLING TIME

(the technical level password is necessary)

Extra Filling Sec.
< 000 >

If the pressure of the water in the supply main is below the minimum value, the filling time of the solenoid valve on the inlet can be increased so as to prevent an alarm from appearing. Enter the required time in seconds (using the **+** and **-** keys). However, it is advisable to restore the right pressure in the water supply system if the times and costs of each cycle are to be optimized.

5.3 HOW TO FLUSH OUT THE WATER TRAP

Flushing out the water trap is a function that rinses the internal pump and hydraulic circuits between one phase and the next. This completely removes the additive used for the previous wash and/or eliminates any mains water residues when rinses with demineralized water are required.

The “water trap flushing” function:

- improves the quality of the washing action

but:

- increases the amount of water used
- lengthens the washing cycle by several minutes

Users can therefore choose whether to program the function or to inhibit it in the following way:

Water Trap Wash

< NO >

5.4 REGENERATION

HARDNESS IN °F <XX>

HARDNESS IN °T <XX>

HARDNESS IN °I <XX>

NUMBER PHASES <XX>

If the user himself enters the number of phases after which regeneration must take place, this value will be given priority over the other items.

Conversion table			
Unit of measurement	°F	°T	°I
1°F	1	0,56	0,7
1°T	1,79	1	1,25
1°I	1,43	0,8	1

The following table gives an initial estimation as to the correspondence between degrees of hardness and number of phases to conduct before regeneration need be done.

Hardness in French degrees [°F]	Hardness in German degrees [°T]	N° Phases
<10	5.6	16
15	8.4	15
20	11.2	14
25	14	14
30	16.8	11
35	19.7	8
40	22.5	6
45	25	5
50	28	4

5.5 COOLER DRAINING

Drain Water Temp <degrees> 0 = NO <00>
--

A minimum temperature of about 75°C can be reached by selecting this option, i.e. by entering a temperature with the + and - keys.

5.6 DRAINED WATER RECOVERY

(the technical level password is necessary)

Recover Drain Water <NO>

Enter “Yes“ if a double draining valve is installed, so that the liquid is transferred to the appropriate drain (for example: a recovery or filter drain) .

5.7 CONDENSER ACTIVATION

(the technical level password is necessary)

Condenser On Temp <degrees> 0 = NO <65>

The steam condenser can be inhibited by entering a temperature of 0°C. If a temperature is entered, that will be the temperature at which the condenser starts operating.

5.8 PERISTALTIC DETERGENT PUMPS

(the technical level password is necessary)

P1 Deter <Yes>
P2 Neutr <Yes>
P3 Disinf <Yes>
P1 Lubri <Yes>

P5 Det N <Yes>

The peristaltic pumps can be activated in these two masks.
Up to 5 pumps can be installed and, thus, activated.

5.9 DETERGENT FLOWMETERS

(the technical level password is necessary)

Flowmeter P1 <Yes>
Flowmeter P2 <Yes>
Flowmeter P3 <Yes>
Flowmeter P4 <Yes>

Flowmeter P5 <Yes>

Each pump can either have or not have its own flowmeter. The relative flowmeter can be selected (or not selected) in this mask.

5.10 CONDUCTIVITY SENSOR

(the technical level password is necessary)

Sensor installed <NO>
Activate meas. k <NO>

The conductivity sensor (optional) checks to make sure that detergent has actually been dispensed. This optional is activated in this mask.

5.11 DWP ACTIVATION

(the technical level password is necessary)

DWP activation
<Yes>

Smeg allows an auxiliary pump ("DWP" pump for demineralized water) to be activated if the pressure of the demineralized water is less than 1 bar.

5.12 TEMPERATURE PROBE SETTINGS

(the technical level password is necessary)

Probe TL1 <0>	Probe TA2 <0>
Probe TL2 <0>	Probe TB <0>
Probe TCL <0>	
Probe TA1 <0>	

As time goes by, the settings of the temperature probes may need to be adjusted. The tenths of a degree required to raise or lower the settings of each individual probe can be entered in this mask. Use the '+' and '-' keys to make the adjustments.

TL1 = probe that monitors the temperature in the washing chamber

TL2 = auxiliary probe (optional) probe that monitors the temperature in the washing chamber

TCL = probe that independently monitors the temperature in the washing chamber (not installed in this model)

TA1 = probe that monitors the temperature of the drying air

TA2 = probe that monitors the temperature of the drying air (not installed in this model)

TB = probe that monitors the temperature of the external boiler (not installed in this model)

6. UTILITIES

The basic parameters concern the washing cycle, but there is another series of parameters concerning the general machine settings and that can be selected by means of various menus.

Access the initial mask

PR: PLASTIC 70°C (000003)
N° Phases: 5
Duration: 1:10:00 <Dryer ON>
05/02/06 08:30:17

then press ESC to access the program selection menu:

<PREWASH>
<DEMI PREWASH >
<PLASTIC 70°C>
<BASES 70°C>
<BASES 80°C>
<NURS.BOT. 93°C>

press ESC again to reach the menu where the various sub-menus can be selected:

<PROGRAMS>
<PROGRAMMING>
<UTILITIES>
<WASH PARAMETERS>

Press ENTER on the “UTILITIES” item to access the menu with the basic settings, which include the 12 items listed below:

<About>	(can only be accessed with a SUPERUSER password)
<Change PSW>	(can also be accessed with the user password. Can only be changed with the operator's personal password)
<Select Language>	(can only be accessed with a SUPERUSER password)
<Calendar>	(can only be accessed with a SUPERUSER password)
<Archive>	(can only be accessed with a SUPERUSER password)
<Print Run Time>	(can only be accessed with a SUPERUSER password)
<User Control >	(can only be accessed with a SUPERUSER password)
<Communication>	(can only be accessed with a SUPERUSER password)
<Config LAN>	(can only be accessed with a SUPERUSER password)
<Default Values>	(can only be accessed with a SUPERUSER password)
<Machine Code>	(can only be accessed with a SUPERUSER password)
<Printer Setup>	(can only be accessed with a SUPERUSER password)
<Maintenance>	(can only be accessed with a TECHNICAL password)
<Output Status>	(can only be accessed with a SUPERUSER password)
<Input Status>	(can only be accessed with a SUPERUSER password)

6.1.1 ABOUT

The software version installed on the machine can be read in this menu. The information regarding the three microprocessors can also be read: the master and slave on the motherboard and the console on the console card. The boot version of the system is provided as well (a low level routine is loaded onto each microprocessor).

```
Master (xx.xx)
      ver.yy.yy
Slave (xx.xx)
      ver.yy.yy
Console (xx.xx)
      ver.yy.yy
```

xx.xx indicates the relative boot version
yy.yy indicates the software version

6.2 HOW TO CHANGE THE PASSWORD/S

```
<Change User >

<USER 1>

<USER 2>

<USER 3>
```

```
<USER 4 >

<USER 5>

<SUPERUS>

<TECHNICIAN>
```

The passwords can be changed by means of these two masks; only a personal password or a lower level password can be changed (if authorization is provided)

Press ENTER after having gone to the selected item. The following window will appear:

```
Enter Password
-----
```

If, for example, <USER 1> is selected, the screen appears as shown to the side. After writing the new password, confirm it by pressing the ENTER key

6.3 SELECTION OF THE MENU LANGUAGES

```
<ITALIANO>

<ENGLISH>
```

6.4 **CALENDAR**

Yy/mm/dd hh:mm:ss
06/02/06 19:05:02

The machine is equipped with a clock/date indicator. The current time and date can be entered in this menu.

6.5 **ARCHIVE**

Overwrite
<Yes>

The memory in which all the cycle data are stored is of limited entity. If “YES” is entered in this menu, the machine will continue to file the data even when there is no more space in the memory.
If “NO” is entered, the relative warning message will appear on the display: the data must now be downloaded by linking to the serial port.

6.6 **CONTINUOUS PRINTOUT OF THE CYCLES IN PROGRESS IN REAL TIME**

Print run time
<NO>

Select this option and the most significant data will be printed on paper as the cycle advances.

6.7 **ACTIVATION OF THE USER CONTROL OPTION**

User Control
<NO>

If this option is selected, the program will ask for the password before beginning a cycle.
You are advised to disable this option unless there are effective security problems and you need to exercise rigid controls over the operators who use the machine.

6.8 INSTALLATION OF THE MACHINE'S DEFAULT PARAMETERS

MODALITA' <LOC>

TIPO <COM>

The control system can be accessed directly via the serial port located under the front panel (GW4090 model) or in the rear of the machine (GW3060).

LOC mode must be used when the need arises to update the software (for example, when a new version of the control system software is released by Smeg). In this case (operation must only be performed by an authorized technician), LOC must be set and select COM TYPE.

If communication is through an optional card installed at the rear of the machine (and therefore by using the SMEG WDTRACE® software), the MODE must be selected as REM (or switch off the machine, if the connection was through a serial port) and TYPE according to the following:

- . select USB when the connection is via USB;
- . select COM when the connection is via RS-232 or via LAN.

The software can also be updated via a LAN connection.



The software updating is a critic procedure even though WDTRACE® allows to do it easily. Make sure to update the software very carefully.

In no case may SMEG be held liable if this procedure is executed incorrectly or files other than those sent by Smeg are downloaded into the control system.



The machine must be turned off when the physical connection via USB or RS232 is being carried out.

In no case may SMEG be held liable if this procedure is executed incorrectly and the data archive goes lost or more serious damages occur.

6.9 DOWNLOAD PROCEDURE OF THE UPDATED FIRMWARE VIA FRONT END

Under the front panel there is a serial connection through which it is possible to download the new firmware or check the machine status.

1. Connect your PC to the RS-232 port;



2. Turn ON the machine and wait for a few seconds.
3. Launch the program WDTRACE® from the PC and follow the steps shown.

6.10 COMMUNICATION CONFIGURATION VIA LAN

IP ADDRESS
0000 0000 0000 0000

Subnet Mask
0000 0000 0000 0000

IP Gateway
0000 0000 0000 0000

The network data for the machine can be set (once the WD-LAN card is installed) on this screen (LAN config.). To confirm, press ENTER and then ESCAPE to save the setting. However, the machine must also be shut off and turned on again. The network settings configuration is only permitted using the keyboard. It cannot be done using WDTRACE®.

6.11 THE MACHINE'S DEFAULT PARAMETERS

Enter default parameters and values

<NO>

Select "Yes" and then proceed with updating, and the operating parameters memorized in the factory will be restored

6.12 SERIAL NUMBER OF THE MACHINE

Serial number

The machine's serial number is indicated in this mask.

6.13 PRINTER SETTINGS

Local <Yes>

Remote <No>

The machine can either have a printer on board (other models) or can be connected to a table printer (GW3060/GW4090 model); in this last case, REM must be set.

6.14 MAINTENANCE

Maintenance

<000000>

<No>

Change filter

<0000000>

<No>

This mask shows the number of cycles accomplished by the machine and can also be used to enter a reminder as to when the Technical Assistance Center must be called for the inspections required.

A filter replacement reminder can also be activated for the air filter of the drying system.



Saving the changes

Every time parameters are changed, the modifications are saved in two cases:

If the confirmation key (ENTER) can be pressed on the modification menu, the changes become permanent immediately. Use the WDTRACE management software to verify this

Instantly

If the changes cannot be confirmed by pressing the ENTER key on the menu, the changes will become permanent only after the machine is turned off and turned on again.

6.15 STATUS OF THE OUTPUTS

00:0	07:0	14:0	21:0
01:0	08:0	15:0	22:0
02:0	09:0	16:0	23:0
03:0	10:0	17:0	
04:0	11:0	18:0	
05:0	12:0	19:0	
06:0	13:0	20:0	

All outputs can be activated by means of the “+” key. Move with the “**rh arrow**” and “**lf arrow**” keys to select the required output.

Note: If the operator exits the screen with an output activated, the output remains this way!

- .00 = (EK) EVCC condenser filling solenoid valve
- .01 = (EA) EVSA alternative draining solenoid valve (optional)
- .02 = (ER) EVR regeneration solenoid valve
- .03 = (UC) EVUC condenser outlet solenoid valve
- .04 = (EU) EVU outlet solenoid valve (not available in this model)
- .05 = (ED) EVD demineralized water solenoid valve
- .06 = (EC) EVC hot water solenoid valve
- .07 = (EF) EVF cold water solenoid valve
- .08 = (BP) BPE1 door lock 1 solenoid valve
- .09 = (SC) PSC condenser drain pump
- .10 = (A1) MA1 drying motor 1
- .11 = (PD) MD demineralized water pump motor (optional)
- .12 = MS drain pump motor
- .13 = P5 peristaltic pump P5
- .14 = P4 peristaltic pump P4
- .15 = P3 peristaltic pump P3
- .16 = P2 peristaltic pump P2
- .17 = P1 peristaltic pump P1
- .18 = (RA) RA1 drying heating element 1
- .19 = (M2) ML2 washing motor 2
- .20 = (M1) ML1 washing motor 1
- .21 = (R3) CR3 heating element 3 relay actuation
- .22 = (R2) CR2 heating element 2 relay actuation
- .23 = (R1) CR1 heating element 1 relay actuation

6.16 STATUS OF THE INPUTS

00:0	07:0	14:0	21:0
01:0	08:0	15:0	22:0
02:0	09:0	16:0	23:0
03:0	10:0	17:0	
04:0	11:0	18:0	
05:0	12:0	19:0	
06:0	13:0	20:0	

1 will appear instead of 0 if an INPUT is activated.

It may be useful to consult this mask for checking the status of the inputs after an alarm or “abnormal” situation.

- .00 = N.o. (door micro 2)
- .01 = (M1) MCM1 door lock 1 mechanical opening microswitch
- .02 = (C1) MCP1 door lock 1 electrical opening microswitch
- .03 = (AS) MCA acquastop microswitch (not installed in this model)
- .04 = (H1) PAP1 washing pump 1 motor high pressure switch
- .05 = (M2) MCM2 door lock 2 mechanical opening microswitch (not installed in this model)
- .06 = PS safety level pressure switch
- .07 = PL washing level pressure switch
- .08 = SS softener salt sensor
- .09 = T1 temperature sensor
- .10 = T2 temperature sensor (optional)
- .11 = TL washing temperature sensor (T3)
- .12 = TB boiler temperature sensor (T4) (optional)
- .13 = (H2) PAP2 washing pump 2 motor high pressure sensor (not installed in this model)
- .14 = (P4) SP4 can 4 level sensor
- .15 = (SF) SFS flow sensor on drain (not installed in this model)
- .16 = (P3) SP3 can 3 level sensor
- .17 = (SC) SSC condenser safety level sensor
- .18 = (P2) SP2 can 2 level sensor
- .19 = (LC) SSL condenser work level sensor
- .20 = (P1) SP1 can 1 level sensor
- .21 = (P5) SP5 can 5 level sensor
- .22 = IR INTERF_REM
- .23 = IL INTERF_LOC

The state of the inputs and outputs is shown in the figure below when the machine is closed and in stand-by mode.

Stato Macchina

1 PRELAV

Buffer stato Macchina

Stato Lavaggio 1- RIPOSO

Stato Fase 1- RIPOSO

Fase Corrente: 1

Tempo estensione (sec):

Tempo Progressivo(sec): 0:00:00

Qta totale acqua caricata (dl) 0

Programma Abortito

Asciugatura non terminata

Contatore Lavaggi 1

Alarmi

Stato Disabilitazione Allarme

1- NO_RISC_ACQUA

4- TL1 SOVRAT

5- TL1_APERTA

7- TA1_APERTA

11- MANCA_Acq_FREDDA

12- MANCA_Acq_CALDA

13- MANCA_Acq_DEMI

14- ...

15- PRESSIONE_Acq_CALDA

16- PRESSIONE_Acq_DEMI

17- TMO_CARICO_FREDDA

18- TMO_CARICO_CALDA

19- TMO_CARICO_DEMI

20- AVARIA_CARICO_ACQUA

21- COLL_FM_INVERTITI

22- AVARIA_FM_FREDDA

23- ACQ_INSUFF

24- NO_Acq

25- SCHIUMA_IN_VASCA

26- PERDITA_EVF

27- PERDITA_EVC

28- PERDITA_EVD

29- NO_SCARICO

30- LIVELLO_SICUREZZA

31- AVARIA_LIV_SICUREZZA

32- ACQUA_IN_VASCA

33- NO_ACQUA_CV

34- NO_SCARICO_CV

Output

H1-EVF

H11-RA1

H3-EVD

H4-EVU

H5-P3

H6-EVR

H7-EVSA

H8-EVCC

H12-P4

H14-MS

H15-P2

H16-MA1

H17-P1

H26-PSC

L7-BPE1

TEST_AD1

up-PS_LED_1

L8-BPE2

TEST_AD2

H9-CR1

H10-ML1

H13-RELE K1

H18-MD

H19-

H20-

H23-

H24-

H25-

H9-CR2

H9-CR3

H27-Sol1

CMD-RST

CMD-REM

LED-1

Input

in220_1

116-MCA

117-MCP1

118-MCM1

119-MCP2

120-MCM2

121-PL

122-PS

123-PAP1

124-PAP2

128-TRA1

129-TRA2

130-TRL

131-TRB

11-SS

12-SP1

13-SLC

14-SP2

15-SSC

16-SP3

17-SFS

18-SP4

19-SP5

INTERF_REM

INTERF_LOC

69- TANICA_VUOTA_P2

70- TANICA_VUOTA_P3

71- TANICA_VUOTA_P4

72- TANICA_VUOTA_P5

73- ARCH_ERR

74- PERDITA_VASCA

75- MANCA_SALE

77- TEMP_OVER_45

78- RESTORE_FAILS

79- PGM_NON_CONG

80- AVARIA_FM_CALDA

81- AVARIA_FM_DEMI

82- AVARIA_SOLENODEI

84- TASC1_SOVRAT

86- P1_BLOCCATA

87- P2_BLOCCATA

88- P3_BLOCCATA

89- P4_BLOCCATA

91- ARCH_FULL

92- SERV_FILTRO

93- SERV_MANUT

94- TEMP_NON_STAB

35- AVARIA_SLC_CV

36- AVARIA_PS_CV

37- NO_TEMP_SC_MISC

38- AVARIA_SC_MISC

39- MISURA_SK

41- NO_IMMIS_P1

42- NO_IMMIS_P2

43- NO_IMMIS_P3

44- NO_IMMIS_P4

46- OSTRUZ_P1

47- OSTRUZ_P2

48- OSTRUZ_P3

49- OSTRUZ_P4

51- COLLEG_FM_DET

52- PORTA1_APERTA_EL

53- PORTA2_APERTA_EL

54- PORTA1_APERTA_MEC

55- PORTA2_APERTA_MEC

56- AVARIA_BLOCCOPORTA1

57- AVARIA_BLOCCOPORTA2

58- NO_RISC_ASC1

59- NO_RISC_ASC2

67- NO_RAFFRED_ASC

68- TANICA_VUOTA_P1

Note: In this case I1-SS is flagged because of the lack of salt (that therefore must be added).

7. PROGRAMMING AND PROGRAM EDITING

SMEG's GW4090 instrument washer has an extremely flexible programmable control that allows highly sophisticated washing cycles to be programmed and all the parameters that affect the quality and type of washing process to be entered.

Thanks to this flexible system, users can enter the washing programs that suit their requirements and that take the least possible time.

The machine's memory contains 20 programs that should cover the washing requirements of hospitals.

We advise users to process the materials with one of the basic programs, which should be chosen to suit the case in question.

If it is necessary to edit an existing program or create a new one, we advise you to only modify one parameter at a time (e.g. to shorten the washing time or change the washing temperature) and to then test the program with dirty materials in order to find out whether there have been any improvements.

If several parameters are changed at the same time, it often becomes difficult to ascertain which parameter has worsened or improved the quality of the wash obtained.

There are essentially two types of programs:

1. SMEG programs: only the percentage of detergents and drying parameters can be changed in these programs.
2. NON-SMEG programs: all the data can be changed in these programs.

7.1 PROGRAMMING MENU

There are 20 programs in the memory of the microprocessor (20 for operating theaters).

These programs reside in a permanent memory.

20 programs (operating theater) are installed in the operating memory (similar to the hard disk of a personal computer) by means of the installation procedure described below.

Besides the 20 positions filled by standard programs, this operating memory also contains further vacant positions where another 10 programs can be memorized.

In order to be edited or executed, the washing programs must be transferred to a work memory (similar to the RAM of a PC) one at a time.

This operation takes place automatically.

7.2 MAIN PROGRAMMING MENU

<LOAD PROG>
<SELECT PROG>

Select "LOAD PROG" to access the following window:

<HOSPITAL>
<DENTAL SURGERY >
<LABORATORY>

from which Smeg's programs for various purposes can be "loaded" (only "HOSPITAL" programs can be selected for this model, GW4090).

Select "SELECT PROG" to access the window below:

<NEW>
<PREWASH>
<DEMI PREWASH >
<PLASTIC 70°C>

A new program can be created or an existing one edited in this mask, depending on the program selected.

7.3 HOW TO CREATE A NEW PROGRAM

If “NEW” has been selected in the previous menu, you will access the first menu for formulating your new program:

<p><PROGRAM NAME ></p> <p><DRYING></p> <p><PHASES></p> <p><Confirm Update></p>
--

The first information to enter is the name of the new program. Select the “PROGRAM NAME” option to access the mask below:

<p>PROGRAM NAME</p> <p>-----</p>

Create the name of the new program using the +, -, Rh Ar and Lh Ar keys.
Confirm with ENTER to go back to the previous mask.

<p><PROGRAM NAME ></p> <p><DRYING></p> <p><PHASES></p> <p><Confirm Update></p>
--

Select the “DRYING” item to access the menu where this phase can be entered:

<p>Drying <NO></p> <p>Temperature <000></p> <p>Extension <0000></p>

Confirm the entered data to go back to the previous menu.

Having entered the name of the program and the drying phase, you can now proceed by editing the true phases of the washing and/or disinfecting program.

7.3.1 How to select the program phases

Each washing program consists of up to 6 different phases.

The actions required can be programmed for each of the 6 phases, i.e.:

1. define the type of water used in the phase (cold, hot or demineralized)
2. define the type of detergent to add to the washing water by activating one or more of the dispenser pumps.
3. define whether to activate an auxiliary function or not, e.g. activation of an alternative drain valve
4. define whether the washing water must be heated or not, and at what temperature
5. define how long the materials are to be washed at the selected temperature

These data must be defined for each of the 6 phases.

Access the main menu for creating the new program and select the “PHASES” item:

<PROGRAM NAME >

<DRYING>

<PHASES>

<Confirm Update>

Now proceed by creating a new phase:

<NEW>

Press “ENTER” to access:

PHASE 00

Confirm?

Confirm, to access the setup menu of the individual phase.

7.3.2 How to edit each individual phase

Cold WATER <NO>

Hot WATER <NO>

Demi WATER <NO>

Quantity <09.0>

Mask for selecting the type and amount of water (expressed in liters)

FIRST TARGET

Heating <NO>

Temperature <00>

Extension <000>

“FIRST TARGET” means the temperature required for the washing/thermal disinfection phase. A temperature can be fixed, or the “Heating” option selected without a precise temperature. Heating will last for the time entered for the “Extension” item (expressed in seconds).

SECOND TARGET

Heating <NO>

Temperature <00>

Extension <000>

“SECOND TARGET” means the temperature required for the washing/thermal disinfection phase. A temperature can be fixed, or the “Heating” option selected without a precise temperature. Heating will last for the time entered for the “Extension” item (expressed in seconds).

P1 Deter <NO> <00.2>

P2 Neutr <NO> <00.2>

P3 Disinf <NO> <00.2>

P4 Lubri <NO> <00.2>

P5 Det N <NO> <00.2>

Pay close attention to the measuring unit: the value are given in ml/lt. Take into account that the machine loads 10 lt/phase.

Phase with drained water recovery

<NO>

7.4 HOW TO EDIT AN EXISTING PROGRAM

If, for example, the “PREWASH” program is selected from the program selection menu (“PROGRAMMING” sub-menu) instead of “NEW”,

<NEW>
<PREWASH>
<DEMI PREWASH >
<PLASTIC 70°C>

you will access the following mask:

<CHANGE>
<COPY>
<DELETE>

There are 3 options to choose from.

7.4.1 How to modify an existing program

<CHANGE>
<COPY>
<DELETE>

Once you have confirmed the selection made, you will access the following menu:

PROGR NAME
<DRYING>
<PHASES>
<CONFIRM UPDATE>

Apart from the name of the program, you can change the drying parameters and the various phases. However, in relation to these latter, only the quantities of detergent can be modified (expressed in percentage form).

7.4.2 How to copy an existing program

<CHANGE>
 <COPY>
 <DELETE>

Confirm, to access the following menu, where the copy of the program about to be copied can be renamed.

PROGRAM NAME

7.4.3 How to delete an existing program

<CHANGE>
 <COPY>
 <DELETE>

The programs installed by Smeg cannot be deleted.
Only new programs created by the customer can be cancelled.

Once the DELETE option has been selected, it must be confirmed on the following screen.

Confirm

 Prog Deletion

 <Y>

The selection is confirmed by selecting the response with the + and – keys and pressing ENTER

8. SUPPLEMENTARY MAINTENANCE

8.1 *PUMP BLOCKED AFTER A LONG IDLE PERIOD*

During long periods at a standstill, the motor's rotating carbon seal may block owing to lack of lubrication or simply because it has remained without water.

This will be evident when the first cycle starts as the motor will fail to function.

WARNING:

DO NOT ALLOW THE MOTOR TO FUNCTION IN THE BLOCKED CONDITION FOR LONG PERIODS AS IT MAY BURN OUT.

Attempt to release the motor in the following way:

- Remove the coarse filter
- You will see 2 holes at the bottom of the trap:
 - One is vertical in the central position
 - The other hole is horizontal, towards the bottom of the machine
- The pump vanes can be seen through the hole
- Try to turn the vanes without forcing with a finger
- If they fail to turn by the simple force exercised by a finger, use a large wide screwdriver to delicately lever (take care not to break the vanes) in one direction or the other, shifting the vanes step by step until the pump releases.
- Attempt to start the machine again.

9. TABLE OF PROGRAMS

Key to the tables

C = cold water

H = hot water

D = demineralized water


P1, P2, P3, P4, P5 = detergent dispenser pumps

EXT = target temperature upkeep time


A₀ = degree of intensity of the thermal disinfection phase: the higher this is, the more efficacious thermal disinfection will be (however, the temperature must be at least 70°C)

10. HYDRAULIC CIRCUIT


10.1 SYSTEM OVERVIEW




16.2 MAIN HYDRAULIC CIRCUIT - CONNECTIONS




10.3 DETERGENT DISPENSING SYSTEM



10.4 DRYING SYSTEM




11. WIRING DIAGRAMS



11.1 RACK BEHIND THE FRONT LOWER DOOR


TERMINAL
BOARD

POWER SUPPLY
LFM60S240



11.2 TERMINAL BOARD

center-star configuration (the picture below is for reference only, but the real situation is other)



Grey terminal board (GW3060 wire configuration):

Front side: from right to left

NR1: neutral for washing heating element R1

NR2: neutral for washing heating element R2

NR3: neutral for washing heating element R3

Mors – N: neutral from terminal board

T phase: phase for washing heating element R3

S phase: phase for washing heating element R2

Mors- R (red wire): phase for electronics and controls


R phase: phase for washing heating element R1

Back side:

Supply cable entry


(In case of GW4090 model, there are two more wires: blue and brown (neutral + phase) for the cabinet

11.3 ELECTRONICS



The rack has two electronic boards joined by connectors: motherboard and backpanel.

11.4 ELECTRONIC BOARD - layout



The list of connectors is given below with an indication as to the signals:

. P1

1. Neutral (blue)
2. Phase (brown)

. P2

1. P3 (peristaltic pump)
2. P4 (peristaltic pump)
6. P1(peristaltic pump)/Detergent Dispenser
7. P2 (peristaltic pump) / Acid batcher
8. EVU output solenoid valve: not installed in this model
9. EVSA alternative draining solenoid valve: optional
10. MA1 drying motor 1
11. PSC condenser draining pump
12. EVCC condenser filling solenoid valve

. P3

1. EVF cold water filling solenoid valve
2. MS draining pump motor
3. ..
4. ..
5. EVR regeneration solenoid valve
6. BPE1 door lock 1
7. EVD demineralized water filling solenoid valve
8. EVC hot water filling solenoid valve

. P4

1. SS salt sensor
2. common contact
3. ..
4. ..
5. ..
6. ..
7. ..
8. PAP1 high pressure switch on was pump
9. MCA acquastop microswitch (not installed in this model)
10. MCP1 door lock electric opening microswitch
11. MCM1 door lock mechanical opening microswitch

. P5

1. PL washing level pressure switch
2. common contact
3. PS safety level pressure switch

. P6

(not used in this model)

. P7

1. TA1 drying temperature probe
- 2.
3. TL1 temperature probe of the washing chamber
- 4.
5. TB temperature probe in external boiler (not installed)
6. TL2 additional temperature probe of the washing chamber (optional)

. P8 (not used in this model)

1. IN
2. OUT
3. +24
4. GND
5. GND
6. ..

. P9

1. SP1 can 1 level probe
2. SP4 can 4 level probe
3. GND ground reference signal
4. SP3 can 3 level probe
5. SP5 can 5 level probe
6. SLC operating level probe of condenser
7. SP2 can 2 level probe

. P10

..

. P11

- 1.. ..
- 2.. ..
3. FMD demineralized water flowmeter
4. FMF cold/hot water flowmeter
- 5.. ..
6. GND ground reference signal

. P12 (not used in this model)

. P13 (not used in this model)

. P14

1. --
2. +24V
3. GND

List of relays on the backpanel

K1 = maximum speed MA1
K2 = RA1
K3 = R3
K4 = booster pump for demi water
K5 = R1
K6 = R2

The list of fuses on the backpanel board is given below:

- . BF1 = not present
- . BF2 = **T2A, 250V, 5x20 mm** (power supply +24V)
- . BF3 = **T2A, 250V, 5x20 mm** (PPSC)
- . BF4 = **T2A, 250V, 5x20 mm** (MA1)
- . BF5 = **T2A, 250V, 5x20 mm** (MS)
- . BF6 = **T2A, 250V, 5x20 mm** (P1/Detergent dispenser)

12. TROUBLESHOOTING – Failure search and solutions

12.1 PROBLEMS or FAILURES NOT ASSOCIATED WITH A FAULT

The GW4090 is a complex system equipped with many sensors used to control the machine and make searching for failures easier. One or more faults correspond to each sensor. The appearance of these faults help the search for the cause of the failure.

In any case, “malfunctions” may appear that are not associated with any fault or that present the user with problems that cannot be solved immediately.

In this paragraph, we will try to identify the possible causes at the heart of the most common malfunctions.

12.1.1 The machine does not turn on

Check, in the following order:

- . connection to the electrical cabinet
- . fuses in the main terminal board (see §11.2) or clamps that have not been locked
- . green light on the power supply in the lower electrical box
- . defective console card

12.1.2 Water leaks

If water is leaking from the front of the machine, it probably is temporary. In any case, make sure that:

- . the water load is correct
- . there is no water in the tank at the end of the drain pipe. If this is not the case, clean the water trap filters and check the drain pipe position
- . the pressure switch trips (in any case, the relative alarm should intervene), using WDTRACE®
- . the problem could be due to water leakage from the water trap or from the various holes present on the bottom of the tub. In this case, incline the machine and check where traces of limestone are found or humid parts. Then, replace the component/gasket.

If water is leaking from the back of the machine, remove the rear guard and make sure that:

- . the condenser is not leaking
- . the condenser is draining correctly

12.1.3 Problems regarding drying

The most important parameters that influence drying are temperature and duration. The higher these values are, the more effective is drying. In any case, we recommend not exceeding 110°C.

If the material is still wet when removed from the machine after the drying phase has completed, do the following:

- . make sure that the material is not placed in a manner that blocks water downflow onto the surfaces due to gravitational forces
- . try to increase the time and temperature of the drying phase (without exceeding 110°C)
- . leave the material in the washing chamber for 10 minutes, if possible, with the door open so that the external air helps the evaporation process of the residual water drops present on the material

12.1.4 Problems with the quality of washing results

If a customer has a problem with bad washing results, first make sure that the following factors were taken into consideration:

- . the material should be arranged in the washing trolley in the best way possible for optimal washing results
- . the quality and quantity of the water used
- . the water pressure in the washing ducts
- . the temperature of the machine during washing
- . detergents used and correct dosage

12.1.4.1 Arrangement of the instruments in the washing trolley

Take the following precautions:

- . do not fill the trolleys with too much material at the same time
- . make sure that the instruments are positioned so that they do not overlap each other
- . if the instruments have articulated parts, open the articulations as far as they will go
- . do not place instruments covered with blood or other substances that have been allowed to dry, into the trolleys
- . bulky instruments must be put in the baskets in a manner that prevents "shadow areas" from being created, as these areas obstruct cleaning of other instruments
- . instruments with cavities (flexible tubes, respiration systems, etc.) must have the possibility of being completely rinsed, even internally. Therefore, inserts with washing devices designed specifically for treating this type of instrumentation must be used.
- . position the instruments so that they are unable to obstruct the spraying arms

12.1.4.2 Quality and quantity of water

It is essential that good quality water is used. Have the water in the water supply system periodically checked. Connect the demineralized water supply when possible (conductivity < 8-10 µS).

It may be advisable to increase the amount of water used for each phase (by means of the relative menu) according to the selected washing trolley and program. This parameter may not be changed as regards Smeg's programs. Be careful when using this parameter and do not exceed 12-13 liters.

12.1.4.3 Water pressure

The correct water pressure in the machine's hydraulic circuit is assured by correct operation of the washing pump and the entire hydraulic circuit. If the pressure is not correct, a fault is activated on the display. However, the sprayers must be checked to see if its holes are blocked/clogged. Also check that the filters present in the water trap are clean and that the pump suction device in the water trap is not clogged.

12.1.4.4 Washing temperature

The optimum temperature at which washing takes place depends on the type of detergent used and the instruments to process. However, washing temperatures must usually exceed 45°C, otherwise the chemicals will not be effective. However, the temperature must never exceed 65°C for detergents recommended by Smeg.

12.1.4.5 Detergents used

The detergents used and the time they are in action represent the most important factors for obtaining an optimum washing result. Smeg guarantees that optimum washing results are obtained so long as the customer uses the products tested along with the appliance. This element is essential to reach the desired washing results.

A list of the products for which Smeg guarantees the correct washing results follows:
for programs 3 to 11, which are the ones normally used, Smeg recommends:

alkaline detergents: Smeg's DETERLIQUID C and C2
Smeg's DETERLIQUID D and D2

slightly alkaline detergents Smeg's EMODET (for programs from 15 to 20)

acid neutralizers: Smeg's ACIDGLASS C and C2 (recommended)
Smeg's ACIDGLASS P and P2 (very strong, use only if absolutely necessary)

disinfectant: Smeg's STREPTOBAT (for chemical disinfection of thermolabile instruments)

lubricant: Smeg's LUBMILK (introduced in the last phase. It is used primarily for making instrument drying easier and protecting surface conditions)

12.1.4.6 Spots on the instruments

The possible causes may be:

- . instruments arranged incorrectly
- . bad quality mains water used in the instrument washer
- . demineralized water or good quality softened water missing from the final washing phase. This usually represents the decisive factor: we recommend always using demineralized water in the final washing phase
- . excessive detergent dosage and washing temperature too high
- . not enough neutralizer added
- . salt residue in the tank due to incorrect sealing or a deformity of the decalcification cap

12.2 FAILURES ASSOCIATED WITH A FAULT

There are two different situations, depending on how serious the fault is:

- . warning messages;
- . true alarm signals.

In the first case, a message will appear on the display allowing the user to ignore the warning and proceed with the selected cycle. In the second case, the user must press the Reset button and comply with the relative procedure in order to resolve the situation.

By and large, an alarm indicates that the machine is operating in a faulty way and needs to be repaired by a technician. Sometimes, however, the alarm may be caused by a temporary situation.

The possible alarm messages are indicated below.



To release the machine when one of these messages appears, press the RESET key for a couple of seconds in order to activate the RESET cycle.

If a cycle has been stopped during a high temperature phase, remember that the temperature must drop below 45°C before the door can be opened.

Alarm ID	Message	Action
1	WATER HEATING FAILED	<p>This alarm is only active during the heating phase.</p> <p>If a <u>three-phase machine</u> is being used, the alarm is issued if the temperature does not increase by 3°C in 2 min.</p> <p>If a <u>single-phase machine</u> is being used, the alarm is issued if the temperature does not increase by 2°C in 1 min.</p> <p>Check the condition of the safety thermostat (first plunger on the upper wall of the detergent reservoir). Use WDTRACE® to make sure that the resistor drive relay activates. Check the state of the fuses R1, R2, R3 in the main terminal board.</p>
2	TEMP. PROBE TL1-TC	<p>This alarm is always active during cycle execution, except during the sub-phases, namely: loading, unloading or mixed unloading, reset and drying.</p> <p>It intervenes when $TL1 > TC + 2^{\circ}\text{C}$</p> <p>Repeat the cycle: It could be a temporary alarm. Check the state of the two probes in the water trap. Clean their surfaces, if necessary. If the problem persists, replace one of the two on the basis of a comparison made with an external thermometer.</p>
4	OVERTEMPERATURE TL1	<p>This alarm is always active during cycle execution, except during the sub-phases, namely: loading, unloading, reset and drying.</p> <p>It intervenes when $TL1 > T_{MAX} + 5^{\circ}\text{C}$, where T_{MAX} is the max. target temperature value of the phase being executed.</p> <p>Repeat the cycle: It could be a temporary alarm. Otherwise, act as indicated in the previous case (point 2).</p>
5	PROBE TL1 DISCONNECTED	<p>Probe TL1 is damaged or disconnected.</p> <p>Check the connection and replace the probe, if necessary.</p>
7	PROBE TA1 DISCONNECTED	<p>Probe TA1 is damaged or disconnected.</p> <p>Check the connection and replace the probe, if necessary.</p>
9	PROBE TB DISCONNECTED	<p>Probe TB (optional boiler) is damaged or disconnected.</p> <p>Check the connection and replace the probe, if necessary.</p>
10	PROBE TCL DISCONNECTED	<p>Probe TCL is damaged or disconnected.</p> <p>Check the connection and replace the probe, if necessary.</p>
11	LACK OF COLD WATER	<p>This alarm is active during the loading phase.</p> <p>It intervenes when the relative turbine does not count any impulses for 30 min. and the state of the pressure switch does not change.</p> <p>Check the water supply:</p> <ul style="list-style-type: none"> . cock and connection pipe open/closed . solenoid valve input . mains water pressure

		<ul style="list-style-type: none"> . the state of the pressure switch . turbine operation
12	LACK OF WARM WATER	<p>This alarm is active during the loading phase. It intervenes when the relative turbine does not count any impulses for 30 min. and the state of the pressure switch does not change.</p> <p>Check the water supply:</p> <ul style="list-style-type: none"> . cock and connection pipe open/closed . solenoid valve input . mains water pressure . the state of the pressure switch . turbine operation
13	LACK OF DEMI WATER	<p>This alarm is active during the loading phase. It intervenes when the relative turbine does not count any impulses for 30 min. and the state of the pressure switch does not change.</p> <p>Check the demineralized water supply:</p> <ul style="list-style-type: none"> . cock open/closed . tank/can empty . solenoid valve input . when the tank is loaded, check water input into the tub; . connection pipe . connection/program correspondence . turbine
14	COLD WATER PRESSURE	<p>Check:</p> <ul style="list-style-type: none"> . water input pressure . the input pipe should not be restricted in any way
15	WARM WATER PRESSURE	<p>Check:</p> <ul style="list-style-type: none"> . water input pressure . the input pipe should not be restricted in any way
16	DEMI WATER PRESSURE	<p>Check:</p> <ul style="list-style-type: none"> . water input pressure . the input pipe should not be restricted in any way
17	COLD WATER LOAD TIME EXCEEDED	<p>The machine took too much time to load the cold water. Check the input pressure. Increase the loading timeout (see Extra Loading Time menu). It has a maximum of 180 minutes.</p>
18	WARM WATER LOAD TIME EXCEEDED	<p>The machine took too much time to load the hot water. Check the input pressure. Increase the loading timeout (see Extra Loading Time menu). It has a maximum of 180 minutes</p>
19	DEMI WATER LOAD TIME EXCEEDED	<p>The machine took too much time to load the demineralized water. Check the input pressure. Increase the loading timeout (see Extra Loading Time menu). It has a maximum of 180 minutes</p>
20	WATER LOAD SYSTEM FAILURE	<p>This alarm is active in the loading phase with the exception of static loads and water trap washing. It intervenes if: {water load ok & water quantity. >7.5 liters} & {PAP1 = off}, where PAP1 is the high pressure switch. Therefore, check PAP1. PAP1 status is identified only after the load has been completed.</p> <p><u>Note:</u> Do not set static loads with higher water quantities than 7.5 liters.</p>
21	FLOWMETER CONNECT.	<p>The flow meter connections relative to the</p>

	INVERTED	<i>detergents are not correct. Check the correspondence between the number of FM's and the associated peristaltic pump.</i>
22	COLD WATER FLOWMETER FAILURE	<i>Check the water supply (cock open/closed, water pressure, connection pipe, etc.) Repeat the cycle because there may be a temporary malfunction. Check the cold water turbine.</i>
23	INSUFFICIENT WATER IN CHAMBER	<i>There is not enough water in the tub. It intervenes if: $\{PL(ON) \rightarrow PL(OFF)\} \text{ and } (PAP1 = 1)$. That is, the PL level pressure switch closes and opens while PAP1 is closed at the same time (the high pressure switch. The one regarding the hydraulic unit pressure). Check the water supply (cock open/closed, water pressure, connection pipe, etc.) and pressure switch operation.</i>
24	NO WATER IN CHAMBER	<i>The alarm is active from the end of the loading phase to unloading. It intervenes when: $\{PL(ON) \rightarrow PL(OFF)\} \text{ and } (PAP1 = 0)$. That is, there is enough water to trip the level pressure switch, but not enough to trip the high pressure switch. Check the water supply (cock open/closed, water pressure, connection pipe, etc.).</i>
25	NO PRESS IN HYD SYS: FOAM	<i>Check the type of detergent used. Above all, check Lubmilk and Emodet, as they are products that may generate a lot of foam. Be careful with the quantity used. Contact Smeg for any doubts. Make sure the instruments are not placed in the machine before eliminating detergent residue used in previous manual decontamination treatments. Repeat the cycle after cleaning the foam from the tub. It intervenes when: $\{PL(ON) = 0\} \text{ and } (PAP1 = 0)$ for more than 3 min. There is water in the tub, but no pressure in the hydraulic circuit.</i>
26	EVF LEAKAGE (EVF = cold water valve)	<i>Check the water supply (cock open/closed, water pressure, connection pipe, etc.) and solenoid valve operation. Make sure the turbine is operating correctly (by exchanging the connections between the supply pipes).</i>
27	EVC LEAKAGE (EVC = hot water valve)	<i>Check the water supply (cock open/closed, water pressure, connection pipe, etc.) and solenoid valve operation. Make sure the turbine is operating correctly (by exchanging the connections between the supply pipes).</i>

28	EVD LEAKAGE (EVD= demi water valve)	Check the water supply (cock open/closed, water pressure, connection pipe, etc). Check solenoid valve operation. Make sure the turbine is operating correctly (by exchanging the connections between the supply pipes).
29	WASHING CHAMBER DRAINAGE FAILURE	It intervenes when: (PAP1oPL) = 1 after 3' from the start of the unloading phase. Make sure the drain is positioned as shown in the manual. Make sure the drain pump operates and that the pressure switch activates. Check the positions of the instruments: do not insert objects such as basins with the concave side up.
30	SAFETY CHAMBER LEVEL EXCEEDED	Active during the washing cycle. It intervenes if (PS = 0). That is, the safety level has tripped for at least 5 min. Make sure that the upstream solenoid valve is not leaking. Check loading pressure switch operation (work + safety level).
31	SAFETY LEVEL FAILED	This alarm only appears after alarm 32 appears if PS does not reset.
32	SUMP FULL	Active when the loading pressure switch trips, even though there is no cycle being executed. Check the state of the solenoid valve input and make sure that the machine is able to drain.
33	LACK OF WATER IN THE STEAM CONDENSER	Indicates that there is no water in the vapor condenser when there should be due to sprayer nozzle activation. Make sure the solenoid valve input to the condenser activates and operates correctly. Also make sure that the drain pump does not remain activated permanently.
34	CONDENSER DRAINAGE FAILED	Intervenes if SLC (work level sensor of the condenser) does not deactivate after 120 min. from drain pump activation. Make sure the drain pump operates correctly. Also make sure that the condenser drain respects the recommended values and it not obstructed. Replace the level sensor, if necessary.
35	CONDENSER LEVEL SWITCH FAILED	Not active
36	CONDENSER DRAINING PUMP FAILED	Not active
37	NOT REACHED TARGET DRAINAGE AT T=	Alarm active during the mixed draining phase if the target temperature is not reached. The alarm appears at the end of the cycle. This may be a temporary malfunction. However, if this is not the case, make sure cold water is poured into the tub during draining.
38	COOLED DRAINAGE FAILED	Intervenes if TL1 does not start dropping when the mixed draining phase begins. See previous case.
41	DETERGENT 1 INFLOW FAILED	Intervenes during the phase when the detergent is put into the tub. The flow meter associated with

	(flowmeter nr.1 optional)	<i>P1 does not count any impulses for 60 min. Make sure that the suction lift pipes of can P1 are in good condition and not crushed. Make sure P1 activates correctly. Check the delivery pipe path to the tub inlet.</i>
42	DETERGENT 2 INFLOW FAILED (flowmeter nr.2 optional)	<i>See previous case associated to P2.</i>
43	DETERGENT 3 INFLOW FAILED (flowmeter nr.3 optional)	<i>See previous case associated to P3.</i>
44	DETERGENT 4 INFLOW FAILED (flowmeter nr.4 optional)	<i>See previous case associated to P4.</i>
46	PUMP 1 TUBE CLOGGED	<i>Alarm active during P1 activation. It intervenes if FM1 (flow meter associated to P1) does not reach the ml target number within 90 min.</i>
47	PUMP 2 TUBE CLOGGED	<i>See previous case associated to P2.</i>
48	PUMP 3 TUBE CLOGGED	<i>See previous case associated to P3.</i>
49	PUMP 4 TUBE CLOGGED	<i>See previous case associated to P4.</i>
51	VERIFY CONNECTIONS FLOWMETER/PUMPS	<i>Alarm disabled.</i>
52	DOOR 1 ELECTRICALLY OPEN	<i>Alarm active during the entire washing cycle (from the beginning to the end of the cycle). It intervenes if MCP=0. That is, if the microswitch for electric door closing is in an incorrect state. Check MCP operation and make sure the door closes properly.</i>
54	DOOR 1 MECHANICALLY OPEN	<i>Alarm active during the entire washing cycle (from the beginning to the end of the cycle). It intervenes if MCM=0. That is, if the microswitch for mechanical door closing is in an incorrect state. Check MCM operation and make sure the door closes properly. Above all, make sure that door opening/closing is not obstructed at the bottom of the door.</i>
56	DOORLOCK 1 FAILURE	<i>Make sure that door opening/closing is not obstructed at the bottom of the door. Make sure the solenoid valve stroke is not obstructed in any way.</i>
58	DRYING 1 FAILURE	<i>Alarm active during heating of the drying phase. Make sure the temperature measured by the drying probe increases as time passes. Check connection to the probe. Make sure the drying cycle safety thermostat has not tripped.</i>
60	DRYING 1 SAFETY SWITCH ON	<i>It activates if the temperature of the drying resistor exceeds the safety limit value. There may be a serious problem. Therefore, shut the machine off immediately and check the state of the drying system. Make sure, in particular, that the filter at the can inlet is not clogged and that the valve in</i>


		<i>the tub air inlet duct is not blocked. Press the 2nd button on the bottom of the drawer (see §11.1) to restart the machine.</i>
62	HEATING THERMOSTAT ON	<i>It activates if the temperature of the drying resistor exceeds the safety limit value. There may be a serious problem. Therefore, shut the machine off immediately and check the state of the heating elements in the tub. It is not a problem if the surface darkens slightly. Press the 1st button on the bottom of the drawer (see §11.1) and repeat the cycle. If the problem persists, call the Technical Assistance Service.</i>
63	BLOWER 1 FAILED	<i>Not active</i>
65	FAN SENSOR 1 FAILED	<i>Not active</i>
67	COOLING FAILED - WARNING: HIGH TEMP	<i>The drying cycle ends with a cooling phase. If the final temperature is not lower than 80°C, the alarm intervenes. Execute a drying phase. Make sure the motor operates and check the temperature variations on the display. Check the temperature probe.</i>
68	TANK 1 EMPTY	<i>Fill the relative can.</i>
69	TANK 2 EMPTY	<i>Fill the relative can.</i>
70	TANK 3 EMPTY	<i>Fill the relative can.</i>
71	TANK 4 EMPTY	<i>Fill the relative can.</i>
73	ARCHIVE ERROR	<i>Error regarding archive management of the processes saved in the microprocessor. It may be a temporary problem. If the problem remains, re-download the firmware on the master processor (motherboard). Replace the motherboard if necessary.</i>
74	CHAMBER LEAKAGE	<i>Not active</i>
75	LACK OF SALT	<i>Fill the softener's salt reservoir (plug in tub).</i>
76	BOILER SAFETY SWITCH ON	<i>Not active (present on another machine model).</i>
77	TEMPERATURE > 45°C	<i>Intervenes (once the relative option is selected) if the water inlet temperature is higher than 45°C during the first phase. If a cycle has just ended, leave the tub door open for 10 minutes to cool it (especially if a drying cycle was not performed).</i>
78	RECOVERY FAILED	<i>Internal microprocessor error. This may occur when firmware download is completed. This is a temporary alarm. Therefore, perform the operation again. If the problem persists, contact headquarters.</i>
79	PROGRAM NOT CONGRUENT	<i>Non-congruent program error. Make sure the settings regarding the desired cycle to be executed are not in contrast with setting the washing parameters.</i>
80	WARM WATER FLOWMETER	<i>Check the water supply (cock open/closed, water</i>

	FAILURE	<i>pressure, connection pipe, etc).</i> <i>Repeat the cycle as the malfunction may be temporary.</i> <i>Check the hot water turbine.</i>
81	DEMI WATER FLOWMETER FAILURE	<i>Check the water supply (cock open/closed, water pressure, connection pipe, etc).</i> <i>Repeat the cycle as the malfunction may be temporary.</i> <i>Check the demineralized water turbine.</i>
82	SYSTEM FAILURE SOL.1	<i>Alarm active during the open door phase.</i> <i>It intervenes if the door remains closed after a series of opening impulses have been provided.</i> <i>Make sure the solenoid valve stroke is not blocked and that the opening button on board the tub trips correctly.</i>
84	OVERTEMPERATURE TA1	<i>It intervenes if the drying probe measures a temperature above 140°C for over 5 min.</i> <i>Make sure the probe is inserted correctly in the air duct and that nothing is blocking passage through the tub inlet.</i> <i>The relative safety thermostat may trip after this alarm appears.</i>
86	PUMP P1 BLOCKED	<i>Peristaltic pump blocked. Check and replace it if necessary.</i>
87	PUMP P2 BLOCKED	<i>Peristaltic pump blocked. Check and replace it if necessary.</i>
88	PUMP P3 BLOCKED	<i>Peristaltic pump blocked. Check and replace it if necessary.</i>
89	PUMP P4 BLOCKED	<i>Peristaltic pump blocked. Check and replace it if necessary.</i>
91	ARCHIVE FULL	<i>Indicates the archive of saved cycles is full.</i> <i>Delete the archive once the data is unloaded.</i>
92	CHANGE FILTER	<i>Replace the air filter.</i>
93	MAINTENANCE REQUEST	<i>Periodic maintenance.</i>

13. INSTALLATION LAYOUTS

13.1 MODEL GW4090

KEY	
E- YELLOW/G REEN	GROUND (YELLOW GREEN)
N- BLUE	NEUTRAL
- BROWN -	PHASE 1- BROWN
- BLACK -	PHASE 2 - BLACK
- BLACK -	PHASE 3 - BLACK
CW	COLD WATER CONNECTION



ELECTRICAL AND HYDRAULIC CONNECTIONS REQUIRED	
TYPE OF VOLTAGE PERMITTED ± 10 %	400V 3 PHASES WITH NEUTRAL - 50 Hz
POWER RATING	7.0 kW
ELECTRICAL CONNECTION	MAGNETOTHERMIC CIRCUIT-BREAKER 3P+N In=20A
WET CONNECTIONS	CW cold HW hot
TYPE OF WATER	standard threaded 3/4" gas
TYPE OF CONNECTION	DRINKING WATER
QUALITY OF THE WATER	DRINKING WATER
PRESSURE IN BAR MIN - MAX	3-10
MAX° FHARDNESS total max (CONDUCTIVITY)	42°F
iron ppm Fe max	< 0.5 ppm
WATER DRAIN PIPE	FLOOR LEVEL
height	minimum height 650
diameter	minimum 30 mm
	ON WALL
	minimum height 650 - maximum height 1,000 mm
	hose nipple 21 mm
	UNDER SINK
	minimum height 650 mm
	raccordo T con portagomma 21 mm

