




# Camlab Water Model CW/6120

Conductivity Meter - Instruction Manual  
Ref: CW/72.23.00



## **CONTENTS**

1. GENERAL DESCRIPTION
2. SPECIFICATION
3. FUNCTIONAL DESCRIPTION
4. TAKING MEASUREMENTS
5. CALIBRATION
6. CHANGING THE BATTERY
7. ACCESSORIES

### **1. GENERAL DESCRIPTION**

\* General purpose conductivity meter with broad application including AQUARIA & FISH HATCHERIES, FOOD & BEVERAGE PROCESSING, PHOTOGRAPHY, LABORATORY, PAPER INDUSTRY, PLATING INDUSTRY, QUALITY CONTROL, EDUCATION, SWIMMING POOLS & WATER CONDITIONING

\* High quality, compact unit with a separate electrode that is designed for easy operation

\* Water resistant front panel with easy to read LCD display and rubberised function keys

## Unpacking



Please check that the shipment includes the following items:

- Camlab Water CW6120 Meter
- Conductivity Sensor
- 9V Battery
- Protective Cover
- 7 Screws
- Screwdriver
- Plastic Cover Calibration Screws
- Instruction Manual

## 2. SPECIFICATION

|                             |   |
|-----------------------------|---|
| Display                     | LCD, 21.5 mm (0.7") digit height<br>Maximum display count no. 1999                  |
| Measurement Range           | 2 ranges: 0 - 1.999 mS/cm,<br>0 - 19.99 mS/cm                                       |
| Resolution                  | 0.001 mS/cm for 0 - 1.999 mS/cm<br>range<br>0.01 mS/cm for 0 - 19.99 mS/cm<br>range |
| Accuracy (23 ± 5°C)         | 3% full scale + 1 digit   |
| Sample Time                 | Approx. 0.4 seconds   |
| Over Range Indicator        | Display shows „1“   |
| Data Hold                   | Freezes the conductivity value on the<br>display                                    |
| Temperature<br>Compensation | Automatic, 0 - 50°C (32 - 122°F)  |
| Operating Temperature       | 0 - 100°C (32 - 212°F)  |
| Operating Humidity          | Maximum 80% relative humidity   |

|               |  |
|---------------|--|
| Power Supply  | 006P DC 9V battery (heavy-duty type),<br>MN1604 (PP3) or equivalent                  |
| Power Current | Approx. DC 5 mA  |
| Weight        | 380 g  |
| Dimensions    | Meter: 208 x 110 x 34 mm (L x W x H)<br>Electrode: 22 mm diameter x 120 mm<br>length |



### 3. FUNCTIONAL DESCRIPTION




Figure 1

- 1 Display
- 2 Power ON Button
- 3 Power OFF Button
- 4 Data Hold Button
- 5 2 mS/cm Range Button
- 6 20 mS/cm Range Button
- 7 Electrode Input Socket
- 8 Calibration Adjust (HR, LR)
- 9 Battery compartment/cover
- 10 Stand
- 11 Electrode Handle
- 12 Conductivity Electrode
- 13 Protective Cover

## 4. TAKING MEASUREMENTS

### Battery installation

Prior to first use take the instrument out of the protective cover open the battery compartment and insert the 9V battery.

Ensure polarity is correct.

### Protective Cover

The instrument is equipped as standard with the protective cover.

Prior to the first measurement please fix the cover by using the enclosed screws. The protective cover ensures reliable operation even in harsh environments.


When the instrument is used without the protective cover please use the enclosed plastic cover to protect the calibration screws against dust and humidity.

- 1) Connect the Conductivity Electrode (Figure 1, 12) to the Electrode Input Socket (Figure 1, 7).
- 2) Power on the instrument by pressing the Power ON Button.
- 3) Select the 2 mS/cm or 20 mS/cm range by pressing the appropriate range button (Figure 1, 5 or 6).
- 4) Holding the electrode handle (Figure 1, 11), immerse the conductivity electrode completely in the sample.

Shake the electrode several times to remove air bubbles from the electrode and therefore ensure stable readings. The instrument will display the conductivity value in mS/cm.

Note: „1“ in the display indicates an out-of-range measurement. When operating in the 20 mS/cm range, if the value obtained has one or more zeros after the decimal point, change to the 2 mS/cm range for improved accuracy.


## 5. CALIBRATION




To calibrate the meter:

- i) Prepare a 1.413 mS/cm calibration solution (Ref RN/CSKCL) or similar.
- ii) Select the 2 mS/cm range (Figure 1, 5).
- iii) Holding the conductivity electrode by its handle (Figure 1, 11), immerse it completely in the calibration solution. Shake the electrode to disperse air bubbles on it and allow measurements to stabilise. Adjust „Calibration Adjust LR“ (Figure 1, 8) until the display reads exactly 1.413 mS/cm.

## 6. CHANGING THE BATTERY



- i) „“ in the left corner of the display indicates that it is necessary to replace the battery. However, accurate measurements may still be made for several hours after the “Low Battery” indicator first appears.
- ii) Remove the protective cover before replacing the battery. To replace the battery, remove the Battery Compartment Cover (Figure 1, 9) on the rear of the meter.
- iii) Remove the battery, install a replacement one (006P DC 9V battery (heavy duty type), MN1604 (PP3) or equivalent) and replace the battery compartment cover.

## 7. Replacement Parts

Calibration Solution 1413  $\mu\text{s}/\text{cm}$  500 ml




Ref: RN/CSKCL

Conductivity Measuring Cell




Ref: TT/72.22.25






**Notes**



**Notes**



**Notes**



Choice, Service and Trust

**camlab**

Camlab Limited , Camlab House, Norman Way, Over  
Cambridge CB4 5WE

Tel: 01954 233110 • Fax: 01954 233101

[www.camlab.co.uk](http://www.camlab.co.uk)

Technical changes may occur without notice