



SODIUM HYDROXIDE KIT

DIRECT READING TITRATOR, 0-10%

CODE 7516-DR-02

QUANTITY	CONTENTS	CODE
30 g	*Barium Chloride Powder	*6073-G
2 x 30 mL	*Hydrochloric Acid, 2.5N	*6251DR-G
60 mL	Deionized Water	5115PT-H
15 mL	*Phenolphthalein Indicator, 1%	*2246-E
1	Pipet, 1.0 mL, plastic	0354
1	Test Tube, 5-10-15-20-25 mL, plastic, w/cap	0715-DRT
1	Direct Reading Titrator, 0 - 10 Range	0377
1	Spoon, 0.5 g, plastic	0698

*WARNING: Reagents marked with an * are considered to be potential health hazards. To view or print a Material Safety Data Sheet (MSDS) for these reagents go to www.lamotte.com. To obtain a printed copy, contact LaMotte by e-mail, phone or fax.

To order individual reagents or test kit components, use the specified code number.

Read the LaMotte Direct Reading Titrator Manual before proceeding. The Titrator is calibrated in % Sodium Hydroxide. Each minor division equals 0.2%.

PROCEDURE

1. Use the 1.0 mL pipet (0354) to add 1.0 mL of sample to test tube (0715-DRT).
2. Fill the test tube to 10 mL line with Deionized Water (5115).
3. Use the 0.5 g spoon (0698) to add one level measure of *Barium Chloride Powder (6073). A white precipitate will form if carbonates are present.
4. Add two drops of *Phenolphthalein Indicator, 1% (2246). Cap and mix. Solution will turn pink.
5. Fill the Direct Reading Titrator (0377) with the *Hydrochloric Acid, 2.5N (6251). Insert the tip of the Titrator into the center hole of the test tube cap.
6. While gently swirling the tube, slowly press the plunger to titrate until the pink color disappears.
7. Read the concentration of the test sample where the plunger meets Titrator scale. Record as % Sodium Hydroxide.

EXAMPLE: Titrator reading is 3 minor divisions below line 7.

$$7 + (3 \text{ divisions} \times 0.2) = 7.6\% \text{ NaOH}$$

LaMOTTE COMPANY

Helping People Solve Analytical Challenges®

PO Box 329 • Chestertown • Maryland • 21620 • USA
800-344-3100 • 410-778-3100 (Outside U.S.A.) • Fax 410-778-6394
Visit us on the web at www.lamotte.com