**T1**

**800 mg/l:** HCOO<sup>-</sup>

**500 mg/l:** CH<sub>3</sub>COO<sup>-</sup>

**60 mg/l:** SO<sub>3</sub><sup>2-</sup>

**20 mg/l:** S<sup>2-</sup>

**6 mg/l:** NO<sub>2</sub>-N

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**LCK 388**

**Carbonate total / carbon dioxide (CO<sub>2</sub>)**

**Please check the "Edition Date" (see data table).  
Safety advice and expiry date on package.**

**Principle**

The reaction converts all carbonates and dissolved carbon dioxide into gaseous carbon dioxide, which passes through a membrane into the indicator cuvette. The change of colour of the indicator is photometrically evaluated.

**Range of Application**

Waste water, surface water, carbonated beverages

**Interferences**

The ions listed in **T1** have been individually checked up to the given concentrations and do not cause interference. We have not determined cumulative effects and the influence of other ions.

Higher concentrations of these ions cause high-bias results. Be careful to use only double-distilled water, which contains no carbon dioxide, to dilute the sample.

The measurement results must be subjected to plausibility checks (dilute and/or spike the sample).

**pH** sample ..... 4 – 10  
**Temperature** sample/reagents ..... 15 – 25°C

**Special note**

- Contamination by ambient air: Never** leave cuvette open, because carbon dioxide in the ambient air can cause high-bias results. Cuvettes must only be opened when necessary (e.g. to add sample) and must be closed or further processed **immediately** afterwards.
- Thermostat:** Heat the thermostat to **100°C** (check the temperature – higher temperatures lead to dangerous excess pressure). When this temperature has been reached, insert the cuvette combinations and start the reaction time (**60 min**) again. Insert cuvette combinations only in the small shafts in the thermostat. Do not insert them in the large shaft with reducing sleeves.
- Disposal:** Do **not** screw the cuvette combinations apart when the analysis has been completed, but press them back into the blister pack.

**NB:**

**Be sure to set the required temperature to 100°C (at 148°C the cuvette combinations may break apart).**

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
**Attention:**


Carry out each CO<sub>2</sub>-determination in sequence in order to avoid contamination by ambient air.

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
Close cuvette with original cap and invert a few times.

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
Close indicator cuvette **very tightly** with membrane double-cap. (**NB:** The barcode label must point towards the indicator cuvette.)

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
Pipette **1.0 ml** sample into the digestion cuvette.

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
**Immediately** close the digestion cuvette **tightly** with the prepared indicator cuvette. Hold cuvette combination vertically. Do **not** invert.

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
Heat in the preheated thermostat at **100°C** for **60 min**. Allow to cool to room temperature.

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**Tighten** the cuvette combination again **before inverting** it.

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Invert cuvette combination. Thoroughly clean the outside of the indicator cuvette and evaluate.

	<b>Indikatorküvette</b> <b>Cuve indicatrice</b> <b>Indicator cuvette</b>
 Barcode <sup>1)</sup>	✓

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