**T1**

**1000 mg/l:** Cl<sup>-</sup>, SO<sub>4</sub><sup>2-</sup>

**500 mg/l:** NO<sub>3</sub><sup>-</sup>, NO<sub>2</sub><sup>-</sup>, Na<sup>+</sup>, Ca<sup>2+</sup>, K<sup>+</sup>, Mg<sup>2+</sup>, NH<sub>4</sub><sup>+</sup>,  
CO<sub>3</sub><sup>2-</sup>, PO<sub>4</sub><sup>3-</sup>

**50 mg/l:** Zn<sup>2+</sup>, Cd<sup>2+</sup>, Ni<sup>2+</sup>, Pb<sup>2+</sup>

**5 mg/l:** Hg<sup>2+</sup>, Sn<sup>2+</sup>, Cr<sup>3+</sup>, Cr<sup>6+</sup>

**3 mg/l:** Fe<sup>3+</sup>, Fe<sup>2+</sup>

**Cu = Koper / Copper**

**Cu 902 = Koper na gebruik van de  
Crack-Set LCW 902**

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**GB**

**LCK 529 Copper**

**Please check the "Edition Date"  
(see data table).**

**Safety advice and expiry date on  
package.**

**Principle**

Copper(I) ions form an orange-coloured complex with the disodium salt of bathocuproine disulphonic acid. Any copper(II) ions present in the water sample are reduced to copper(I) ions by ascorbic acid before the complex is formed.

**Range of Application**

Drinking water, waste water, process analysis

**Interferences**

The ions listed in **T1** have been individually checked up to the given concentrations and do not cause interference. We have not determined cumulative effects and the influence of other ions.

Higher amounts of iron and chromium cause high-bias results.

Undissolved copper and copper contained in complexes can only be determined after digestion with Crack-Set LCW 902.

The measurement results must be subjected to plausibility checks (dilute and/or spike the sample).

**For Special attention**

For exact evaluation it is very important that there are no air bubbles in the beam path (lower half of the cuvette). To fill the cuvette in a bubble-free manner hold the pipette against the narrower face of the slightly tilted cuvette and allow the liquid to run slowly into the cuvette down the side wall. If any air bubbles should adhere to the cuvette walls they can be removed by gentle shaking or tapping the cuvette.

**pH** sample .....2.5 – 11  
pH values outside this range cause very low-bias results.

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