

# neodisher<sup>®</sup> Entschäumer S

Liquid component with anti-foaming action  
for use in special washing machines



## xMain fields of application:

Use as an additive to suppress or prevent foam in special washing machines in medical laboratories (blood investigation, serology), the foodstuffs industry, the cosmetics industry, and industrial laboratories as well as in special washing machines for operating tables.

## Characteristics:

neodisher Entschäumer S prevents the development of foam during the removal of foam-generating residues such as surfactants, proteins, or soaps. neodisher Entschäumer S can be used as a component with all neodisher cleaning agents for the preliminary or main cleaning, and it is effective even at low temperatures.

## Application and dosage:

0.05 - 0.1 ml/l added manually or using suitable dosing devices.

Avoid overdosing, otherwise precipitation and deposition of silicone residues on the surfaces of the items being cleaned may result.

The addition should be made either during or immediately after water intake in washing machines with water change. In the case of tank machines, addition is effected continuously via a suitable metering pump.

Do not mix with other products.  
Only for professional use.

## Technical data:

Specific gravity (20 °C): 1.0 g/cm<sup>3</sup>

pH-range: neutral

## Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:

< 5 % non-ionic surfactants

also preservation agent (Benzylhemiformal)

## CE-marks:

neodisher Entschäumer S conforms to the requirements of the legislation for medical devices according to the European directive 93/42/EC.

## Storage information:

Storage temperatures should be between 5 °C and 30 °C.

## Hazard warnings and safety advice:

neodisher Entschäumer S has not to be classified according to the Dangerous Preparation Directive 1999/45/EC.

Dispose only when container is empty and closed. For disposal of product residues, refer to Material Safety Data Sheet.

For further safety information see EC-Safety data sheets. These are available e.g. under [www.drweigert.de](http://www.drweigert.de) under the headline "Service".

MB 4301 GB/2-2  
02/06



Chemische Fabrik Dr. Weigert GmbH & Co. KG  
Mühlenhagen 85, D – 20539 Hamburg

Telefon: (040) 789 60 - 0  
Telefax: (040) 789 60 - 120

E-Mail: [info@drweigert.de](mailto:info@drweigert.de)  
Internet: [www.drweigert.de](http://www.drweigert.de)

Die Angaben dieses Merkblattes beruhen auf unseren derzeitigen Kenntnissen und Erfahrungen. Sie befreien den Verwender nicht von eigenen Prüfungen und Versuchen. Eine rechtlich verbindliche Zusicherung bestimmter Eigenschaften kann hieraus nicht abgeleitet werden.

With the above information, which is appropriate to our current knowledge we describe our product regarding possible safety necessities, but we do not involve any quality description or promise certain properties.