Shaking Water Bath GLS Aqua Plus

Operating Manual

EN

FR

DE

IT

ES

Grant Instruments, based near Cambridge, England, is an independent, privately owned company, founded in 1951 by Peter Ward and Cecil Chapman.

Grant is a world leader in the manufacture and design of equipment for sample preparation, scientific analysis, data acquisition and data analysis providing solutions to the global scientific and industrial markets.

Standards Compliance and Quality

Grants' brand and reputation are based around quality, reliability and accuracy. We ensure our products stringently meet all necessary international safety standards.

We pay particular attention to the safety testing of products and remain at the forefront of the product safety standard for laboratory equipment IEC 61010-1. The company is committed to operating its safety test laboratory in accordance with the requirements of ISO 17025.

Grant operates a Quality Management System that complies with the requirements of BS EN ISO 9001:2008.

Beyond compliance to the standard, Grant is committed to continually improving in everything we do; with particular emphasis on understanding what matters to our customers and suppliers, and designing our systems and work to meet their needs.

If you have any feedback on Grant's products or services we would like to hear from you. Please send all feedback to:

Quality Manager
Grant Instruments (Cambridge) Ltd
Shepreth
Cambridgeshire
SG8 6GB
England

Tel: +44 (0) 1763 260 811
Fax: +44 (0) 1763 262 410
E-mail: feedback@grant.co.uk

Contents

1.0	Use of products	3
2.0	How to use this operating manual	3
3.0	Safety information	3
3.1	Safety compliance	3
3.2	Safety symbols	4
3.3	Safety warnings	4
4.0	Operating instructions	5
4.1	Unpacking instructions	5
4.2	Assembly of the equipment and components	5
4.3	Installation	6
4.4	Electrical supply	6
5.0	Operating procedures	7
5.1	Operation	7
5.1.1	Water level	7
5.1.2	Operation above 60°C	7
5.1.3	Use as a thermostatic bath without shaking	7
5.1.4	Emptying the baths	7
5.2	Using the GLS Aqua Plus	8
5.2.1	Bath controls	8
5.2.2	Setting the control temperature	8
5.2.3	Setting the shaking speed	9
5.2.4	Starting & stopping shaking	9
5.2.5	Setting the sample protection thermostat	9
5.2.6	Configuring a bath preset	10
5.2.7	Running a bath preset	10
5.2.8	Setting a countdown period	10
5.2.9	Cancelling the timer	11
5.2.10	Locking & unlocking the control panel	11
5.2.11	Single point calibration	11
5.2.12	Twin point calibration	11
6.0	Technical specifications	13
7.0	Technical Tips	13
7.1	Which water should you use in your bath?	13
7.2	How to prevent rust in water baths	14
7.3	How to prevent algae and bacteria?	14
8.0	Warranty information	15
9.0	Maintenance and service	15
9.1	Cleaning	15
9.2	Fuses	15
9.2.1	Replacing fuses on GLS Aqua Plus	15
9.2.2	Routine safety tests	16
9.3	Service	16
10.0	Optional accessories	16
11.0	Troubleshooting	16
12.0	Contact Grant Instruments	17
13.0	Compliance	17

1.0 Use of products

The following products are covered by this operating manual:

- GLS Aqua (12 and 18) Plus

The products listed above are a general purpose series of thermostatically controlled shaking water baths designed for indoor laboratory use by a professional user.

2.0 How to use this operating manual

This operating manual will allow you to unpack, set-up and operate this water bath correctly and safely. Important safety information, symbols and warnings are listed below and should be read carefully. Section 4 gives information about how to unpack and install the product correctly. Section 5 gives generic operating information in the first section followed by specific operating details for the GLS Aqua Plus baths. Product technical specifications and tips are provided in the sections 6 and 7. The warranty for this water bath is detailed in section 8 but to register you should return the warranty card or complete the on-line registration form at www.grant.co.uk.

If there is a technical matter that this operating manual does not address, or any other question concerning this product please contact Grant Instruments or your local distributor who will be able to provide any additional information.

A laminated quick start guide is provided with the baths as quick reference guide but should not be used until the full user manual has been read.

3.0 Safety information

3.1 Safety compliance

Grant water baths meet the requirements of international safety standard IEC 61010-2-10 and IEC 61010-2-51 and national standards based on them including:

UL 61010A-2-010

UL 61010-2-051-04

CAN/CSA-C22.2 NO. 61010-2-010-04

CAN/CSA-C22.2 No. 61010-2-051

3.2 Safety symbols

The symbols below are marked on the equipment to indicate:

Caution: Surfaces and water can be hot during and after use

Read this manual before using the bath

Important safety warning

3.3 Safety warnings

Read the whole of these instructions. Safety may be impaired if they are not followed.

If the equipment has been transported or stored in cold or humid conditions, condensation may form inside it. If that could have happened, allow time (at least 2 hours at room temperature) for the condensation to evaporate before using the equipment.

Do not use the bath to heat any material that could cause a fire or any other kind of hazard.

Do not use the equipment in an area where there are hazardous or explosive chemical mixtures or atmospheres.

If a potentially hazardous liquid is spilt onto the equipment, disconnect it from the power supply and have it checked by a competent person. It is the user's responsibility to carry out appropriate decontamination if hazardous material is spilt on the equipment.

Before emptying a bath, allow the water temperature to fall to a safe level. For 12 and 18 litre baths, empty the bath before moving it.

The shaking speed will revert to previously set value when reconnected after power interruption.

The protection provided by the equipment may be impaired if it is used with accessories not provided or approved by Grant.

The bath is for use only with water as the bath liquid. Make sure that it cannot become contaminated by other liquids. We recommend the use of de-ionised water.

Before first switching on the bath please remember to fill the bath with water. Switching the bath on dry will damage the heater and could invalidate the product warranty.

4.0 Operating instructions

4.1 Unpacking instructions

Standard equipment includes:

- Thermostatic bath
- Mains cord with plug
- Gabled polycarbonate lid
- Drive magnet
- Shaking trolley
- Operating manual
- Quick start guide


Remove packing materials carefully, and retain them for future shipment or storage of the equipment.

4.2 Assembly of the equipment and components

The water bath has four main components, the bath, the lid, the tray (supplied separately) and the shaking mechanism. The lid should only be lifted by the handle provided as other parts can become hot during use. It also has a vent/thermometer hole - this hole should not be sealed as pressure could build up inside the bath.

The shaking mechanism and tray should be fitted as follows

1. Put the drive magnet into the hole in the base of the tank.
2. Place the shaking trolley in the tank with the axles parallel to the front panel, and with the peg of the drive arm in the hole in the top of the drive magnet.
3. Place the shaking tray on the shaking trolley (see diagram below showing the universal plain tray).

4.3 Installation

Place the water bath on a level, non-combustible surface. Ensure that the mains plug and the switch are easily accessible.

4.4 Electrical supply

Check that the supply voltage marked on the serial number label, and the type of mains plug, are correct for your mains supply outlet, which must have a ground connector.

To disconnect the equipment from the mains supply, remove the mains plug from the mains supply outlet.

5.0 Operating procedures

5.1 Operation

5.1.1 Water level

Ensure that the baths are used with the shaking tray in place and that the water always covers the tray. The maximum water level should not be higher than the swage line around the top of the bath (around 2.5cm from the top). These criteria apply both when there are no vessels in the bath and with the maximum contents. Note that water level may have to be reduced at higher shaking speeds to prevent spillage and stalling of the mechanism.

Avoid letting the bath run dry. In the event this does happen, the safety cutout will trip and disconnect the heater. If this does happen, unplug the bath and have the cut-out reset by a competent person.

5.1.2 Operation above 60°C

The lid must be used above 60°C to maintain proper temperature control and to ensure that the water temperature reaches the set point

The lid will also prevent excessive evaporation that requires the bath to be filled more often and will save energy

5.1.3 Use as a thermostatic bath without shaking

If the shaking trolley and drive magnet are lifted out, a shaking bath can be used as a simple thermostatic bath. When used in this way, if flat bottomed vessels or objects are to be placed in the bath, always use a raised base tray (optional accessory) to avoid damage to the under-tank heater.

5.1.4 Emptying the baths


Before emptying any bath allow the water temperature to fall to a safe level and take reasonable precautions to prevent accidental spillage.

Drain taps are included in both baths to allow convenient emptying. These baths should be emptied using the drain tap prior to moving. To empty the bath using the drain tap, push the supplied drain insert into the drain tap (see picture below). Note that the water will begin to empty from the bath as soon as the drain insert is fully engaged. A length of hose can be added to the barbed end of the drain insert if required.

5.2 Using the GLS Aqua Plus

5.2.1 Bath controls

5.2.2 Setting the control temperature ()

The water temperature of the bath can be set using the control keys.

1. Turn the knob of the sample protection thermostat to maximum.
2. Whilst the display is showing the bath temperature press either the '+' or '-' key. This will cause the display to flash indicating that it can be set.
3. Use the '+' or '-' to set the desired temperature.
If no key is pressed for 4 seconds then the display will revert back to showing the bath temperature and the set temperature will remain at its original value.
4. Press the enter key to store the requested value and the display will revert to showing the bath temperature.
The water temperature will change to the new set value - during heating or cooling to the set point the last digit will change to a rising or falling bar until the actual temperature is within one degree of the set value.
5. When the temperature has stabilised measure the temperature with a thermometer and adjust the set temperature again if necessary.
6. Set the sample protection thermostat as detailed below if needed.

5.2.3 Setting the shaking speed

The shaking speed of the bath can be set using the control keys

1. Press the enter key to show ' ' .
2. Use the '+' or '-' keys to navigate to the animated tray option: 
3. Press the enter key to select the shaking speed; the display will show the present shaking speed.
4. Use the '+' or '-' keys to set the shaking speed in strokes/min.
The shaking will respond as the value changes.
5. Press the enter key to set the shaking speed to the new value.
If no key is pressed for 4 seconds then the display will revert back to showing the bath temperature and the shaking speed will remain at its original value.

The maximum speed for a particular application will depend on number and size of vessels, and the depth of the water. The shaking trolley is driven by a drive magnet under the tank, magnetically coupled to the trolley through the tank. As the speed is increased a point is reached where the magnetic coupling de-latches, causing the trolley either to move erratically or to stop.

5.2.4 Starting & stopping shaking

The shaking can be started or stopped at any time using the '↑' key. The stop shaking function stops the tray immediately whereas the start shaking function ramps up the shaking slowly to the desired speed to prevent stalling. When the '↑' key is pressed to start shaking the current speed setting is displayed briefly on the display.

5.2.5 Setting the sample protection thermostat

The bath is fitted with an adjustable sample protection thermostat to protect the sample. This device is not a safety feature. The sample protection thermostat can be set to be a few degrees above the bath set temperature.

1. Turn the knob of the sample protection thermostat to maximum.
2. Set the control temperature to 2°C above the desired operating temperature and wait for the temperature to stabilise.
3. Turn the knob of the sample protection thermostat slowly anti-clockwise until a click is heard and ' ' is displayed.
The audible alarm will sound continuously to warn you that the bath is controlling using the sample protection thermostat.
4. Re-set the control temperature to the desired temperature using the '+' and '-' keys.
The audible alarm will sound continuously until the primary control system is re-activated.

In the unlikely event of failure of the primary temperature control system, the sample protection thermostat will maintain the water in the bath at a temperature a few degrees above the set temperature and with greater fluctuations. The display will cycle between ' ' and the actual bath temperature and the audible alarm will sound intermittently. If this happens, first check that the sample protection thermostat is not set to a lower temperature than the primary control system. If that is **not** the reason, the bath can continue to be used without compromising the safety of persons or the surroundings until current work is completed. However, it is recommended to have the bath checked by a competent person as soon as conveniently possible.

5.2.6 Configuring a bath preset

Each bath contains two presets which can be configured to different set temperatures and shaking speeds to allow the bath to be conveniently run under frequently used conditions. Use the method below to configure preset 1; other presets can be set in a similar manner.

1. Press preset key '1' to activate the preset temperature setting. The display will flash ' ' and then the current preset set temperature.
2. Use the '+' or '-' to set the desired temperature.
If no key is pressed for 4 seconds then the display will revert back to showing the bath temperature and the set temperature will remain at its original value.
3. Press preset key '1' to save the preset temperature.
4. Press preset key '1' twice to activate the preset shaking speed setting. The display will flash ' '; then  and finally the current preset shaking speed
5. Use the '+' or '-' to set the desired shaking speed.
If no key is pressed for 4 seconds then the display will revert back to showing the bath temperature and the set temperature will remain at its original value.
6. Press preset key '1' to save the preset shaking speed.

5.2.7 Running a bath preset

To run bath preset 1:

1. Press preset key '1' to activate the preset. The display will flash ' ' and then the current preset set temperature.
2. Press the enter key to run the preset; the bath will heat (or cool) to the requested temperature.

5.2.8 Setting a countdown period ()

The bath includes a countdown timer in minutes which triggers an audible alarm on completion.

1. Press the enter key to show ' '.
2. Use the '+' or '-' keys to navigate to ' '.
3. Press enter key to select the countdown timer; the display will flash ' '.
4. Use the '+' or '-' keys to set the countdown timer in minutes.
5. Press the enter key to confirm the countdown timer.
If a timer has been set then a dot in the bottom right corner of the display will flash during the countdown.

When the countdown timer reaches zero the audible alarm will sound and the display will alternate between the water temperature and zero. The alarm can be accepted by pressing the enter key.

5.2.9 Cancelling the timer

The countdown timer can be cancelled by using the instructions above but using the '-' key to reset the value to zero. Use the enter key to confirm that you wish to clear the timer.

5.2.10 Locking & unlocking the control panel

The keys on the front panel can be locked to prevent settings being accidentally changed or modified. The front panel will display ' ' if any buttons are pressed and the control panel is locked.

The control panel can be locked by pressing the enter and '+' keys simultaneously for at least 3 seconds. The display will show ' ' to confirm that the keypad has locked. The control panel can be unlocked by pressing the enter and '-' keys simultaneously for at least 3 seconds. The display will show ' ' to confirm that the keypad is operational.

5.2.11 Single point calibration ()

The display of set temperature and actual temperature is accurate to within 1°C at 37°C but may be up to 5°C out at maximum temperature. Single point calibration allows the bath to be adjusted to be accurate to $\pm 0.1^\circ\text{C}$ at a single temperature.

1. Press the enter key to enter the menu. The display will show ' '.
2. Use the '+' or '-' keys to navigate to ' '.
3. Press enter key to select calibration; the display will flash with the current offset value.
4. Use the '+' or '-' keys to set an offset correction value (positive or negative) as measured by a thermometer.
If no key is pressed for 8 seconds then the display will revert back to showing the bath temperature and the calibration will remain at its original value.
5. Press the enter key to store the new corrected set temperature.
The display will show the corrected bath temperature.

5.2.12 Twin point calibration ()

The bath can also be calibrated using a twin point method to give you a linear correction over the full temperature range. The process consists of three steps:

1. Set low point calibration value.
2. Set high point calibration value.
3. Confirm calibration settings.

The bath is not calibrated until all three steps have been completed. The greater the difference between the low and high points, the better the calibration will be.

Low calibration point ()

1. Set the bath to a temperature between 20°C and 40°C and allow the bath to stabilise, indicated by the display not fluctuating by more than 0.2°C.
2. Measure the temperature of the water with a calibrated thermometer and make a note of this reading.
3. Press the enter key. The display will show ' '.
4. Use the '+' or '-' keys to navigate to ' ' and press the enter key.
5. Use the '+' or '-' keys to navigate to ' ' and press the enter key.
6. Using the '+' or '-' keys adjust the display until the value displayed is the same as the reading taken earlier.

7. Press the enter key to confirm the value and the bath will display ' '. Wait for the bath to return to normal mode.

High calibration point ()

1. Set the bath to a temperature between 70°C and 90°C, allow the bath to stabilise, indicated by the display not fluctuating by more than 0.2°C.
2. Measure the temperature of the water with a calibrated thermometer and make a note of this reading.
3. Press the enter key. The display will show ' '.
4. Use the '+' or '-' keys to navigate to ' ' and press the enter key.
5. Use the '+' or '-' keys to navigate to ' ' and press the enter key.
6. Using the '+' or '-' keys adjust the display until the value display is the same as the reading taken earlier.
7. Press the enter key to confirm the value and the bath will display ' '. Wait for the bath to return to normal mode.

Setting the calibration ()

1. Press the enter key. The display will show ' '.
2. Use the '+' or '-' keys to navigate to ' ' and press the enter key.
3. Use the '+' or '-' keys to navigate to ' ' and press the enter key to apply the new low and high calibration values.

The bath is now calibrated and the display will read the corrected value.

Restoring factory calibration ()

1. Press the enter key. The display will show ' '.
2. Use the '+' or '-' keys to navigate to ' ' and press the enter key.
3. Use the '+' or '-' keys to navigate to ' ' and press the enter key to restore the factory settings.

6.0 Technical specifications

Operating conditions

Ambient Temperature	5 to 40°C
Maximum relative humidity	80% R.H. in room temperatures up to 31°C decreasing linearly to 50 % R.H. at 40°C
Altitude above sea level	Up to 2,000 m (6,500 ft)
Operating Environment	Indoor use only

Electrical details

Mains supply: 230V @ 50/60 Hz or 120V @ 50/60 Hz

Pollution degree: 2

Installation Category: II

Note: Mains supply voltage fluctuations are not to exceed $\pm 10\%$ of the nominal supply voltage.

Models	Capacity (L)	Current Rating (A)	
		120V	230V
GLS Aqua 12 Plus	12	5.8	2.8
GLS Aqua 18 Plus	18	8.5	5.4

GLS Aqua Plus bath performance

Temperature range	5°C above ambient to 99°C
Shaking speed range	40 to 400 strokes/min
Display (also used for setting)	10.0 to 99.0°C in steps of 0.1°C
Temperature stability to DIN 12876	$\pm 0.2^\circ\text{C}$

The maximum shaking speed for a particular application will depend on the number and size of vessels in the shaking tray, and the depth of the water.

7.0 Technical Tips

7.1 Which water should you use in your bath?

For the long-term reliability of water baths it is important to use oxygenated water that is free from ions and minerals that can cause corrosion of stainless steel. We recommend the use of distilled water and de-ionised water from modern ion exchange systems that do not use salt back flushing to regenerate the ion-exchange cartridges.

Stainless steel is protected from corrosion by a layer of chromium oxide. If the layer is damaged, oxygen present in water can reform the oxide layer. If the water is still or de-oxygenated, and the oxide layer is damaged, ions can corrode the stainless steel tank. If a water bath has been unused for some time, or water boiled, we recommend changing to fresh distilled water or correct de-ionised water.

Water normally contains calcium or magnesium ions. De-ionised water has most ions removed as indicated by its conductivity level; the purer the water the lower the conductivity. It is important to use only de-ionised water from an ion exchange system with replaceable cartridges. Do not use de-ionised water generated from an ion-

exchange system that incorporates a salt back-flush system to regenerate the ion-exchange resin as this can leave sodium ions that are very corrosive to stainless steel.

7.2 How to prevent rust in water baths

Most Grant tanks, as well as immersed parts, are made from type 304 stainless steel, an extremely versatile general purpose grade of stainless steel. It is the excellent forming characteristic that has made this grade dominant in the manufacture of laboratory and industrial water baths, as well as domestic sinks and saucepans. Type 304 stainless steel is highly suitable for applications where hygiene is important; it exhibits good heat resistance and excellent resistance to corrosion.

However, despite resistance to general surface corrosion, stainless steel is susceptible to specific types of corrosion, in particular pitting (small pin hole style corrosion) and stress corrosion cracking. It can also undergo general corrosion in specific environments, such as one containing hydrochloric or sulphuric acids.

Stainless steel is protected by its high content of alloying elements, primarily chromium and nickel. Chromium is the most important with respect to corrosion resistance, although the nickel assists in allowing the chromium to do its job. The chromium forms an oxide layer on the surface of the steel, which inhibits further oxidation. This layer adheres extremely well to the metal substrate, but it is essential that it remains intact, and must be protected from various forms of damage.

If the surface chromium oxide layer becomes damaged, oxygen present in water can partially reform the oxide layer, so it is advisable to ensure that water is always fresh and well oxygenated. Baths that will be out of use for an extended period should be emptied, and all moisture should be wiped from the bottom of the tank.

In some cases a brown layer may appear on the surface of a stainless steel tank. In most of these cases this is not rust, but it may be a surface deposit of minerals from the local water supply, or ferrous particles or salts that have fallen into the tank. These surface deposits can usually be removed by using a household cleaner such as Duraglit or Silvo metal polish.

7.3 How to prevent algae and bacteria?

Water baths provide the ideal environment for the growth of micro-organisms. If left uncontrolled the growth of these organisms can result in a range of serious problems and health risks from pathogenic bacteria.

The growth of algae on the surface of parts will cause biofouling which can reduce performance.

Micro-organisms that produce acidic metabolic by-products can cause bio-corrosion by depolarisation of metal surfaces.

There are a number of biocides available on the market.

8.0 Warranty information

When used in laboratory conditions and according to these working instructions, this product is guaranteed for THREE YEARS against faulty materials or workmanship.

Extended warranty for years four and five can be purchased by contacting our sales department at labsales@grant.co.uk.

9.0 Maintenance and service

No routine maintenance is required except for cleaning. There are no user serviceable parts inside the unit.

9.1 Cleaning

Clean the outside of the equipment with a damp cloth, using water only. Do not use chemical cleaning agents. Before using any other cleaning or decontamination method, check with Grant Instruments or your local representative to make sure that the proposed method will not damage the equipment.

Scale on immersed parts can be removed using chemical de-scaling products designed for use on kitchen equipment that have metal parts. De-scaling products may be toxic and manufacturer's instructions should always be followed.

9.2 Fuses

Fuses are Littelfuse 3AB 314 series, fast-acting, high breaking current (max breaking current at least 750 A); dimensions are 1.25 inch long, 0.25 inch diameter. Replace fuses only by the same type and rating (250V).

Models	Fuse Rating (A)	
	120V	230V
GLS Aqua 12 Plus	10	10
GLS Aqua 18 Plus	15	15

9.2.1 Replacing fuses on GLS Aqua Plus

Disconnect the unit from the power supply

1. Remove the mains lead from the socket at the back of the bath.
2. Press down the fuse drawer catch and pull out the fuse drawer. Check the fuse(s) and replace if necessary using the specified fuse type and rating.
3. Push back the drawer and replace the mains lead.

9.2.2 Routine safety tests

If routine tests are to be made, we recommend a test of the integrity of the protective earth conductor and an insulation test at 500 V DC. Routine flash tests are **not** recommended for any electrical equipment, because repeated high voltage tests degrade insulation materials.

9.3 Service

If service is required, switch off the unit and contact Grant Instruments or your local representative for repairs.

Service Department
Grant Instruments (Cambridge) Ltd
Shepreth
Cambridgeshire
SG8 6GB
England.

Tel: +44 (0) 1763 260 811
Fax: +44 (0) 1763 262 410
E-mail: labservice@grant.co.uk

10.0 Optional accessories

A full listing of product accessories and options is available in the Grant Scientific Reference Catalogue (a copy of which is available upon request) and on the Grant website at www.grant.co.uk.

11.0 Troubleshooting

No lights or display - check power source, mains switch and fuses.

Water is not reaching temperature - check that the set temperature and the sample protection thermostat are set correctly.

Display shows ' ' - check that the sample protection thermostat is set correctly.

Display temperature is different from actual temperature - use single () or twin point calibration () to correct this.

There is no last digit on the display - allow the bath to stabilise.

Temperature of water fluctuates - check that the sample protection thermostat is set correctly.

EN

12.0 Contact Grant Instruments

At Grant we are continuously trying to improve the performance we offer our customers. If you have any feedback on Grant's products or services we would like to hear from you. Please send all feedback to:

Quality Manager
Grant Instruments (Cambridge) Ltd
Shepreth
Cambridgeshire
SG8 6GB
England

Tel: +44 (0) 1763 260 811
Fax: +44 (0) 1763 262 410
E-mail: feedback@grant.co.uk

13.0 Compliance

WEEE directive

Grant Instruments complies fully with the Waste Electrical & Electronic Equipment (WEEE) regulations 2006. We are a member of the B2B compliance scheme (Scheme Approval Number WEE/MP3338PT/SCH), which handle our WEEE obligations on our behalf. Grant Instruments have been issued with a unique registration number by the Environmental Agency, this reference number is WEE/GA0048TZ.

For information regarding WEEE collections in the UK please contact our B2B Compliance Scheme directly on 01691 676 124
For other countries please contact your equipment supplier.

For General WEEE information please visit: www.b2bcompliance.org.uk

RoHS directive

All the products covered by this manual comply with the requirements of the RoHS Directive (Directive 2002/95/EC).

Electrical safety and electromagnetic compatibility

All the products covered by this manual comply with the requirements of the Low Voltage Directive (2006/95/EC) for electrical safety and the EMC directive (2004/108/EC) for electromagnetic compatibility. See the Declaration of Conformity on the inside back page.

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Quality Manager
Grant Instruments (Cambridge) Ltd
Shepreth
Cambridgeshire
SG8 6GB
Angleterre

Tél. : +44 (0) 1763 260 811
Fax : +44 (0) 1763 262 410
E-mail : feedback@grant.co.uk

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Declaration of Conformity

Manufacturer	GRANT INSTRUMENTS (CAMBRIDGE) LTD, Shepreth, Cambridgeshire SG8 6GB
Equipment Name	GLS AQUA PLUS: GLS AQUA 12, GLS AQUA 18, GLS AQUA 12L, GLS AQUA 18L
Equipment Description	Electronically controlled thermostatic baths
Directives	EMC Directive 2004/108/EC LVD Directive 2006/95/EC
Including Accessories	Mains cord-set, plastic lid and tray
CE mark first applied	2010

Applied Standards:	
Harmonised Standards:	BS EN 61326-1:2006 Electrical Equipment for measurement, control and laboratory use - EMC requirements- Part 1: General requirements BS EN61010 Part 2-010:2003 Safety requirements for electrical equipment for measurement, control and laboratory use; particular requirements for laboratory equipment for the heating of materials BS EN61010 Part 2-051:2003 Safety requirements for electrical equipment for measurement, control and laboratory use; particular requirements for laboratory equipment for mixing and stirring

The product complies with the above Directives when used with the supplied mains cord-set; it may not comply if an alternative cord-set is used.

I confirm that this apparatus conforms to the requirements of the above Directive(s)

.....

Ludo Chapman
Managing Director
Grant Instruments (Cambridge) Ltd.

Dated.....30/4/2010.....

USA

This device complies with Part 15 of the FCC rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

EN **IMPORTANT!** Before using this water bath for the first time, please clean the inside tank with general household cleaner, rinse with de-ionised water, then dry and wipe clean thoroughly. This will help prevent bath corrosion. Please refer to “How to prevent rust in water baths” section of the operating manual for more detail.

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