



2019 Chemical Filtration Pump

robust PTFE diaphragm pump

- **Suitable for basic and acid media**
- **Reliable and compact due to design**
- **Low maintenance even in aggressive applications**

2019C-02

Chemical Duty Filtration Pump

Economical Model 2019 vacuum pump is durable and chemical resistant. Use Model 2019 for vacuum filtration, dessication, aspiration, degassing, and other mid-range vacuum applications. Free air flow 31 l/min; maximum vacuum 200 mbar. Ideal for single funnel or up to 6 holder manifold filtrations – for organic and light acid/base (e.g. acetic acid) solutions.

Pump can also be used to pressurize up to 18 PSIG. Comes standard with UK and Schuko line cord, plug and convenient power switch.



Design Features:

- PTFE wetted parts
- Low noise emission
- Direct drive
- Continuous operation
- No oil contamination

Applications:

- Filtration
- Vacuum drying
- Dessicators
- Concentration

Scope of Delivery:

Pump with Schuko, and UK plug leads

Certification:

CE marked, UL/CSA approved

Accessories:

Silicone Hose 6 mm ID
CAT NO. 828374

Particulate Filter
CAT NO. 112555-04

Technical Features/Ordering Information:

Model	Pumping Speed	Ultimate Pressure	Tubing Needed	Dimensions W/D/H	Weight	CAT. No.
2019	31 l/min	200 mbar	6 mm	264/127/222 mm	4.6 Kg	2019C-02