

# ACCESSORIES



## Angle rotor, 12-place



∠35°  
n = 5,000 min<sup>-1</sup>  
max. RCF 2,879

**Cat. No. 1133**

capacity in ml	5	6	7
∅ x L in mm	12x75	12x82	12x100
<b>Cat. No.</b>	<b>0553</b> <sup>1)</sup>	<b>0501</b> <sup>1)</sup>	<b>0578</b> <sup>1)</sup>



**rotor**  
**Cat. No. 1133**

<b>Cat. No.</b>	-
boring ∅ x L in mm	12.5 x 66
tubes per rotor	12
max. RCF <sup>1)</sup>	2,879
radius in mm	103
run-up in sec	8
run-down in sec, braked	10


## Angle rotor, 6-place



∠35°  
n = 6,000 min<sup>-1</sup>  
max. RCF 4,025

**Cat. No. 1137**

capacity in ml	7	15	25	50	9–10	10	1.6–5	4–7	15	50	30	50
∅ x L in mm	12x100	17x100	24x100	34x100	16x92	15x102	13x75	13x100	17x120	29x115	26x95	29x107
<b>Cat. No.</b>	<b>0578</b> <sup>1)</sup>	<b>0518</b> <sup>1)</sup>	<b>0519</b> <sup>1)</sup>	<b>0521</b> <sup>1)</sup>	<b>blood collection / urine tubes</b>				<b>0509</b>	<b>0513</b>	<b>0545</b>	<b>0546</b>



**rotor**  
**Cat. No. 1137**

<b>Cat. No.</b>	<b>1632</b>	<b>1635</b>	<b>1633</b>	-	<b>1635</b>				<b>1631</b>	<b>1641</b>	<b>1633</b>	<b>1634</b>
boring ∅ x L in mm	13x92	17.5x95	26x88	35x96	17.5x95				17x98	30x98	26x88	29x95
tubes per rotor	18	6			6				6	3	6	
max. RCF <sup>1)</sup>	3,944	3,783	3,703	4,025	3,783		2,978	3,783	3,824	3,824	3,703	3,904
radius in mm	98	94	92	100	94		74	94	95	95	92	97
run-up in sec	20											
run-down in sec, braked	17											

<sup>1)</sup> Please note that the RCF values indicated refer only to rotor performance. The max. permissible RCF of tubes used should be verified with the individual manufacturers. The max. RCF for glass tubes annotated with footnote <sup>1)</sup> is 4,000.