



The Dune series offers all of the necessary features for quick and easy weighing. Every Dune is equipped with an AC adapter, auto calibration through the keypad, low battery indicator, and overload indicator. With its alternate battery operation, splash-proof keypad, and sturdy plastic construction throughout, the Dune delivers accurate readings in field, classroom, workplace or lab, a truly versatile and portable product.

Streamlined design for basic weighing in the lab, classroom or field




## Features

- High-contrast display is easily visible from a distance
- Colour-coded keys facilitate quick recognition of the most frequently used buttons
- Lightweight design simplifies transport and portability
- Overload protection helps prevent damage to internal components
- Sealed keypad protects against dirt and spills
- Large platform removes easily for speedy cleaning
- Auto power-off to save energy
- Zero-tracking feature ensures display returns to zero reading
- Operates with AC adapter (included) or batteries (not included)

## Applications

- Weighing



DCT 302

## About ADAM:

For more than 40 years, professionals worldwide have relied on Adam Equipment for an extensive selection of dependable and affordable weighing equipment. Adam provides the right balance of speed, performance and value.

Model	DCT 302	DCT 201	DCT 601	DCT 2001	DCT 2000	DCT 5000
<b>Capacity</b>	300g	200g	600g	2000g	2000g	5000g
<b>Readability</b>	0.01g	0.1g	0.1g	0.1g	1g	2g
<b>Repeatability (S.D.)</b>	0.01g	0.1g	0.1g	0.1g	1g	2g
<b>Linearity (+/-)</b>	0.02g	0.2g	0.2g	0.2g	2g	4g
<b>Pan Size</b>	100mm ø	145x145mm				
<b>Weighing Units</b>	g, oz		g, lb, oz			
<b>Tare Range</b>	Full range					
<b>Stabilisation Time (sec)</b>	2					
<b>Calibration</b>	External calibration					
<b>Display</b>	LCD with 15mm -high digits					
<b>Power Supply</b>	12VAC 50/60Hz 150mA adapter / 6 x AA batteries					
<b>Operating Temperature</b>	0° to 40°C					
<b>Housing</b>	ABS plastic					
<b>Draught Shield</b>	Square 145x148x120mm	N/A				
<b>Overall Dim.</b>	147x208x150mm (wx dxh)					
<b>Net Weight</b>	0.5kg					

## Accessories

Item number	Description
700660290 .....	Calibration certificate
700100012.....	F1 200g Calibration Weight
700100010 .....	F1 2kg Calibration Weight
700100030 .....	M1 2kg Calibration Weight
700100060 .....	M1 500g Calibration Weight
700100072.....	M1 200g Calibration Weight

Distributed and Supported in the UK by:



24 Norman Way Industrial Estate, Over, Cambridge, CB24 5WE  
 Tel: +44(0)1954 233 100 Fax: +44(0)1954 233 101  
 Email: sales@camlab.co.uk Web: www.camlab.co.uk