



RAFR05203 Advanced Under Counter Freezer

Part Number	RAFR05203
Capacity (Litres)	150
Style	Under Counter
Exterior Dimensions (HxWxD)	820 x 600 x 640mm
Energy Consumption	1.3 kWh/24h
Operating Temperature	-20°C
Temperature Set Point Range	-18°C to -22°C
Interior Construction	Stainless Steel
Exterior Construction	Stainless Steel
Shelves	2
Shelf dimensions (W xD)	490 x 430mm
Fixed Plate Shelves	✓
Door Hinging	Right hand only
Door Style	Solid
Rear Castors	✓
Door Lock	✓
Access Port	Optional Extra
Remote Alarm Contacts	Optional Extra
Digital Temperature Control	✓
Digital Temperature Display	✓
High Temperature alarm	✓
Low Temperature alarm	✓
Door Open alarm	✓
Manual defrost	✓
Environmentally friendly insulation	✓
Eco Friendly Refrigerant	✓
Ammonia Free	✓

Robustly constructed from Stainless Steel this freezer from the **Labcold** Advanced range is designed to provide close temperature control for temperature sensitive products.



Insulated with Zero ODP Injected Polyurethane Insulation, the refrigerant is also ammonia, HCFC and CFC free.

Fitted with a sturdy door lock, this Advanced freezer is ideal for demanding environments such as busy laboratories and features a manual defrost putting you in control.

Remote alarm contacts and an access port can be fitted at time of manufacture at an additional cost.