

Compatibility Chart UV-Cuvettes micro

Status as of: 1107

	Center height mm	Compatibility	Remarks
Amersham Biosciences/Biochrom/Pharmacia			
GeneQuant	15	●	
Ultrospec 1100 pro	15	●	
Ultrospec 2000	15	●	
Ultrospec 2100 pro	15	●	
Ultrospec 3000	15	●	
Ultrospec 3100 pro	15	●	
Ultrospec 3300 pro	15	●	
Ultrospec III	15	●	
Analytik Jena AG			
SPECORD 40/30	8.5 / 15	●	Please use adjustable cuvette holder!
SPECORD 50	8.5 / 15	●	
SPECORD 200	8.5 / 15	●	
SPECORD 205/210	8.5 / 15	●	
SPECORD S 100/S 100 Bio	8.5 / 15	●	
SPEKOL 1100/1200	8.5 / 15	●	
Biochrom Ltd./WPA (see Amersham)			
Lightwave S200	15	●	
Lightwave/Biowave S2100	15	●	
UV 1101	15	●	
UV 1101T	15	●	
Eppendorf AG			
BioPhotometer	8.5	●	
Hitachi			
all types	-	-	
JASCO			
V-550	8.5	●	
Krüss AG			
UV6500	8.5	●	
Perkin Elmer Instruments			
Lambda 2	15	●	Please use the bottom attachment for UV-cuvettes micro! For further informa- tion, please contact info@brand.de
Lambda 10	15	●	
Lambda 11	15	●	
Lambda 12	15	●	
Lambda 14	15	●	
Lambda 20	15	●	
Lambda 25	15	●	
Lambda 35	15	●	
Lambda 40/Lambda 40 P	15	●	
Lambda 45	15	●	
Lambda 650	15	●	
Lambda 750	15	●	
Lambda 800	15	●	
Lambda 850	15	●	
Lambda 900	15	●	
Lambda 950	15	●	
Shimadzu			
UV-120-02	15	●	
UV-160 A	15	●	
UV-1601	15	●	
ThermoSpectronic			
AquaMate	15	●	Recommended accessory parts are available from the photometer's manufacturer.
BioMate 3	8.5	●	
BioMate 5	15	●	
Evolution 60	8.5	●	
Evolution 100	15	●	
Evolution 300	8.5	●	
Evolution 600	8.5	●	
Genesys Serie	8.5	●	
Helios Serie	15	●	
UV Serie	15	●	
UV1	15	●	
Varian			
Cary 50	15	●	

One glance at the chart will tell you which BRAND UV-cuvette micro is compatible with your photometer. The chart serves as a guideline. It has been carefully proven and it's up-to-date on what we know by now. Photometers not mentioned have not been tested yet.

We are continuously updating the table with information from manufacturers and feedback from our customers. If you are working with a photometer not mentioned, ask for free samples of our UV-cuvettes micro with no obligation (www.brand.de) to check compatibility with your instrument. A short feedback from you (to info@brand.de) about the results will help us to complete the chart.

To thank you for your effort, we will send you free of charge a pack of 100 BRAND micro cuvettes of your choice!

● = applicable

- = not applicable



BRAND