



## ***Basic Weighing Capabilities at an Economical Price***

The Compass CR series provides basic weighing functionality and competitive performance at an economical price. Features such as a slim, stackable design, and an enlarged weighing platform, make the CR series perfect for a variety of portable weighing applications.

### ***Standard Features:***

- **Basic Weighing at a Competitive Price**

The Compass CR series is designed for basic weighing at a competitive price point, and is ideal for use in kitchens, classrooms, offices and more.

- **Space-Saving Stackable Design for Convenient Storage**

The Compass CR series is specially designed for convenient storage. Its slim, space saving design allows users to stack up to 4 scales, saving valuable storage space when the scale is not in use.

- **Enlarged Weighing Platform Allows for Versatile Use**

With an enlarged weighing platform, the Compass CR series can adequately handle a wide range of sample sizes. You will only need one scale on your counter for basic weighing with the CR series.

# COMPASS™ CR SERIES Compact Scales

| Model                                | CR221  | CR621 | CR2200 | CR5200 |
|--------------------------------------|--|-------|--------|--------|
| Capacity (g)                         | 220  | 620   | 2200   | 5200   |
| Readability (g)                      | 0.1  |       | 1      |        |
| Span Calibration Mass (Not Included) | 200 g  | 500 g | 1 kg   | 2 kg   |
| Stabilization Time (s)               | 2  |       |        |        |
| Construction                         | ABS housing and pan  |       |        |        |
| Calibration                          | User-selectable external span / Digital with external weight |       |        |        |
| Tare Range                           | Full capacity by subtraction                                 |       |        |        |
| Weighing Units*                      | g, ozt, dwt, lb:oz   |       |        |        |
| Application Modes                    | Weighing   |       |        |        |
| Power Requirement                    | 3 AA batteries (included) or AC adapter (not included)       |       |        |        |
| Typical Battery Life                 | 300 hours  |       |        |        |
| Display Type                         | Liquid crystal display (LCD)                                 |       |        |        |
| Display Size                         | 15 mm digits   |       |        |        |
| Pan Size (W × D)                     | 125 × 132 mm   |       |        |        |
| Balance Dimensions (W × D × H)       | 140 × 205 × 41 mm  |       |        |        |
| Shipping Dimensions (W × D × H)      | 170 × 235 × 55 mm  |       |        |        |
| Net Weight                           | 0.4 kg   |       |        |        |
| Shipping Weight                      | 0.6 kg   |       |        |        |

\* Availability is dependent on model.

## Other Standard Features and Equipment

ABS housing and pan, slip-resistant feet, stability indicator, mechanical and software overload/underload protection, auto shut-off


## Compliance

- **Product Safety:** IEC/EN 61010-1
- **Electromagnetic Compatibility:** IEC/EN 61326-1 Class B, Basic environments; FCC Part 15 Class B; Canada ICES-003 Class B
- **Compliance Marks:** CE; RCM

## Accessories

|   |          |
|---|----------|
| Carrying Case .....                       | 30467763 |
| Bowl, 1050 mL .....                       | 80850075 |
| Adapter with plugs (EU-US-UK-AU-KR) ..... | 30467987 |
| Adapter, 12V 0.5A, fixed UK plug .....    | 30043166 |

## Dimensions



### OHAUS Europe GmbH

Im Langacher 44  
8606 Greifensee  
Switzerland

e-mail: [ssc@ohaus.com](mailto:ssc@ohaus.com)  
Tel: 0041 22 567 53 19  
e-mail: [tsc@ohaus.com](mailto:tsc@ohaus.com)  
Tel: 0041 22 567 53 20

[www.ohaus.com](http://www.ohaus.com)

*The management system governing the manufacture of this product is ISO 9001:2015 certified.*

Europe Literature #

