

## GW2145



### LABORATORY GLASSWARE WASHER

Basic line - glassware washer 45 cm

[Download PDF](#)

- washing up to 95°C - antibacterial thermal disinfection
- programs: 9 default + 6 free
- microprocessor programmable electronic control
- serial port RS232 for connection to external printer or PC
- active thermodynamic drying
- alkaline powder detergent dispenser
- peristaltic pump for acid liquid neutralizer
- possibility of installing up to 1 additional peristaltic pump
- washing chamber and internal door in AISI 316L stainless steel
- external panels made in AISI 304 stainless steel
- 2 independent washing levels through telescopic guides
- external dimensions: LxDxH: 450 x 620 x 850 mm
- internal working dimensions LxDxH: 370 x 480 x 570 mm
- incorporated water softener with salt sensor
- cold water and demineralized water connections
- electronic door lock system for user safety
- electrical connection: 1/N/PE 230 V - 50 Hz - 3,3 kW max
- 60 Hz version available
- Schuko plug 16 A - 250 V included
- upper and lower trolleys D-CS1 and D-CS2 included

Smeg S.p.A. Via Leonardo Da Vinci, 4 - 42016 Guastalla (RE) Italy - P. IVA 01555030350

Phone +39 0522 8211 - Fax +39 0522 821453 - [service.instruments@smeg.it](mailto:service.instruments@smeg.it)