

GW7010M

HIGH CAPACITY LABORATORY GLASSWARE WASHER (266 LT)
PROFESSIONAL LINE - LARGE DIMENSION GLASSWARE WASHERS
[Download PDF](#)

- ✓ washing up to 95°C - antibacterial thermal disinfection
- ✓ capacity up to 5 adjustable levels
- ✓ programs: 20 default + 20 custom
- ✓ multi-language TFT colour touch screen
- ✓ functioning protected by different levels password
- ✓ integrated built-in panel printer (optional)
- ✓ USB port for downloading washing cycles log
- ✓ ethernet port for LAN connection (optional)
- ✓ sliding down full glass door
- ✓ electronic door lock system for user safety
- ✓ high pressure hot air drying system
- ✓ high efficiency Ecoslim steam condenser (optional)
- ✓ peristaltic pump for liquid detergent
- ✓ peristaltic pump for liquid neutralizer
- ✓ possibility of installing up to 2 additional peristaltic pumps
- ✓ detergent cabinet with removable tray for storing up to 3 x 5 liters cans + 1 x 1 liter can
- ✓ washing chamber made of AISI 316L stainless steel
- ✓ external panels made of AISI 304 stainless steel
- ✓ external dimensions LxDxH: 780 x 800 x 1950 mm
- ✓ internal working dimensions LxDxH: 550 x 620 x 685 mm
- ✓ cold water, hot water and demineralized water connections
- ✓ floor drain connection
- ✓ drain pump (optional)
- ✓ electrical connection: 3/N/PE 400 V - 50 Hz - 13 kW max
- ✓ 60 Hz version available
- ✓ without racks

Smeg S.p.A. Via Leonardo Da Vinci, 4 - 42016 Guastalla (RE) Italy - P. IVA 01555030350
Phone +39 0522 8211 - fax +39 0522 821453 - service.instruments@smeg.it