

**GENO™ Type II Water Systems 10-50 litres per hour**

<b>Camlab Code</b>	<b>Part Number</b>	<b>Description</b>
1223777	GENO10	Geno™ - Type II Water System - 10L/hr
1227225	GENO20	Geno™ - Type II Water System - 20L/hr
1227226	GENO50	Geno™ - Type II Water System - 50L/hr
<b>Please select a reservoir for use with Geno and state position of tap to be centre or base</b>		
<b>Reservoirs</b>		
1213584	TANK30	30 Lit Storage Reservoir
1213585	TANK60	60 Lit Storage Reservoir
1213586	TANK100	100 Lit Storage Reservoir
1213587	TANK100PUMP	100 Lit Storage Reservoir with pump(8l/min)
<b>Accessories</b>		
1227415	AV004	Valve Kit for Washing Machine Feed (Single POU)
1227418	AV006	Wall Mounting Kit (Alto™/Duo™/Geno™/Puro™)
1227417	AV007	External Pre-filter Housing Kit (Excludes Filter)
1227416	AV008	Wall Mounting Kit (TANK30)
1227438	AV009	Wall Mounting Kit (TANK60)
1227439	AV012	UV Upgrade Kit (TANK30/TANK60)
1227442	AV013	UV Upgrade Kit (TANK100)
1213588	TC006	RO Membrane - 20L/hr
1227388	TC013	RO Membrane - 10L/hr
(1 x TC013 required for Geno 10, 2 x TC013 required for Geno 20, 3 x TC006 required for Geno 50)		
<b>Geno Consumables</b>		
1213590	TC001	Cartridge - Pre-treatment - Carbon
1227393	TC002	Endure Purification Pack 1
1226843	TC005	Filter - Air Vent (for any TANK)
1227389	TC008	Replacement Tank UV
1227385	TC012	Filter - Air Vent with CO <sub>2</sub> Trap (for any TANK)
1227403	TC038	Filter - 10" Spun Element (for AV007)
1227402	TC039	Filter - 10" Pleated Element (for AV007)
1227401	TC040	Filter - 10" Carbon Element (For AV007)
1227386	TC017	Annual Consumable Pack (Geno™) (1xTC005 + 2xTC001 + 2xTC002) If TC017 purchased with a Service Contract