

Labcold Flake Ice Makers



LITV-IQ50/IQ85C

Labcold Flake Ice Makers produce crushed ice ideal for filling packs and moulding to fit any item or object. The ice produced is both hygienic and clear so ideal for medical and laboratory applications.

Robustly constructed from stainless steel with well-insulated plastic collection bins and a lime scale reducing spray system.

Labcold Ice Makers require both a water and electricity supply and a floor level drain.



LITV-IQ135

Part Number	LITV-IQ50	LITV-IQ85C	LITV-IQ135
Daily Ice Production	43kg per day	100kg per day	163kg per day
Ice Storage	14kg	20kg	60kg
Style	Freestanding	Freestanding	Freestanding
Exterior Dimensions (HxWxD)	750 x 400 x 520mm	800 x 470 x 600mm	1360 x 520 x 550mm
Interior Construction	Plastic Bins	Plastic Bins	Plastic Bins
Exterior Construction	Stainless Steel	Stainless Steel	Stainless Steel
Automatic Production	✓	✓	✓
Door Hinging	Bottom inward folding	Bottom inward folding	Bottom inward folding
Door Style	Solid	Solid	Solid
Feet	4	4	4