

JB Nova - 5L, 12L, 18L and 26L

**NEW**

- PERFORMANCE › Set and Forget™ technology - fast heat-up, reliable temperature control
- INTUITIVE › Simple to use controls, clear bright display
- QUALITY › 3 year warranty, UK design and manufacture
- PRACTICAL › Display lock - disables front panel controls
- CONVENIENT › Drain tap, tray and lid included. Suitable for use with heat transfer beads
- VALUE › Advanced dry start and run dry protection - prevents costly service repairs







### Applications

- Education - school/college university teaching lab
- Healthcare - sample warming
- Industry - QA/QC testing

Ideal choice for schools and colleges

## Specifications

CE

	JBN5	JBN12	JBN18	JBN26
	 3 kg l: 215 mm w: 335 mm h: 200 mm	 6 kg l: 380 mm w: 360 mm h: 225 mm	 9 kg l: 590 mm w: 335 mm h: 275 mm	 9 kg l: 590 mm w: 335 mm h: 275 mm
Capacity	5L	12L	18L	26L
Temperature range	ambient +5 to 95°C			
Temperature display and setting resolution	0.5°C			
Stability (DIN 12876-3 @70 °C)	±0.5°C			
Working area (l x w) mm	131 x 281	281 x 306	485 x 281	481 x 278
User defined calibration points	1	1	1	1
Drain tap	-	•	•	•
Heater power 230V/120V (kW)	0.35 / 0.35	0.8 / 0.8	1.4 / 1.05	1.4 / 1.05
Supply voltage (V)	120 or 230 (50-60Hz)			

## Accessories (full range available on website)

### Flat lid with ring sets

Variable hole diameter to accommodate tall vessels whilst reducing evaporation



### Polypropylene spheres

Useful alternative to a lid, minimises evaporation and heat loss. Easy access to vessels in the bath; particularly useful for tall vessels



### Test tube racks

Stainless steel - range of sizes available for 0.5 to 1.5 ml microtubes or 10 to 30 mm Ø tubes



### Raised shelves

Stainless steel - shelf covers half area of bath. Reversible, allows two shelf depths



### Drain stopper

Recommended when using heat transfer beads, to prevent beads entering drain hole



## Alternative water bath ranges:




### SUB Aqua Pro - advanced range

- Temperature range amb. +5 to 99°C
- Stability ±0.2°C
- 8 bath sizes: 2L, 2L shallow, 5L, 12L, 18L, 26L, 5L & 12L dual bath
- 3 temperature presets
- User defined high temperature alarm and cutout

### JB Nova - basic

- Simple and practical, great value
- 3 bath sizes: 5L, 12L, 18L



  
www.camlab.co.uk  
**camlab**

Call us:  
+ 44 (0) 1954 233 110

Email:  
sales@camlab.co.uk