

## High-performance wiper solutions for the laboratory environment

### KIMTECH SCIENCE\* Wipers

#### **KIMTECH SCIENCE\* Delicate Task Wipers**

- Unique combination of low-lint tissue softness and multi-ply (2-ply) absorbency that's perfect for a multitude of delicate tasks.
- Available in two box quantity, 200 wipes and 100 wipes.

Designed specially to handle a variety of delicate wiping tasks in a laboratory such as surface/parts cleaning, glassware cleaning and instrument cleaning. KIMTECH SCIENCE\* Wipers are ideal to wipe up liquids and dust. They are supplied in handy pop-up boxes equipped with an anti-static dispensing to reduce lint and antistatic discharge.

#### **KIMTECH SCIENCE\* Precision Wipers**

- Light duty, single ply wiper for a variety of precision wiping tasks.
- Available in two sizes: 30.4cm x 30.4cm and 11.5cm x 21.5cm.



# High-performance wiper solutions for the laboratory environment

## KIMTECH SCIENCE\* Wipers

KIMBERLY-CLARK PROFESSIONAL\* has long been recognised for the development of technologically advanced products that consistently meet the quality of laboratories worldwide.

Count on KIMBERLY-CLARK PROFESSIONAL\* for everything you need to maintain your laboratory environment, from technical data to product quality testing and exceptional service.



### KIMTECH SCIENCE\* Delicate Task Wipers

Code	Delivery System	Colour	Material	Size WxL	Wipers/Unit	Case Content
7557	Pop-up	white	paper	20cm x 21cm	100	100 x 24 = 2400
7558	Pop-up	white	paper	20cm x 20,5cm	200	200 x 15 = 3000

### KIMTECH SCIENCE\* Precision Wipers

Code	Delivery System	Colour	Material	Size WxL	Wipers/Unit	Case Content
7551	Pop-up	white	paper	30.5cm x 30.5cm	196	196 x 15 = 2940
7552	Pop-up	white	paper	11.5cm x 21cm	280	280 x 30 = 8400

#### INFORMATION SERVICE

For technical enquiries please email [infofax@kcc.com](mailto:infofax@kcc.com)  
For sales enquiries please email [kimtech.support@kcc.com](mailto:kimtech.support@kcc.com)

**[www.kcprofessional.com](http://www.kcprofessional.com)**

For more information on the KIMTECH\* brand visit our website:

**[www.kimtech.eu](http://www.kimtech.eu)**