

PRODUCT INFORMATION

Model: DF3517L

Type: Labcold ilShin Ultralow Temperature Freezer (-35 to -45°C)

For laboratory and general use where storage at temperatures down to -45°C is required. This range offers low noise, a high level of construction, robustness and features coupled with great reliability and a comprehensive five year warranty just for good measure.



**Shown with
optional
temperature chart
recorder**

Capacity: 484 litres (16.9 cu.ft.)

Construction: Scratch resistant epoxy exterior; White enamel interior; (stainless steel optional)

Cabinet mounted on heavy duty castors

Inner compartments: 4; with insulated inner doors

Exterior dimensions (mm): 1980 high x 978 wide x 1010 deep

Features:

Microprocessor temperature control with digital display

Audible and visible hi/lo temperature alarms with remote contacts

Lockable handle

Exceptionally low noise: approx 50dBa

Low frost build-up due to door design and frame heaters

Two access ports (1 rear, 1 side)

CFC and HCFC free

Temperature set-point range: -35 to -45°C

Maximum ambient temperature: 32°C

5 year full parts and labour warranty

Order reference: DF 3517L