

PRODUCT INFORMATION

Model: DF9007

Type: Labcold ilShin Ultralow Temperature Freezer (-85°C)

For laboratory and general use where storage at temperatures down to -85°C is required. This range offers low noise, a high level of construction, robustness and features coupled with great reliability and a comprehensive five year warranty just for good measure.



Shown with
optional
temperature chart
recorder

Capacity: 191 litres (6.7 cu.ft.)

Construction: Scratch resistant epoxy exterior; White enamel interior; (stainless steel optional)

Cabinet mounted on heavy duty castors

Inner compartments: 1; with insulated inner lids

Exterior dimensions (mm): 1010 high x 1280 wide x 790 deep

Features:

Microprocessor temperature control with digital display

Audible and visible hi/lo temperature alarms with remote contacts

Lockable lid (with suitable padlock, not supplied)

Exceptionally low noise: approx 50dBa

Low frost build-up due to lid design and frame heaters

2 Access ports

CFC and HCFC free

Temperature set-point range: -50 to -85°C

Maximum ambient temperature: 32°C

5 year full parts and labour warranty

Order reference: DF 9007