

## PRODUCT INFORMATION

**Model: RAFR44042**

**Type: Labcold Advanced Laboratory Refrigerator**

For laboratory and general use where close temperature control and uniformity is required together with a high level of construction, robustness and features, and where there is no risk that volatile, flammable materials or components might be stored.



**Capacity: 1245 litres (gross)**

**Construction: White and blue enamel exterior, stainless steel interior.**  
Cabinet mounted on castors, lockable.

**Shelves Supplied: 8; Plastic coated wire mesh**

**Exterior Dimensions (mm): 1980 high x 1440 wide x 820 deep**

**Features:**

Microprocessor temperature control with digital display

Audible and visible hi/lo temperature alarms

Fan circulated interior

Door Lock

CFC and HCFC free

Approximate controllable range: +2 to +12°C

Accuracy: at 4°C ± 3°C, including defrost

Automatic defrost with condensate removal

Maximum ambient temperature: +32°C

**Order Reference: RAFR 44042**

**Optional Accessories: Remote Alarm Contacts: REMCONADV**

**32mm/50mm Access Ports: RLAP00\***

**Wall Mount Chart Recorder: RECW5050**