


PRODUCT INFORMATION

Model: RECW RANGE

Type: Labcold Electronic Chart Recorder

Labcold chart recorders provide permanent records of storage temperatures in a wide variety of units. Recorders have a thermocouple temperature sensor and are supplied for wall mounting (fitting not included).



Temperature Ranges:

RECW0051: 0°C to +50°C
RECW5050: -50°C to +50°C
RECW5001: -50°C to 0°C

Sensor Type:

Type K, Thermocouple

Power:

9V power supply. 230/240V, 50Hz power adapter included

Battery back up:

9V Battery

Pen Life:

Approx. 4-6 months

Exterior Dimensions (mm):

190 x 190 x 80

Chart Size:

153mm (6") diameter

Recording Period:

7 days per chart

Order Reference: RECW****