

## **PRODUCT INFORMATION**

**Model: RLBF1740 and RLBF1740/LK (Lockable)**

**Type: Labcold Sparkfree Laboratory Superfreezer (-40°C)**

For laboratory and general use where storage at below normal freezing temperatures down to -40°C and/or where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.



**Capacity: 50 Litres (gross)**

**Construction: White enamel exterior, moulded plastic interior**

**Shelves supplied: 1, wire mesh**

**Exterior dimensions (mm): 580 high x 550 wide x 620 deep**

### **Features:**

**Audible and visible high temperature alarm**

**Digital temperature controller with display**

**Sparkfree interior**

**Counter-top use**

**CFC and HCFC free**

**Temperature set-point range: -30 to -40°C**

**Operator controlled defrost**

**Remote alarm connections available as optional extra**

**Order reference: RLBF 1740**

**RLBF 1740/LK**