

PRODUCT INFORMATION

Model: RLCF1520

**Type: Labcold Sparkfree Laboratory Freezer (-20°C)
Top opening, Chest style**

For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future. Chest configuration provides improved temperature uniformity and stability with much reduced temperature rise following door openings.



Capacity: 447 litres (gross)

Construction: White enamel exterior, Aluminium interior

Baskets supplied: 2, wire mesh

Exterior dimensions (mm): 860 high x1300 wide x 730 deep (incl. handle and hinges)

Features:

Sparkfree interior

Door lock

Easy-grip door handle

CFC and HCFC free

Temperature set-point range: -18 to -22°C

Visual high temperature alarm

Operator controlled defrost

Order reference: RLCF 1520