

PRODUCT INFORMATION

Model: RLCH01501

Type: Labcold Refrigerated (cooled) Incubator (Glass door)

For maintaining samples at temperatures accurately controlled within the range 10 to 40°C. Temperatures can be cycled where appropriate accessories are fitted.



Capacity: 145 litres (gross)

Construction: White enamel exterior, moulded plastic interior, glass door.

Shelves supplied: 2; Plastic coated wire mesh

Exterior dimensions (mm): 860 x 600 x 650

Features:

Microprocessor temperature control with digital display

Temperature set-point range: 10 to 40°C

Temperature precision: better than +/- 0.5°C at 5°C

Fan circulated interior

CFC and HCFC free

Manual defrost with condensate drain

Maximum ambient temperature: 32°C

Order reference: RLCH 01501