

PRODUCT INFORMATION

Model: RLDF02042

Type: Labcold Pharmacy/Vaccine Refrigerator

For laboratory, surgery and general use for the storage of drugs and vaccines. Conforms in full with NEW UK Department of Health "Green Book" requirements.



Capacity: 60 litres

Construction: White enamel exterior, fully moulded plastic interior

Shelves supplied: 2 wire mesh

Exterior dimensions (mm): 572H x 485W x 535D

Features:

Compressor refrigeration (does not use ammonia)

Digital display of temperature with max/min record

Audible alarm with mute facility

Door lock

CFC, HCFC, ammonia free

Accurate temperature control: +2°C to +8°C

Order reference: RLDF 02042

Labcold Limited, The Loddon Centre, Wade Road, Basingstoke, RG24 8FL, UK

MARCH 2006