

PRODUCT INFORMATION

Model: RLDF1804

Type: Labcold Pharmacy/Vaccine Refrigerator

For laboratory, surgery and general use for the storage of drugs and vaccines. Conforms in full with New UK Department of Health "Green Book" requirements.



Capacity: 490 litres (gross)

Construction: White enamel exterior and interior

Shelves supplied: 6, plastic coated wire mesh

Exterior dimensions (mm): 1775H x 760W x 720D

Features:

Digital display of temperature with max/min record

Audible alarm with mute facility

Fan assisted air circulation

Door lock

CFC, HCFC and Ammonia free

Accurate temperature control: +2°C to +8°C

Automatic defrost

Order reference: RLDF 1804