

PRODUCT INFORMATION

Model: RLFR1804

Type: Labcold Fan Circulated Economy Refrigerator

For laboratory and general use where large capacity combined with value for money is required and there is no risk that volatile, flammable materials or components might be stored.



Capacity: 505 litres (gross)

Construction: White enamelled steel exterior and interior

Shelves supplied: 4; Plastic coated wire mesh

Exterior dimensions (mm): 1710 high x 765 wide x 725 deep

Features:

Fan circulated interior

Door lock

CFC and HCFC free

Temperature set-point range: +2 to +5°C

Analogue temperature display

Automatic defrost

Order reference: RLFR 1804