

## **PRODUCT INFORMATION**

**Model: RLHD14301**

**Type: Labcold Laboratory Superfreezer (-30°C)**

For laboratory and general use where storage at below normal freezing temperatures down to -30°C is required together with a high level of construction, robustness and features, and where there is no risk that volatile, flammable materials or components might be stored.



**Capacity: 385 litres (gross)**

**Construction: White enamel exterior, stainless steel interior. Cabinet mounted on castors.**

**Shelves supplied: 4; Plastic coated wire mesh**

**Exterior dimensions (mm): 1995 high x 600 wide x 700 deep**

**Features:**

**Microprocessor temperature control with digital display**

**Audible and visible hi/lo temperature alarms**

**Fan circulated interior**

**Door lock**

**CFC and HCFC free**

**Temperature set-point range: -20 to -30°C**

**Automatic defrost with condensate removal**

**Maximum ambient temperature: 32°C**

**Order reference: RLHD 14301**

N.B. Remote alarm contacts no longer supplied as standard but available at additional cost – order ref. REMCONADV