

PRODUCT INFORMATION

Model: RLHE0845

**Type: Labcold Sparkfree Laboratory Superfreezer (-45°C)
Top opening, Chest style**

For laboratory and general use where there is a need to store at lower than normal freezing temperatures and/or there is concern that volatile, flammable materials or components might be stored now, or at some point in the future. Chest configuration provides improved temperature uniformity and stability with much reduced temperature rise following door openings.



Capacity: 237 litres (gross)

Construction: White enamel exterior, Aluminium interior, Additional insulation

Baskets supplied: 1, wire mesh

Exterior dimensions (mm): 860 high x 1050 wide x 730 deep (incl. handle and hinges)

Features:

Sparkfree interior

Door lock

Easy-grip door handle

CFC and HCFC free

Temperature set-point range: -20 to -45°C

Audible high temperature

Operator controlled defrost

Remote alarm connections available as optional extra

Order reference: RLHE 0845