

## **PRODUCT INFORMATION**

**Model: RLHE1145**

**Type: Labcold Sparkfree Laboratory Superfreezer (-45°C)  
Top opening, Chest style**

For laboratory and general use where there is a need to store at lower than normal freezing temperatures and/or there is concern that volatile, flammable materials or components might be stored now, or at some point in the future. Chest configuration provides improved temperature uniformity and stability with much reduced temperature rise following door openings.



**Capacity: 314 litres (gross)**

**Construction: White enamel exterior, Aluminium interior, Additional insulation**

**Baskets supplied: 2, wire mesh**

**Exterior dimensions (mm): 860 high x 1305 wide x 730 deep (incl. handle and hinges)**

### **Features:**

**Digital display of temperature**

**Sparkfree interior**

**Door lock**

**Easy-grip door handle**

**CFC and HCFC free**

**Temperature set-point range: -20 to -45°C**

**Audible high temperature**

**Operator controlled defrost**

**Remote alarm connections available as optional extra**

**Order reference: RLHE 1145**