

## **PRODUCT INFORMATION**

**Model: RLPF05042**

**Type: Labcold Fan-Assisted Laboratory Refrigerator**

For laboratory and general use where improved temperature uniformity is required and there is no risk that volatile, flammable materials or components might be stored.



**Capacity: 160 litres (gross)**

**Construction: White enamel exterior, fully moulded plastic interior**

**Shelves supplied: 3, toughened glass**

**Exterior dimensions (mm): 850 high x 550 wide x 620 deep**

**Features:**

**Fan assisted air circulation**

**Removable door pockets**

**Door lock**

**Easy-grip door handle**

**CFC and HCFC free**

**Temperature set-point range: +2 to +8°C**

**Low energy operation**

**Automatic defrost**

**Order reference: RLPF 05042**