

PRODUCT INFORMATION

Model: RLPR05042

Type: Labcold Sparkfree Laboratory Refrigerator

For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.



Capacity: 160 litres (gross)

Construction: White enamel exterior, fully moulded plastic interior

Shelves supplied: 3, toughened glass

Exterior dimensions (mm): 850 high x 550 wide x 620 deep

Features:

Sparkfree interior

Removable door pockets

Door lock

Easy-grip door handle

CFC and HCFC free

Temperature set-point range: +2 to +8°C

Low energy operation

Automatic defrost

Order reference: RLPR 05042