

PRODUCT INFORMATION

Model: RLPR06042 and RLPR06042/LK (Lockable)

Type: Labcold Sparkfree Laboratory Refrigerator

For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future. This unit features a removable top panel to facilitate under counter installation.



Capacity: 160 Litres (gross)

Construction: White enamel exterior, fully moulded plastic interior

Shelves supplied: 3, glass

Exterior dimensions (mm): 850 high (820 without top) x 550 wide x 620 deep

Features:

Sparkfree interior

Removable top panel to reduce height

CFC and HCFC free

Temperature set-point range: +2 to +8°C

Automatic defrost

Order reference: RLPR 06042

RLPR 06042/LK