

## **PRODUCT INFORMATION**

**Model: RLSP1804**

**Type: Labcold Sparkfree Laboratory Refrigerator  
with Fan Air Circulation**

For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future. Fan circulation of chamber air provides much improved temperature uniformity.



**Capacity: 490 litres (gross)**

**Construction: White enamel exterior and interior**

**Shelves supplied: 6; Plastic coated wire mesh**

**Exterior dimensions (mm): 1775 high x 760 wide x 720 deep**

**Features:**

**Sparkfree interior**

**Internal air circulation powered by external motor**

**Door lock**

**CFC and HCFC free**

**Temperature set-point range: 0 to +8°C**

**Digital temperature display with audible alarm**

**Automatic defrost**

**Order reference: RLSP 1804**