

PRODUCT INFORMATION

Model: RPBB0006

Type: Labcold Portable Digital Transporter

For the temperature controlled transport, and short-term holding of samples and specimens at a fixed temperature of $4^{\circ}\text{C} \pm 2$.



Capacity: 12 litres

Weight: 6.8kg

Construction: Fully moulded plastic exterior, aluminium interior

Exterior dimensions (mm): 290H x 415W x 285D

Interior dimensions (mm): 200H x 335W x 160D

230 volt A.C. and 12 volt D.C. operation (cables supplied)

Features:

Digital temperature display

Solid state refrigeration system

CFC, HCFC, and ammonia free

Temperature: $4^{\circ}\text{C} \pm 2$ (in 30°C ambient)

Removable wire-mesh basket

Carry strap and internal lid supplied

Order reference: RPBB 0006