

## PRODUCT INFORMATION

**Model: RPDF0012**

**Type: Labcold Portable Vaccine and Sample Carrier**

For laboratory, surgery and general use for the transport of drugs, vaccines and other temperature sensitive materials. Not suitable for the transport of blood and blood components.



**Capacity: 12 litres**

**Weight: 6.2kg**

**Construction: Fully moulded plastic exterior, aluminium interior**

**Exterior dimensions (mm): 290H x 415W x 285D**

**Interior dimensions (mm): 200H x 335W x 160D**

**230 volt A.C. and 12 volt D.C. operation (cables supplied)**

**Features:**

**Solid state refrigeration system**

**CFC, HCFC, ammonia free**

**Temperature set-point range: -5 to +60°C (in 20°C ambient)**

**Carry strap and internal lid supplied**

**Order reference: RPDF 0012**