

PRODUCT INFORMATION

Model: RPFR14041

Type: Labcold Advanced Pharmacy Refrigerator

For pharmacy, laboratory and general use where close temperature control and uniformity is required together with a high level of construction, robustness and features, and where there is no risk that volatile, flammable materials or components might be stored.



Capacity: 385 litres (gross)

Construction: White enamel exterior, stainless steel interior. Cabinet mounted on castors.

Shelves supplied: 4; Plastic coated wire mesh

Exterior dimensions (mm): 1995 high x 600 wide x 700 deep

Features:

Microprocessor temperature control with digital display

Audible and visible hi/lo temperature alarms

Fan circulated interior

Door lock

CFC and HCFC free

Temperature set-point range: +2 to +12°C

Automatic defrost with condensate removal

Maximum ambient temperature: 43°C

Order reference: RPFR 14041

N.B: Remote alarm contacts no longer supplied as standard but available at additional cost – order ref. REMCONADV