

## **PRODUCT INFORMATION**

**Model: RTBB0006**

**Type: Labcold Insulated Transport Box**

A professional, low cost, reusable insulated box for the transport of samples, specimens and biological materials between hospital departments within a 30 minute period.



**Capacity: 6 Litres**

**Construction: Strong, wipeable, red plastic with white insulated interior. Label provided.**

**Exterior Dimensions (mm): 180 high, 336 wide, 165 deep.**

### **Features:**

- Extremely lightweight – 400g empty**
- Manufactured from wipeable plastics**
- Stick-on label with space for writing**
- Easy-grip carrying handle**
- Supplied flat-packed in packs of 10 for ease of storage**

**NOTE: This product is not for holding blood for extended periods of time.  
Recommended maximum hold time of half an hour.**

**Order Reference: RTBB0006**