

MSE

Medium Capacity Bench Top Centrifuges

- Simple - Easy operation
- Fast - Up to 18,000rpm
- Powerful - Brushless induction motor
- Flexible - 10 acceleration and brake rates
- Practical - Low height for easy loading
- Safe - Strong lid inter-lock
- Functional - Up to 48 blood tubes per run
- Dependable - Quality British build

The **Harrier 15/80** and refrigerated **Harrier 18/80R** from MSE



Quality Laboratory Equipment

MSE UK Ltd

Distributed by Camlab Ltd

Camlab Ltd, 24 Norman Way Industrial Estate,
Over, Cambridge CB24 5WE, UK

Tel: +44 1954 233 110
 Fax: +44 1954 233 101
www.camlab.co.uk
sales@camlab.co.uk



Harrier 15/80 and 18/80R Technical Data

	Harrier 15/80	Harrier 18/80R
Height (mm)	355	355
Width (mm)	470	670
Depth (mm)	580	580
Weight (kg)	41	65
Maximum speed: swing out rotor	4,500 rpm	6,000 rpm
Maximum RCF: swing out rotor	3,637 g	6,466 g
Maximum speed: angle rotor	15,000 rpm	18,000 rpm
Maximum RCF: angle rotor	20,125 g	28,980 g
Maximum capacity	4 x 200 ml	4 x 200 ml
Number of acceleration/brake rates	10	10
Timer	0-59 mins 59 secs, 0-99hrs 59 minutes and hold	
Speed display	Digital LED \pm 10 rpm	Digital LED \pm 10 rpm
Selectable temperature range	N/A	-9°C to + 40°C
Catalogue number: 220/230V 50Hz	MSB080.CX1.5	MSB020.CR1.K
Catalogue number: 120V 60Hz	MSB080.CX1.1	MSB080.CR1.H

43124-134 Swing Out Rotor with 43551-1215 Buckets

Tube or Sample Type	Adapter Part Number	Maximum Rotor Capacity
2-5ml blood tube	34159-306A	48
5-7ml blood tube	34159-3055	36
10ml flat bottomed tube	34159-3036	16
10-15ml blood tubes	34159-3040	28
15ml flat bottomed tube	34159-308	28
15ml Falcon tube	34159-310	12
50ml Falcon tube	34159-309	4
50ml universal tube	34159-3017	4
100ml universal tube	34159-3021	4
200ml Polypropylene bottle	Bottle p/n 43551-103A	4

46124-134 Swing Out System for 50ml Falcon Tubes (15/80)

Tube or Sample size	Trunnion	Bucket	Adapter	Rotor Capacity	Maximum G-force
50ml Falcon tube	34136-110	43156-603	34142-110	8	3,637

Angle Rotor System

Tube or Sample Type	Angle Rotor	Bucket	Adapter	Rotor Capacity	Maximum G-force
2.0ml	43117-612	-	-	24	28,980
15ml	43117-607	43156-601	cushion	12	5,290
50ml universal	43117-608	46156-602	cushion	6	5,574

Sealed Angle Rotor System

Tube or Sample size	Angle Rotor	Bucket	Adapter	Rotor Capacity	Maximum G-force
50ml Falcon tube	43117-609	43156-603	34142-110	4	4,812

Microplate Rotor System

Sample type	Rotor part number	Carrier part number	Rotor capacity	Maximum G-force
Microplates	43121-117	included	4 normal / 2 deep well	982

Centrifuges sold subject to our standard terms and conditions of sale, available on request. All information contained in this brochure is correct at time of printing, and may be subject to change at any time. The MSE logo is a registered trademark of MSE (UK) Ltd, registered address, Worsley Bridge Road, London, SE26 5AZ.