


**MSE**

## Small Capacity Bench Top Centrifuges

- Simple - Easy operation
- Fast - Up to 6,000 rpm / 3,340 g
- Quiet - Low noise output
- Flexible - carries 2 - 15ml blood tubes
- Practical - Just 5 kilogrammes
- Safe - Lid locking system
- Functional - Up to 6 blood tubes per run
- Dependable - Reliable

The **Kestrel** from MSE



Quality Laboratory Equipment

### MSE UK Ltd

**Distributed by Camlab Ltd**

Camlab Ltd, 24 Norman Way Industrial Estate,  
Over, Cambridge CB24 5WE, UK


Tel: +44 1954 233 110  
Fax: +44 1954 233 101  
www.camlab.co.uk  
sales@camlab.co.uk



## Kestrel Technical Data

Height (mm)	180
Width (mm)	220
Depth (mm)	270
Weight (kg)	5
Maximum speed	6,000 rpm
Maximum RCF	3,340 g
Maximum capacity	6 x 15ml
Timer	1-60 minutes or continuous
Time Increments	15 seconds
Speed display	Digital LED $\pm$ 10 rpm
Selectable temperature range	N/A
Catalogue number: 220/230V 50Hz	MSB015.CX1.5
Catalogue number: 120V 60Hz	MSB015.CX1.1

Tube or Sample Type	Extra Cushion Required?	Maximum Rotor Capacity
8-15ml blood tube (110mm L)	No	6
15ml BD Falcon© (conical) tube	No	3
10ml flat bottomed tube	No	6
5-10ml Sarstedt Monovette©	No	6
5-7ml blood tube (110mm L)	No	6
2-4ml blood tube (80mm L)	Yes	6



The initials 'BD' and 'Falcon' are registered trademarks of Beckton Dickenson. The words 'Sarstedt' and 'monovette' are registered trademarks of Sarstedt. Please verify maximum permissible RCF with tube manufacturers prior to centrifugation.

Centrifuges sold subject to our standard terms and conditions of sale, available on request. All information contained in this brochure is correct at time of printing, and may be subject to change at any time. The MSE logo is a registered trademark of MSE (UK) Ltd, registered address, Worsley Bridge Road, London, SE26 5AZ.