

# Multi-Parameter Meters at a glance

Enter this number into the search box at [www.camlab.co.uk](http://www.camlab.co.uk) to see more information, prices and live stock data

	Our Barcode	Product	Temperature	pH	Conductivity	Total Dissolved Solids (TDS)	Salinity	Oxidation Reduction Potential (ORP) (aka Redox)	Dissolved Oxygen (DO)
Handheld / Portable	1137070	EC/TDS/ Salt Tracer Pocketester	✓		✓	✓	✓		
	1187525	Sension+ MM150 multi pH Redox Conductivity TDS and Temp meter. Instrument Only	✓	✓	✓	✓		✓	
	1134146	HQ30D Flexi Meter With Gel filled pH Probe 1M and LDO Probe 1M		✓					✓
	1160642	HQ40D Multi MeterField Kit 1M Gel Filled pH and 1M Conductivity Probe		✓	✓				
	1134154	HQ40D Multi Meter 1M Gel Filled pH probe and 1M LDO probe		✓					✓
	1135990	Model 430 pH/Conductivity meter supplied in carry case with probe		✓	✓				
	1130373	MyronL Ultrameter 6PII Bundle		✓	✓			✓	
	1137074	pH /Conductivity Tracer Pocketester	✓	✓	✓	✓	✓		
	1134454	Portable Blue pH & conductivity meter with HI-1285 probe		✓	✓	✓			
	1134031	Sension 156 pH/DO/Conductivity Meter With Standard Electrodes		✓	✓				✓
	1160391	Techpro II Meter TDS, conductivity, temperature	✓		✓	✓			
	1160582	Techpro II Meter TPH1 pH, Conductivity, TDS and Temperature	✓	✓	✓	✓			
	1139362	Ultrameter 4PII Waterproof Conductivity TDS Resistivity Meter			✓	✓			
	1173075	Ultrameter 6PIIc Conductivity TDS Resistivity pH ORP Temperature	✓	✓	✓	✓		✓	
1177478	Waterproof pH/Conductivity/TDS Meter		✓	✓	✓				
Benchtop	1135984	Model 3540 pH/conductivity meter with dual graphics display		✓	✓				
	1188376	HQ430d, Digital Multimeter, benchtop, One-channel pH/ORP/Conductivity/DO(LDO) Meter without probes		✓	✓			✓	✓
	1191678	HQ440d, Digital Multimeter, benchtop, two-channel pH/ORP/Conductivity/DO(LDO) Meter without probes		✓	✓			✓	✓