




COMPARE PRODUCTS

[Home](#) / [Products](#) / [Compare Products](#)

[print this page](#)

[Back to the Product Selector](#)

COMPARE PRODUCTS

			
	Navigator® NV	Navigator® XL	Navigator® XT
CAPACITY	210 g - 5100 g	510 g - 20000 g	1600 g - 16000 g
READABILITY	.01 g, .1 g, .05 ct, .5 ct, 5 ct, 1 g, .002 lb, .2 g, .5 g, .0005 lb, 1 ct, 2.5 ct	.1 g, .5 ct, 5 ct, 1 g, .5 g, 2.5 ct	.1 g, .5 ct, 5 ct, 1 g, .002 lb, .2 g, .5 g, .0005 lb, 1 ct, 2.5 ct, 2 g, .0002 lb, .001 lb, 10 ct
PLATFORM SIZE	120mm(4.7in), 190x138mm(7.5x5.4in)	194x203mm	230x174mm
UNITS	Gram, Kilogram, Pound, Ounce (decimal or fractional), Pound:Ounce (decimal or fractional oz), Carat, Pennyweight, Troy Ounce, Grain, Newton, Hong Kong Tael, Singapore Tael, Taiwan Tael, Tical, Tola	Gram, Kilogram, Pound, Ounce (decimal or fractional), Pound:Ounce (decimal or fractional oz), Carat, Pennyweight, Troy Ounce, Grain, Newton, Hong Kong Tael, Singapore Tael, Taiwan Tael, Tical, Tola	Gram, Kilogram, Pound, Ounce (decimal or fractional), Pound:Ounce (decimal or fractional oz), Carat, Pennyweight, Troy Ounce, Grain, Newton, Hong Kong Tael, Singapore Tael, Taiwan Tael, Tical, Tola
APPLICATION MODES	Percent Weighing, CheckWeighing, Parts Counting, Accumulation/Totalization, Display Hold	Percent Weighing, CheckWeighing, Parts Counting, Accumulation/Totalization, Display Hold	Percent Weighing, CheckWeighing, Parts Counting, Accumulation/Totalization, Display Hold
DISPLAY	6 full segment LCD white auto backlight or Red LED, both with multiple icons and indicators	6 full segment LCD white auto backlight or Red LED, both with multiple icons and indicators	6 full segment LCD white auto backlight or Red LED, both with multiple icons and indicators
COMMUNICATION	Optional RS232 or USB or Ethernet Interface Kits with Integrated Cables	Optional RS232 or USB or Ethernet Interface Kits with Integrated Cables	Optional RS232 or USB or Ethernet Interface Kits with Integrated Cables
CONSTRUCTION	ABS plastic housing, 304 stainless steel pan	ABS plastic housing, 304 stainless steel pan	ABS plastic housing, 304 stainless steel pan
OPERATION	AC Adapter (included) or 4AA Batteries (not included)	AC Adapter or Internal Rechargeable Battery (not included) or 4 C Batteries (not included)	AC Adapter or Internal Rechargeable Battery (not included) or 4 C Batteries (not included)
DESIGN FEATURES	Up front illuminated level indicator, removable stainless steel weighing platform, menu lockswitch, integral security bracket, integral shipping lock, stability indicator, overload and underload indicators, low battery indicator, auto shut-off	Up front illuminated level indicator, removable stainless steel weighing platform, menu lockswitch, integral security bracket, integral shipping lock, stability indicator, overload and underload indicators, low battery indicator, auto shut-off	Up front illuminated level indicator, removable stainless steel weighing platform, menu lockswitch, integral security bracket, integral shipping lock, stability indicator, overload and underload indicators, low battery indicator, auto shut-off