




Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5816, FC5816R

| | | | | | | |
|---|--|----------------|----------------|--|----------------|-------------------|
| Rotor Photo  | Camlab product Code | 1207248 | Description | Rotor, Angle, 24x1.5/2.0ml, ID, Sealable | | |
| | Compatible with: | FC5714 | FC5718 | FC5718R | FC5816 | FC5816R |
| | Camlab product Code | 1206875 | 1206876 | 1206877 | 1206878 | 1206879 |
| | max. speed | 14000 rpm | 14000 rpm | 15000 rpm | 15000 rpm | 16000 rpm |
| | max. RCF | 18624 xg | 18624 xg | 21379 xg | 21379 xg | 24325 xg |
| Accessories | | | | | | |
| Camlab product Code | Description | | Vol./Tube | Diameter | Tube type | No.of tubes/rotor |
| No Adapter | Rotor, Angle, 24x1.5/2.0ml, ID, Sealable | | 1.5/2.0ml | 11mm | | 24 |
| 1207351 | Adapter, 1x0.2/0.4ml D6mm, 6/pk | | 0.2/0.4ml | 6mm | | 24 |
| 1207352 | Adapter, 1x0.5ml D8mm, 6/pk | | 0.5ml | 8mm | | 24 |


Rotors and Accessories for Ohaus Frontier FC 5714, FC5718, FC5718R

| | | | | | | |
|---|--------------------------------|----------------|----------------|-------------------------------|-----------|--------------------|
|  | Camlab product Code | 1207268 | Description | Rotor, Swing out, 4x100ml, ID | | |
| | Compatible with: | FC5714 | FC5718 | FC5718R | | |
| | Camlab product Code | 1206875 | 1206876 | 1206877 | | |
| | max. speed | 4500 rpm | 4500 rpm | 4500 rpm | | |
| | max. RCF | 3350 xg | 3350 xg | 3350 xg | | |
| Accessories | | | | | | |
| Camlab product Code | Description | | Vol./Tube | Diameter | Tube type | No. of tubes/rotor |
| 1207289 | Rack, 1x100ml D46mm RB, 2/pk | | 100ml | 46mm | RB | 4 |
| 1207290 | Rack, 1x85ml D38mm RB, 2/pk | | 85ml | 38mm | RB | 4 |
| 1207291 | Rack, 1x50ml D34mm RB, 2/pk | | 50ml | 34mm | RB | 4 |
| 1207292 | Rack, 1x50ml D29mm RB, 2/pk | | 50ml | 29mm | RB | 4 |
| 1207293 | Rack, 1x50ml D29mm FA, 2/pk | | 50ml | 29mm | FA | 4 |
| 1207294 | Rack, 1x50ml D29mm Rim, 2/pk | | 50ml | 29mm | Rim | 4 |
| 1207295 | Rack, 1x30ml D25mm Rim, 2/pk | | 30ml | 25mm | Rim | 4 |
| 1207296 | Rack, 1x25ml D25mm RB, 2/pk | | 25ml | 25mm | RB | 4 |
| 1207297 | Rack, 2x15ml D17mm FA, 2/pk | | 15ml | 17mm | FA | 8 |
| 1207298 | Rack, 1x15ml D17mm Rim, 2/pk | | 15ml | 17mm | Rim | 4 |
| 1207299 | Rack, 7x15ml D17mm RB, 2/pk | | 15ml | 17mm | RB | 28 |
| 1207300 | Rack, 9x5-7ml D13mm RB, 2/pk | | 5-7ml | 13mm | RB | 36 |
| 1207301 | Rack, 10x1.5ml D11mm, 2/pk | | 1.5ml | 11mm | | 40 |
| 1207302 | Rack, 9x4.5ml D12mm FB, 2/pk | | 4.5ml | 12mm | FB | 36 |
| 1207303 | Rack, 7x5-7ml D13mm RB, 2/pk | | 5-7ml | 13mm | RB | 28 |
| 1207304 | Rack, 7x5-7ml D13mm Vac, 2/pk | | 5-7ml | 13mm | Vac | 28 |
| 1207305 | Rack, 5x4-10ml D16mm Sar, 2/pk | | 4-10ml | 16mm | Sar | 20 |


Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R

| | | | | | | |
|---|---|----------------|----------------|---|-----------|--------------------|
| Rotor Photo  | Camlab product Code | 1207269 | Description | Rotor, Swing out, 4x100ml, ID, Sealable | | |
| | Compatible with: | FC5714 | FC5718 | FC5718R | | |
| | Camlab product Code | 1206875 | 1206876 | 1207877 | | |
| | max. speed | 4000 rpm | 5000 rpm | 5000 rpm | | |
| | max. RCF | 2486 xg | 3885 xg | 3885 xg | | |
| Accessories | | | | | | |
| Camlab product Code | Description | | Vol./Tube | Diameter | Tube type | No. of tubes/rotor |
| 1207283 | Bucket, 100ml w/o Cap, 2/pk | | | | | |
| 1207284 | Bucket, 100ml w/ Cap, Sealable, 2/pk | | | | | |
| 1207375 | Screw Cap, 100ml Bucket, Sealable, 2/pk | | | | | |
| 1207376 | O-ring, 100ml Bucket, Sealable, 10/pk | | | | | |
| 1207306 | Rack, 1x100ml D40mm RB, 2/pk | | 100ml | 40mm | RB | 4 |
| 1207307 | Rack, 1x85ml D38mm RB, 2/pk, v2 | | 85ml | 38mm | RB | 4 |
| 1207308 | Rack, 1x50ml D34mm RB, 2/pk, v2 | | 50ml | 34mm | RB | 4 |
| 1207309 | Rack, 1x50ml D29mm RB, 2/pk, v2 | | 50ml | 29mm | RB | 4 |
| 1207310 | Rack, 1x50ml D29mm FA, 2/pk, v2 | | 50ml | 29mm | FA | 4 |
| 1207311 | Rack, 1x50ml D29mm Rim, 2/pk, v2 | | 50ml | 29mm | Rim | 4 |
| 1207312 | Rack, 1x30ml D25mm Rim, 2/pk, v2 | | 30ml | 25mm | Rim | 4 |
| 1207313 | Rack, 1x25ml D25mm RB, 2/pk, v2 | | 25ml | 25mm | RB | 4 |
| 1207314 | Rack, 1x15ml D17mm FA, 2/pk | | 15ml | 17mm | FA | 4 |
| 1207315 | Rack, 2x15ml D17mm FA, 2/pk, v2 | | 15ml | 17mm | FA | 8 |
| 1207316 | Rack, 2x15ml D17mm Rim, 2/pk | | 15ml | 17mm | Rim | 8 |
| 1207317 | Rack, 4x15ml D17mm RB, 2/pk | | 15ml | 17mm | RB | 16 |
| 1207318 | Rack, 7x5-7ml D13mm RB, 2/pk, v2 | | 5-7ml | 13mm | RB | 28 |
| 1207319 | Rack, 5x1.5ml D11mm, 2/pk | | 1.5ml | 11mm | | 20 |
| 1207320 | Rack, 3x10ml D17mm Sar, 2/pk | | 10ml | 17mm | Sar | 12 |
| 1207321 | Rack, 4x5-7ml D13mm RB, 2/pk | | 5-7ml | 13mm | RB | 16 |
| 1207322 | Rack, 4x10ml D17mm Vac, 2/pk | | 10ml | 17mm | Vac | 16 |
| 1207285 | Bucket, 2x50ml D29mm FA, w/ Rack, 2/pk | | 50ml | 29mm | FA | 8 |
| 1207286 | Bucket, 3x15ml D17mm FA, w/ Rack, 2/pk | | 15ml | 17mm | FA | 12 |
| 1207287 | Cyto, Rectangular Bucket, 2/pk | | | | | |
| 1207377 | Cyto, Slide Carrier, 1/pk | | | | | |
| 1207378 | Cyto, Container, 2x1.5ml Cyl., 1/pk | | | | | |
| 1207379 | Cyto, Container, 2x3ml Conical, 1/pk | | | | | |
| 1207380 | Cyto, Container Cover, 1/pk | | | | | |
| 1207381 | Cyto, Filter Card, Cyl., 100/pk | | | | | |
| 1207382 | Cyto, Filter Card, Conical, 100/pk | | | | | |


Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R, FC5816, FC5816R

| | | | | | | |
|---|-----------------------------------|----------------|----------------|------------------------------|----------------|--------------------|
| Rotor Photo  | Camlab product Code | 1207270 | Description | Rotor, Swing out, 2x3MTP, ID | | |
| | Compatible with: | FC5714 | FC5718 | FC5718R | FC5816 | FC5816R |
| | Camlab Product Code | 1206875 | 1206876 | 1207877 | 1206878 | 1206879 |
| | max. speed | 4500 rpm | 4500 rpm | 4500 rpm | 4500 rpm | 4500 rpm |
| | max. RCF | 2716 xg | 2716 xg | 2716 xg | 2716 xg | 2716 xg |
| Accessories | | | | | | |
| Camlab product Code | Description | | Vol./Tube | Diameter | Tube type | No. of tubes/rotor |
| No Adapter | Rotor, Swing out, 2x3MTP, ID | | MTP | | | 6 |
| 1207288 | Bucket, 3xMTP, 2/pk | | MTP | | | 6 |
| 1207356 | Adapter, 48x1.5/2.0ml D11mm, 2/pk | | 1.5/2.0ml | 11mm | | 96 |


Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R

| | | | | | | |
|---|------------------------------------|----------------|----------------|--------------------------------|-----------|--------------------|
| Rotor Photo  | Camlab product Code | 1207254 | Description | Rotor, Angle, 6x50ml RB/FA, ID | | |
| | Compatible with: | FC5714 | FC5718 | FC5718R | | |
| | Camlab product Code | 1206875 | 1206876 | 1206877 | | |
| | max. speed | 6000 rpm | 6000 rpm | 6000 rpm | | |
| | max. RCF | 4427 xg | 4427 xg | 4427 xg | | |
| Accessories | | | | | | |
| Camlab product Code | Description | | Vol./Tube | Diameter | Tube type | No. of tubes/rotor |
| No Adapter | Rotor, Angle, 6x50ml RB/FA, ID | | 50ml | 29mm | RB/FA | 6 |
| 1207353 | Adapter, 1x1.5/2.0ml D11mm, 6/pk | | 1.5/2.0ml | 11mm | | 6 |
| 1207345 | Adapter, 1x7ml D13.5mm RB, 2/pk | | 7ml | 13.5mm | RB | 6 |
| 1207346 | Adapter, 1x5ml D13.5mm, 2/pk | | 5ml | 13.5mm | | 6 |
| 1207347 | Adapter, 1x30ml D26mm, 2/pk | | 30ml | 26mm | | 6 |
| 1207348 | Adapter, 1x16ml D18mm, 2/pk | | 16ml | 18mm | | 6 |
| 1207349 | Adapter, 1x15ml D17mm, 2/pk | | 15ml | 17mm | | 6 |
| 1207350 | Adapter, 1x15ml D17mm RB, 2/pk, v2 | | 15ml | 17mm | RB | 6 |


Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R, FC5816, FC5816R

| | | | | | | |
|---|----------------------------------|----------------|----------------|---------------------------------|----------------|-------------------|
| Rotor Photo  | Camlab product Code | 1207256 | Description | Rotor, Angle, 30x15ml RB/FA, ID | | |
| | Compatible with: | FC5714 | FC5718 | FC5718R | FC5816 | FC5816R |
| | Camlab Product Code | 1206875 | 1206876 | 1206877 | 1206878 | 1206879 |
| | max. speed | 4500 rpm | 4500 rpm | 4500 rpm | 4500 rpm | 4500 rpm |
| | max. RCF | 2830 xg | 2830 xg | 2830 xg | 2830 xg | 2830 xg |
| Accessories | | | | | | |
| Camlab product Code | Description | | Vol./Tube | Diameter | Tube type | No.of tubes/rotor |
| No Adapter | Rotor, Angle, 30x15ml RB/FA, ID | | 15ml | 17mm | RB/FA | 30 |
| 1207353 | Adapter, 1x1.5/2.0ml D11mm, 6/pk | | 1.5/2.0ml | 11mm | | 30 |
| 1207345 | Adapter, 1x7ml D13.5mm RB, 2/pk | | 7ml | 13.5mm | RB | 30 |
| 1207346 | Adapter, 1x5ml D13.5mm, 2/pk | | 5ml | 13.5mm | | 30 |


Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R

| | | | | | | |
|---|----------------------------------|----------------|----------------|------------------------------------|-----------|--------------------|
| Rotor Photo  | Camlab product Code | 1207258 | Description | Rotor, Angle, 12x15ml RB/FA, ID | | |
| | Compatible with: | FC5714 | FC5718 | FC5718R | | |
| | Camlab Product Code | 1206875 | 1206876 | 1206877 | | |
| | max. speed | 6000 rpm | 6000 rpm | 6000 rpm | | |
| | max. RCF | 4427 xg | 4427 xg | 4427 xg | | |
| Accessories | | | | | | |
| Camlab product Code | Description | | Vol./Tube | Diameter | Tube type | No. of tubes/rotor |
| No Adapter | Rotor, Angle, 12x15ml RB/FA, ID | | 15ml | 17mm | RB/FA | 12 |
| 1207353 | Adapter, 1x1.5/2.0ml D11mm, 6/pk | | 1.5/2.0ml | 11mm | | 12 |
| 1207345 | Adapter, 1x7ml D13.5mm RB, 2/pk | | 7ml | 13.5mm | RB | 12 |
| 1207346 | Adapter, 1x5ml D13.5mm, 2/pk | | 5ml | 13.5mm | | 12 |

Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R

| Rotor Photo | Camlab product Code | 1207260 | Description | Rotor, Angle, 30x1.5/2.0ml, ID, Sealable | | |
|---|--|----------------|----------------|--|-----------|--------------------|
|  | Compatible with: | FC5714 | FC5718 | FC5718R | | |
| | Camlab product Code | 1206875 | 1206876 | 1206877 | | |
| | max. speed | 12000 rpm | 13000 rpm | 14000 rpm | | |
| | max. RCF | 15131 xg | 17758 xg | 20595 xg | | |
| Accessories | | | | | | |
| Camlab product Code | Description | | Vol./Tube | Diameter | Tube type | No. of tubes/rotor |
| No Adapter | Rotor, Angle, 30x1.5/2.0ml, ID, Sealable | | 1.5/2.0ml | 11mm | | 30 |
| 1207351 | Adapter, 1x0.2/0.4ml D6mm, 6/pk | | 0.2/0.4ml | 6mm | | 30 |
| 1207352 | Adapter, 1x0.5ml D8mm, 6/pk | | 0.5ml | 8mm | | 30 |

Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R, FC5816, FC5816R

| | | | | | | |
|---|---------------------------------|----------------|----------------|--------------------------------|----------------|-------------------|
| Rotor Photo  | Camlab product Code | 1207261 | Description | Rotor, Angle, 24x1.5/2.0ml, ID | | |
| | Compatible with: | FC5714 | FC5718 | FC5718R | FC5816 | FC5816R |
| | Camlab Product Code | 1206875 | 1206876 | 1206877 | 1206878 | 1206879 |
| | max. speed | 14000 rpm | 14000 rpm | 15000 rpm | 15000 rpm | 16000 rpm |
| | max. RCF | 18624 xg | 18624 xg | 21379 xg | 21379 xg | 24325 xg |
| Accessories | | | | | | |
| Camlab product Code | Description | | Vol./Tube | Diameter | Tube type | No.of tubes/rotor |
| No Adapter | Rotor, Angle, 24x1.5/2.0ml, ID | | 1.5/2.0ml | 11mm | | 24 |
| 1207351 | Adapter, 1x0.2/0.4ml D6mm, 6/pk | | 0.2/0.4ml | 6mm | | 24 |
| 1207352 | Adapter, 1x0.5ml D8mm, 6/pk | | 0.5ml | 8mm | | 24 |