

Dry Block Heaters

Models

Model	Number of Blocks	Control	Temperature Range	Stability @ 37°C	Uniformity Within the Block @ 37°C	Item No.	
HB1DG	1	Digital	Ambient + 5°C to 120°C	± 0.2°C	± 0.2°C	30392061	
HB2DG	2	Digital	Ambient + 5°C to 120°C	± 0.2°C	± 0.2°C	30392082	
HB4DG	4	Digital	Ambient + 5°C to 120°C	± 0.2°C	± 0.2°C	30392089	
HB6DG	6	Digital	Ambient + 5°C to 120°C	± 0.2°C	± 0.2°C	30392096	
HB1AL	1	Analog	Low Range: Ambient +5°C to 100°C High Range: 75°C to 150°C	± 1.5°C	± 0.2°C	30392047	
HB2AL	2	Analog		± 2°C	± 0.2°C	30392054	
HB4AL	4	Analog		± 2.5°C	± 0.2°C	30392068	
HB6AL	6	Analog		± 2.5°C	± 0.2°C	30392075	
HB2DGHL	2	Digital	Ambient + 5°C to 100°C	± 0.2°C	± 0.2°C	30392103	

Accessories

Item No.	Description	HB1DG	HB2DG	HB4DG	HB6DG	HB1AL	HB2AL	HB4AL	HB6AL	HB2DGHL	
30400151	Module Block 10 mm Test Tube	X	X	X	X	X	X	X	X	X	
30400152	Module Block 12/13 mm 20 Hole	X	X	X	X	X	X	X	X	X	
30400153	Module Block 15/16 mm 12 Hole	X	X	X	X	X	X	X	X	X	
30400154	Module Block 20 mm Test Tube 8 Hole	X	X	X	X	X	X	X	X	X	
30400155	Module Block 25 mm Test Tube 8 Hole	X	X	X	X	X	X	X	X	X	
30400156	Module Block Combination	X	X	X	X	X	X	X	X	X	
30400157	Module Block 0.5 mL Microtaper	X	X	X	X	X	X	X	X	X	
30400158	Module Block 6 mm Test Tube	X	X	X	X	X	X	X	X	X	
30400159	Module Block 1.5 mL Microtaper	X	X	X	X	X	X	X	X	X	
30400160	Module Block Solid	X	X	X	X	X	X	X	X	X	
30400161	Module Block Cuvette - 12 Cuvet	X	X	X	X	X	X	X	X	X	
30400162	Module Block 11.5/1.5 mL Microcentrifuge	X	X	X	X	X	X	X	X	X	
30400163	Module Block 250 Microliter/6 mm	X	X	X	X	X	X	X	X	X	
30400165	Module Block 12/13 mm 16 Holes	X	X	X	X	X	X	X	X	X	
30400167	Module Block 4 X 35 mm Test Tube	X	X	X	X	X	X	X	X	X	
30400168	Module Block 50 mL Conical Centrifuge	X	X	X	X	X	X	X	X	X	
30400169	Module Block 0.2 mL Micro 80 Hole	X	X	X	X	X	X	X	X	X	
30400170	Module Block 0.2 mL Micro 64 Hole	X	X	X	X	X	X	X	X	X	
30400172	Module Block 15 mL Conical Centrifuge	X	X	X	X	X	X	X	X	X	
30400164	Module Block Titer Plate		X	X	X		X	X	X	X	
30400166	Module Block Solid Double		X	X	X		X	X	X	X	
30400171	Module Block 0.2 mL 96 Hole Double Block		X	X	X		X	X	X	X	
30400173	Sand Bath, 1 Block Unit	X				X					
30400174	Sand Bath, 2 Block Unit		X				X			X	
30400175	Sand Bath, 4 Block Unit			X				X			
30400176	Sand Bath, 6 Block Unit				X				X		
30400177	Sand, 0.5 kg (1 Lb) Bag	X	X	X	X	X	X	X	X	X	
30400178	Stainless Steel Shot, 0.5 kg (1 Lb)	X	X	X	X	X	X	X	X	X	
30400179	Cover Low-Temperature 2 Block	X	X	X	X	X	X	X	X	X	
30400180	Cover Low-Temperature 4 Block	X	X	X	X	X	X	X	X	X	
30400181	Cover Low-Temperature 6 Block	X	X	X	X	X	X	X	X	X	

Dry Block Heaters

Accessories

Item No.	Description	HB1DG	HB2DG	HB4DG	HB6DG	HB1AL	HB2AL	HB4AL	HB6AL	HB2DGHL	
30400182	Module Block For Vials 12 mm	X	X	X	X	X	X	X	X	X	
30400183	Module Block For Vials 15 mm	X	X	X	X	X	X	X	X	X	
30400184	Module Block For Vials 17 mm	X	X	X	X	X	X	X	X	X	
30400185	Module Block For Vials 19 mm	X	X	X	X	X	X	X	X	X	
30400186	Module Block For Vials 21 mm	X	X	X	X	X	X	X	X	X	
30400187	Module Block For Vials 23 mm	X	X	X	X	X	X	X	X	X	
30400188	Module Block For Vials 25 mm	X	X	X	X	X	X	X	X	X	
30400189	Module Block For Vials 28 mm	X	X	X	X	X	X	X	X	X	
30400190	Module Block For Vials 16 mm	X	X	X	X	X	X	X	X	X	
30400191	Module Block, 2.0 mL Tube	X	X	X	X	X	X	X	X	X	
30400192	Module Block, 2.0 mL Corning Tube	X	X	X	X	X	X	X	X	X	
30400193	Module Block Combo, 50 mL, 15 mL, 1.5 mL	X	X	X	X	X	X	X	X	X	
30400194	Module Block Combo - Micro Tube	X	X	X	X	X	X	X	X	X	
30400195	Module Block 12 Holes 17-18 mm	X	X	X	X	X	X	X	X	X	
30400246	Stainless Steel Probe Kit	X	X	X	X	X	X	X	X	X	