



The Homogenizer Company™



*Ideally suited
for nucleic acid, drug,
protein, enzyme,
and virus extractions.*

Bead Ruptor 4

BEAD MILL HOMOGENIZER

**Compelling performance
at a compelling price**

Bead Ruptor 4

BEAD MILL HOMOGENIZER

ORDER NO. 25-001



Processing Range:
25 µL to 6 mL

Power Rating:
100-230 VAC/50-60 Hz

Dimensions:
8.5" W x 11.5" D x 10" H

Weight: 16.5 kg (54 lbs.)

Run Time: 1 - 5 mins.

Programs: 1

Cycles: 1

Performance Range:
30 g - 150 g in
increments of 1-5 m/s

Efficiently lyse:


- > Human tissues
- > Animal tissues
- > Plant & food
- > Micro-organisms

Omni Bead Ruptor 4 is ideally suited for:

- > Homogenizing up to 4 samples in seconds
- > Reproducible extraction of proteins and nucleic acids
- > Simultaneously processing 4 x 0.5 mL, 4 x 1.5 mL, and 1 x 7 mL samples

*includes Tube Carriage for 1.5 mL microtubes,
2 mL screw cap tubes and 7 mL tubes*

- > Holds 4 x 0.5 mL / 1.5 mL / 2 mL or 1 x 7 mL tubes
- > Includes semi-permanent ring and assembly tools



**Patent
Pending
Snap In
Tube Carriage**



*For best
results, use
with Omni
pre-filled
bead kits*

Distributed and Supported in the Uk by:



24 Norman Way Industrial Estate, Over, Cambridge, CB24 5WE
Tel: +44(0)1954 233 100 Fax: +44(0)1954 233 101
Email: sales@camlab.co.uk Web: www.camlab.co.uk