

# Pharmacy Refrigerators

Porkka pharmacy cabinets are designed for the storage of temperature critical medicines and drugs.

They are energy efficient, safe and reliable. Supplied with three pull out drawers with adjustable dividers and three stainless steel shelves, for the storage of larger items.

Our standard models are supplied with a glass window and LED illumination, allowing item to be pre-selected before the door is opened, reducing energy consumption.

**R290**  
REFRIGERANT



Distributed and Supported in the UK by:



[www.camlab.co.uk](http://www.camlab.co.uk)  
**camlab**

24 Norman Way Industrial Estate, Over, Cambridge, CB24 5WE  
Tel: +44(0)1954 233 100 Fax: +44(0)1954 233 101  
Email: [sales@camlab.co.uk](mailto:sales@camlab.co.uk) Web: [www.camlab.co.uk](http://www.camlab.co.uk)

Designed to have a long life cycle and operating on R290 the most energy efficient refrigerant, Porkka pharmacy cabinets are both environmentally friendly and cost efficient.

## Features

- Digital temperature display with built-in microprocessor controller
- Keyboard tamper protection
- Recording of max./min. temperature
- Visual and audible high/low temperature alarm
- Door open alarm
- Probe failure alarm
- Digital calibration for temperature probe
- Potential free alarm contact for connection to a Building Management System
- TTL-Mod Bus communication port (RS 485 with optional converter)
- Low temperature protection
- Exterior finish in white polyester-coated sheet steel, interior finish in stainless steel grade AISI 304 (all models also available in finish outside/inside stainless steel grade AISI 304)
- Insulation with high density CFC/HCFC free polyurethane
- 3 pcs telescope drawers and 3 pcs full-size stainless steel plate shelves, adjustable in increments of 30 mm
- Door with glass window, equipped with central lock, magnetic gasket, heavy duty spring hinges and LED light
- Efficient and energy saving, hermetic air cooled CFC/HCFC refrigerant R 290 system.



**PC 520 GW**



**PC 720 GW**



**PC 722 GW**



Telescopic stainless steel drawers

## Optional features

- Extra shelves as required
- Extra Drawers as required
- Wire baskets in lieu of drawers or shelves
- Castors in lieu of legs
- Full glass door or full door in stainless steel AISI 304 in lieu of standard door with glass window
- Real-time temperature data logging with download to USB facility, optional with battery back up
- Converter TTL Mod Bus to RS 485
- Coded lock for limited personal access
- Electronic key lock with access notification of whom has opened the cabinet plus time and date
- Additional thermo mass probe connected to the controller
- 8 mm sleeved access port for remote monitoring.

## Specifications

	PC 520 GW	PC 720 GW	PC 722 GW
Temperature	+2...+6 °C	+2...+6 °C	+2...+6 °C
Ambient Temperature	+5...+38 °C	+5...+38 °C	+5...+38 °C
Width (mm)	600	850	730
Depth (mm)	700	700	800
Height (mm)	2030	2030	2030
Usable Net Volume (Ltrs)	291	471	447
Shelves/Drawers	3/3	3/3	3/3
Shelf sizes (width x depth)	530 x 650	650 x 530	530 x 650
Drawer sizes (width x depth)	360 x 485 x 50	610 x 485 x 50	490 x 585 x 50
Gross Weight (kg)	120	140	150
Connection Power (W)	346	346	346
Voltage	230V/50Hz	230V/50Hz	230V/50Hz